



[FDA Home](#)³ [Medical Devices](#)⁴ [Databases](#)⁵

Class 2 Device Recall BD Pyxis Products

⁶ [510\(k\)](#)⁷ [DeNovo](#)⁸ [Registration & Listing](#)⁹ [Adverse Events](#)¹⁰ [Recalls](#)¹¹ [PMA](#)¹² [HDE](#)¹³ [Classification](#)¹⁴ [Standards](#)¹⁵
[CFR Title 21](#)¹⁶ [Radiation-Emitting Products](#)¹⁷ [X-Ray Assembler](#)¹⁸ [Medsun Reports](#)¹⁹ [CLIA](#)²⁰ [TPLC](#)²¹

[New Search](#)

[Back to Search Results](#)

Class 2 Device Recall BD Pyxis Products



Date Initiated by Firm	May 08, 2025
Date Posted	June 24, 2025
Recall Status ¹	Open ³ , Classified
Recall Number	Z-2017-2025
Recall Event ID	96585 ²³
Product Classification	Cabinet, table and tray, anesthesia ²⁴ - Product Code BRY ²⁵
Product	Product Name [REF]: Supply Server CCE Enterprise SW Site License (2-5) BD Pyxis Logistics [806] BD Pyxis MedStation 4000 Main 10885403512629 [303] BD Pyxis MedStation ES Tower 10885403512674 [352] BD Pyxis Anesthesia Station ES 10885403477836 [327] BD Pyxis MedStation ES (Main) 10885403512667 [323] MD BD Pyxis SupplyStation V9 [317] Security Module [139698-01] BD Knowledge Portal for Medication Technologies [136607-01] STOCKSTN V10.X WRD/WL SCANNERS [132-53-02] PRS STOCKSTN V10.X WRD/WL SCANNERS [132-53-01] STOCKSTN V10.X WIRED SCANNER [132-52-02] PRS STOCKSTN V10.X WIRED SCANNER [132-52-01] BD Pyxis CII Safe ES, Desktop PC [1156-00] BD Pyxis SupplyStation Panel PC V11 [1132-00] BD Pyxis SupplyRoller V11 [1127-00] Dell 640 Server [1128-00] BD Pyxis SupplyStation System V11 [1123-00] BD Pyxis CII Safe ES 10885403512605 [1116-00] BD Pyxis Enterprise Server 10885403518348, 10885403520341 [1115-00] BD Pyxis SupplyRoller V10 [1101-00] BD Pyxis CII Safe, V9.X, SGL IM, BIO 10885403512520 [111-221] BD Pyxis CII Safe, V7.X, DBL INTG MAIN, BIO, SOLID DO [111-194] BD Pyxis CII Safe, V8.X Desktop PC, Bio [107-245-01]
Code Information	Each customer was provided a unique list of their affected devices that included computer names and atlas keys. Distributed: April 14, 2014 - March 21, 2025 REF or Name [UDI-DI]: Supply Server [NA], CCE Enterprise SW Site License (2-5) [NA], 806 [NA],

303 [10885403512629],
 352 [10885403512674],
 327 [10885403477836],
 323 [10885403512667],
 317 [NA],
 139698-01 [NA],
 136607-01 [NA],
 132-53-02 [NA],
 132-53-01 [NA],
 132-52-02 [NA],
 132-52-01[NA] ,
 1156-00 [NA],
 1132-00 [NA],
 1127-0 [NA],
 1128-00 [NA],
 1123-00 [NA],
 1116-00 [NA],
 1115-00 [10885403518348, 10885403520341],
 1101-00 [NA],
 111-221 [10885403512520],
 111-194 [NA],
 107-245-01 [NA],

**Recalling Firm/
Manufacturer**

CareFusion 303, Inc.
 10020 Pacific Mesa Blvd
 San Diego CA 92121-4386

**Manufacturer Reason
for Recall**

During automated dispensing cabinet upgrade/installation/reimaging, component manager was left in installed mode, and if a Microsoft patch is then installed this may lead to an error preventing user access to the dispensing software application, which may lead to delay in access to medications/supplies, and delivery/replenishment of patient specific orders; which may lead to delay in patient care.

**FDA Determined
Cause ²**

Software Design Change

Action

On 5/8/2025, correction notices were mailed and emailed to Customers informing them of the following:

Firm has temporarily configured impacted devices to "disabled" mode so the device will not automatically apply validated patches intended to address security vulnerabilities or performance issues.

Device keys can be used to manually access unavailable drawers and must be readily available and retrievable to constitute risk mitigation; medications can also be obtained from other unaffected dispensing devices or pharmacy.

The Pyxis Anesthesia Station unlocks non-secure drawers upon user login, providing immediate access to medications. Non-secure drawers remain unlocked while the user is logged in.

User guides include instructions on how to manually access medications in an emergency: BD Pyxis MedStation ES user guide DME 10000430667 version 00 (page 120), BD Pyxis Anesthesia Station ES user guide DME 10000430730 version 00 (page 84).

Other potential mitigations:

-Healthcare facility system administrators can grant permission to add temporary patients to select staff per user guide DME 10000430667 version 00 (page 25).

-Critical override or temporary Non-RX mode can be configured for impacted profile devices per user guide DME 10000430667 version 00 (pages 74-75).

-Refill reports can be generated at the device to identify items requiring replenishment if inventory information is not communicated to the server.

Share the recall letter with all facilities in your network. Post the recall letter on or near affected products.

Complete and return response form via email to BDRC22@bd.com

Firm will coordinate with customers to obtain a schedule to implement "Managed Mode".

The following firm contacts were provided: North American Regional Complaint Center: 1-844-823-5433, productcomplaints@bd.com; BD Remediation Support: 1-866-583-8783; BD Technical Support: 1-800-727-6102

Quantity in Commerce

6,313

Distribution

Worldwide - US Nationwide distribution including in the states of PR, DE, GU, MO, IA, IL, CA, MD, MA, VA, PA, NY, WA, NJ, AZ, CO, OH, FL, WI, TX, CT, LA, MI, WY, MN, AR, KY, OR, NM, TN, NC, SC, IN, VT, KS, AL, GA, WV, NV, ID, DC, MS, NE, HI, OK, UT, ND, MT, NH, SD, ME, RI, AK, MP and the countries of SA, CA, NZ, CH, FR, IT, BE, AU, BH, QA, AE, GB, EG, AR, JO, DE, BM, TH, LB, BR, IE, FI, MX, TW, ID, SG, BS, IN, MA, GF, ES, HK.

Total Product Life Cycle

[TPLC Device Report](#)²⁶

¹ A record in this database is created when a firm initiates a correction or removal action. The record is updated if the FDA identifies a violation and classifies the action as a recall, and it is updated for a final time when the recall is terminated. Learn more about [medical device recalls](#)²⁷.

² Per FDA policy, recall cause determinations are subject to modification up to the point of termination of the recall.

³ The manufacturer has initiated the recall and not all products have been corrected or removed. This record will be updated as the status changes.

Links on this page:

1. <http://www.addthis.com/bookmark.php?u508=true&v=152&username=fdomain>
2. <http://www.addthis.com/bookmark.php>
3. <https://www.fda.gov/>
4. <http://www.fda.gov/MedicalDevices/default.htm>
5. <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/Databases/default.htm>
6. </scripts/cdrh/devicesatfda/index.cfm>
7. </scripts/cdrh/cfdocs/cfPMN/pmn.cfm>
8. </scripts/cdrh/cfdocs/cfpmn/denovo.cfm>
9. </scripts/cdrh/cfdocs/cfRL/rl.cfm>
10. </scripts/cdrh/cfdocs/cfMAUDE/TextSearch.cfm>
11. </scripts/cdrh/cfdocs/cfRES/res.cfm>
12. </scripts/cdrh/cfdocs/cfPMA/pma.cfm>
13. </scripts/cdrh/cfdocs/cfHDE/hde.cfm>
14. </scripts/cdrh/cfdocs/cfPCD/classification.cfm>
15. </scripts/cdrh/cfdocs/cfStandards/search.cfm>
16. </scripts/cdrh/cfdocs/cfCFR/CFRSearch.cfm>
17. /scripts/cdrh/cfdocs/cfPCD_RH/classification.cfm
18. </scripts/cdrh/cfdocs/cfAssem/assembler.cfm>
19. </scripts/cdrh/cfdocs/Medsun/searchReportText.cfm>
20. </scripts/cdrh/cfdocs/cfClia/Search.cfm>
21. </scripts/cdrh/cfdocs/cfTPLC/tpic.cfm>
22. <https://www.fda.gov/safety/recalls-market-withdrawals-safety-alerts/enforcement-reports>
23. /scripts/cdrh/cfdocs/cfRES/res.cfm?start_search=1&event_id=96585
24. </scripts/cdrh/cfdocs/cfPCD/classification.cfm?ID=BRY>
25. </scripts/cdrh/cfdocs/cfPCD/classification.cfm?ID=BRY>
26. </scripts/cdrh/cfdocs/cfTPLC/tpic.cfm?id=BRY>
27. <https://www.fda.gov/medical-devices/medical-device-recalls/what-medical-device-recall>

Page Last Updated: 07/17/2025

Note: If you need help accessing information in different file formats, see [Instructions for Downloading Viewers and Players](#).

Language Assistance Available: [Español](#) | [繁體中文](#) | [Tiếng Việt](#) | [한국어](#) | [Tagalog](#) | [Русский](#) | [العربية](#) | [Kreyòl Ayisyen](#) | [Français](#) | [Polski](#) |

[Português](#) | [Italiano](#) | [Deutsch](#) | [日本語](#) | [فارسی](#) | [English](#)

[Accessibility](#) [Contact](#) [FDA Careers](#) [FDA Basics](#) [FOIA](#) [No FEAR Act](#) [Nondiscrimination](#) [Website Policies](#) / [Privacy](#)



U.S. Food and Drug Administration
10903 New Hampshire Avenue
Silver Spring, MD 20993
Ph. 1-888-INFO-FDA (1-888-463-6332)
[Contact FDA](#)



[For Government For Press](#)

[Combination Products](#) [Advisory Committees](#) [Science & Research](#) [Regulatory Information](#) [Safety](#) [Emergency](#) [Preparedness](#)
[International Programs](#) [News & Events](#) [Training and Continuing Education](#) [Inspections/Compliance](#) [State & Local Officials](#)
[Consumers](#) [Industry](#) [Health Professionals](#) [FDA Archive](#) [Vulnerability](#) [Disclosure Policy](#)



U.S. Department of **Health & Human Services**

Links on this page:

1. <http://www.addthis.com/bookmark.php?u508=true&v=152&username=fdomain>
2. <http://www.addthis.com/bookmark.php>
3. <https://www.fda.gov/>
4. <http://www.fda.gov/MedicalDevices/default.htm>
5. <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/Databases/default.htm>
6. </scripts/cdrh/devicesatfda/index.cfm>
7. </scripts/cdrh/cfdocs/cfPMN/pmn.cfm>
8. </scripts/cdrh/cfdocs/cfpmn/denovo.cfm>
9. </scripts/cdrh/cfdocs/cfRL/rl.cfm>
10. </scripts/cdrh/cfdocs/cfMAUDE/TextSearch.cfm>
11. </scripts/cdrh/cfdocs/cfRES/res.cfm>
12. </scripts/cdrh/cfdocs/cfPMA/pma.cfm>
13. </scripts/cdrh/cfdocs/cfHDE/hde.cfm>
14. </scripts/cdrh/cfdocs/cfPCD/classification.cfm>
15. </scripts/cdrh/cfdocs/cfStandards/search.cfm>
16. </scripts/cdrh/cfdocs/cfCFR/CFRSearch.cfm>
17. /scripts/cdrh/cfdocs/cfPCD_RH/classification.cfm
18. </scripts/cdrh/cfdocs/cfAssem/assembler.cfm>
19. </scripts/cdrh/cfdocs/Medsun/searchReportText.cfm>
20. </scripts/cdrh/cfdocs/cfClia/Search.cfm>
21. </scripts/cdrh/cfdocs/cfTPLC/tpic.cfm>
22. <https://www.fda.gov/safety/recalls-market-withdrawals-safety-alerts/enforcement-reports>
23. /scripts/cdrh/cfdocs/cfRES/res.cfm?start_search=1&event_id=96585
24. </scripts/cdrh/cfdocs/cfPCD/classification.cfm?ID=BRY>
25. </scripts/cdrh/cfdocs/cfPCD/classification.cfm?ID=BRY>
26. </scripts/cdrh/cfdocs/cfTPLC/tpic.cfm?id=BRY>
27. <https://www.fda.gov/medical-devices/medical-device-recalls/what-medical-device-recall>