

PRODUCT CORRECTION

Cardinal Health
3651 Birchwood Dr
Waukegan, IL 60085
800.292.9332 toll free
614-652-9648 fax



CardinalHealth

Risk Management Department

cardinalhealth.com

SEP 09 2025

August 14, 2025

Received

Description: Determine™ HIV-1/2 Ag/Ab COMBO			
Cardinal Health Catalog Numbers	Abbott Rapid DX Catalog Numbers	Lot Numbers	
B7D2648	7D2648	0000851700	0000930923
B7D2648BD		0000867175	0000939868
B7D2648BJB		0000869145	0000946267
B7D2648IH		0000871699	0000950115
B7D2648JM		0000871729	0000950116
B7D2648MM		0000874643	0000952264
B7D2648NDX		0000927861	0000952265
B7D2648OH		0000927862	0000984187
B7D2648PSJ		0000927863	0000990023
B7D2648SB		0000929961	0000962456
B7D2648UAB		0000930922	0000962457

Event-2025-07216

Dear Recall Coordinator/Materials Management/Purchasing:

Cardinal Health has recently been informed that **Abbott Rapid DX** is initiating a voluntary product correction due to the impacted lots having a higher occurrence of preliminary false reactive complaints due to a subcomponent used in the manufacture of Determine™ HIV-1/2 Ag/Ab COMBO. Please see enclosed **Abbott Rapid DX** Product Correction Notice.

Please take the following action:

1. Count how many kit boxes (25 tests / kit box) of remaining impacted inventory require replacement.
2. Discontinue use of and destroy any remaining inventory of the impacted lots according to your procedures.
3. Complete and return the Customer Acknowledgment Form (Form must be completed, signed and returned to receive replacement product). Even if you do not have impacted inventory in stock. Send completed form to ssc.ardxfieldactions@abbott.com
4. Complete the Customer Request for Replacement Product Form and submit via email ssc.ardxfieldactions@abbott.com to order replacement product.
5. If you have forwarded the product listed above to others in your network then inform them of this Field Product Correction and provide to them a copy of this notice and request they take the necessary action to quarantine and destroy any remaining inventory of impacted lots and RETURN THE FORM TO YOU to coordinate product replacement on their behalf.

PRODUCT CORRECTION

Cardinal Health
3651 Birchwood Dr
Waukegan, IL 60085
800.292.9332 toll free
614-652-9648 fax



CardinalHealth

cardinalhealth.com

6. Retain this letter for your records.
7. All product lots not identified in Appendix I can continue to be used.

We sincerely apologize for any inconvenience this notice may have caused you and your staff. If you or any of the health care providers you serve have questions completing the form or product replacement, please feel free to contact Abbott Rapid Diagnostics Customer Service at 1-877-441-7440, Option 1. (Available 8am – 5pm EST, Monday – Friday)

Sincerely,

Eric Kastory

Sr. Specialist, Quality Compliance



Product Correction Notice

Action Required

Date Issued AUGUST 2025

Product

Product Description	List Number	Lot Number	GTIN
Determine™ HIV-1/2 Ag/Ab COMBO	7D2648	Various – Reference Appendix I	00811877011101

Explanation

Dear Valued Customer,

The purpose of this letter is to inform you of a Field Product Correction impacting select lots of Determine™ HIV-1/2 Ag/Ab COMBO. Abbott has confirmed that the impacted lots identified in Appendix I of this document have a higher occurrence of preliminary false reactive complaints due to a subcomponent used in the manufacture of Determine™ HIV-1/2 Ag/Ab COMBO. The issue has been isolated to specific lots of the subcomponent that was used to manufacture the kit lots listed in Appendix I.

Please review the steps below which provide details on the actions required by you.

Abbott has worked diligently to determine the root cause of this event, and the necessary actions have been taken to prevent recurrence. We remain committed to providing you with the highest quality diagnostic products and support services to meet your needs.

Patient Impact

- There is a potential for unnecessary prophylactic treatment due to preliminary false reactive results.
- Any reactive test results using Determine HIV-1/2 Ag/Ab COMBO result should be confirmed by additional testing using other tests.

Please follow one of the following actions, as applicable.

**Necessary
Actions to be
Taken by
Customer**

If....	Then...
You have impacted inventory in current stock:	• Count how many kit boxes (25 tests / kit box) of remaining impacted inventory require replacement. • Discontinue use of and destroy any remaining inventory of the impacted lots according to your procedures. • Complete and return the Customer Acknowledgment Form (Form must be completed, signed and returned to receive replacement product). • Complete the Customer Request for Replacement Product Form. • Contact Customer Support to order replacement product. • Retain this letter for your records. • All product lots not identified in Appendix I can continue to be used.
You have forwarded the product listed above to others in your network	• Inform them of this Field Product Correction and provide to them a copy of this notice and request they take the necessary action to quarantine and destroy any remaining inventory of impacted lots and <u>RETURN THE FORM TO YOU</u> to coordinate product replacement on their behalf.
You do not have impacted inventory in stock	• Complete and return the Customer Acknowledgment Form acknowledging you do not have affected inventory on hand, in stock. • All product lots not identified in Appendix I can continue to be used.

**Contact
Information**

If you or any of the health care providers you serve have questions completing the form or product replacement, please contact Abbott Rapid Diagnostics Customer Service at 1-877-441-7440, Option 1. (available 8am – 5pm EST, Monday - Friday).

Adverse reactions or quality problems experienced with the use of this product may be reported to the FDA's MedWatch Adverse Event Reporting program either online (<http://www.fda.gov/MedWatch/report.htm>), by mail (<http://www.fda.gov/MedWatch/getforms.htm>), by phone (1-800-332-1088), or by fax (1-800-FDA-0178).

If you have experienced any patient or user injury associated with this Field Product Correction, please immediately report the event to Technical Service at 877-441-7440, Option 2.

Appendix I

Determine™ HIV-1/2 Ag/Ab COMBO				
Part #	Kit Lot #	Device Lot #	Expiry Date	UDI
7D2648	851700	0000850468	2025-12-04	010081187701110110851700
7D2648	867175	0000853194	2026-02-28	010081187701110110867175
7D2648	869145	0000866249	2026-02-28	010081187701110110869145
7D2648	871699	0000853194	2026-02-28	010081187701110110871699
7D2648	871729	0000866252	2026-02-28	010081187701110110871729
7D2648	874643	0000866255	2026-02-28	010081187701110110874643
7D2648	0000927861	0000912600	2026-05-28	0100811877011101100000927861
7D2648	0000927862	0000912601	2026-05-28	0100811877011101100000927862
7D2648	0000927863	0000926277	2026-07-28	0100811877011101100000927863
7D2648	0000929961	0000926278	2026-07-28	0100811877011101100000929961
7D2648	0000930922	0000943284	2026-07-28	0100811877011101100000930922
7D2648	0000930923	0000938702	2026-09-28	0100811877011101100000930923
7D2648	0000939868	0000944760	2026-09-28	0100811877011101100000939868
7D2648	0000946267	0000945409	2026-09-28	0100811877011101100000946267
7D2648	0000950115	0000945412	2026-09-28	0100811877011101100000950115
7D2648	0000950116	0000949678	2026-10-10	0100811877011101100000950116
7D2648	0000952264	0000945414	2026-10-28	0100811877011101100000952264
7D2648	0000952265	0000938700	2026-10-03	0100811877011101100000952265
7D2648	0000984187	0000969031	2026-03-29	0100811877011101100000984187
7D2648	0000990023	0000969046	2026-03-29	0100811877011101100000990023
7D2648	0000962456	0000952951	2026-10-28	0100811877011101100000962456
7D2648	0000962457	0000966007	2026-10-28	0100811877011101100000962457



Product Correction Notice

Customer Acknowledgement Form Appendix II

Determine™ HIV-1/2 Ag/Ab COMBO

Identifier: 2025 07

This response form is to confirm receipt of this notification and to request replacement product. If eligible, please only return Appendix II (pages 4 and 5 of this notice).

1. Customer Details (all fields are mandatory unless otherwise noted)

Account / Customer Number (optional)	
Healthcare Organization Name	
Distributor Name (optional)	
Street	
City	
State / Zip Code	/
Country	United States
Contact Name	
Title or function	
Telephone number	
Fax number	
E-mail	
Shipping Address (if different than above)	

2. Customer action undertaken

Confirm and complete the actions taken (BOX 1 & 2 REQUIRED FOR REPLACEMENT)	
☐	I confirm that we, the Customer, have received, read, and understood this Product Correction Notice for Determine™ HIV-1/2 Ag/Ab Combo
☐	We have taken the necessary actions as directed by this Product Correction Notice and have destroyed the quantity of affected product that we have outlined in the table below where product still existed in our inventory.
Please Select one (1) of the options below:	
☐	I have affected product. Please complete Request for Replacement Product (next page).
☐	I do not have affected product. (Request for Replacement Product is not required) Print and sign below.
Printed Name	Signature / Date



Product Correction Notice

Customer Request for Replacement Product Form

Determine™ HIV-1/2 Ag/Ab COMBO

Identifier: 2025 07

This response form is to confirm receipt of this notification and to request replacement product. If eligible, please only return Appendix II (pages 4 and 5 of this notice).

Kit Name	Part number	Lot numbers	Requested KIT BOX (25 tests/kit box) Replacement Qty
Determine™ HIV-1/2 Ag/Ab Combo NOTE: Please provide quantity of remaining kits on hand. - One kit contains 25 tests. - If you have a partial kit remaining, round to the nearest full kit.	7D2648	851700	_____ kits or ☐ N/A
	7D2648	867175	_____ kits or ☐ N/A
	7D2648	869145	_____ kits or ☐ N/A
	7D2648	871699	_____ kits or ☐ N/A
	7D2648	871729	_____ kits or ☐ N/A
	7D2648	874643	_____ kits or ☐ N/A
	7D2648	0000927861	_____ kits or ☐ N/A
	7D2648	0000927862	_____ kits or ☐ N/A
	7D2648	0000927863	_____ kits or ☐ N/A
	7D2648	0000929961	_____ kits or ☐ N/A
	7D2648	0000930922	_____ kits or ☐ N/A
	7D2648	0000930923	_____ kits or ☐ N/A
	7D2648	0000939868	_____ kits or ☐ N/A
	7D2648	0000946267	_____ kits or ☐ N/A
	7D2648	0000950115	_____ kits or ☐ N/A
	7D2648	0000950116	_____ kits or ☐ N/A
	7D2648	0000952264	_____ kits or ☐ N/A
	7D2648	0000952265	_____ kits or ☐ N/A
	7D2648	0000984187	_____ kits or ☐ N/A
	7D2648	0000990023	_____ kits or ☐ N/A
	7D2648	0000962456	_____ kits or ☐ N/A
	7D2648	0000962457	_____ kits or ☐ N/A
Printed Name		Signature / Date	

3. Return acknowledgement to sender.

Email	ssc.ardxfieldactions@abbott.com
Deadline for returning the customer reply form	Please complete and return this form within 10 business days of receipt.

By signing this Customer Acknowledgement Form, Customer hereby waives, fully releases, and forever discharges Abbott Laboratories and its affiliates, parents, as well as any agents thereof, from any and all claims, liabilities, damages, and causes of action, of any nature whatsoever, arising from or relating to the product identified on this Form.

ORIGIN ID:PHDA (800) 292-9332
REGULATORY MANAGEMENT
CARDINAL HEALTH
3651 BIRCHWOOD

SHIP DATE: 19AUG25
ACTWGT: 0.50 LB MAN
CAD: 0397361/CAFE3908

WAUKEGAN, IL 60086728
UNITED STATES US

BILL THIRD PARTY

TO ATTN RECALL COORDINATOR/MATL MGMT
THE MOUNT SINAI HOSPITAL
20 E 101ST ST

NEW YORK NY 100296538

PO EVENT:2025-07216-ABBOTT

REF: 664



J251024111401uv

TRK#
0201 4629 6253 4449

WED - 20 AUG 10:30A
PRIORITY OVERNIGHT

XG FLUA

10029
NY-US EWR



58CC4/A855/FE2D

