



SECOND NOTICE
URGENT – MEDICAL DEVICE RECALL
3M™ Prevena™ Plus 125 Therapy Units and Kits
Distributed Between March 11, 2024, and December 18, 2024

April 9, 2025

Dear Valued Customer,

KCI USA, Inc., a Solventum company ("Solventum"), values the care of our patients and the quality of our products. This letter is to inform you that Solventum is initiating a voluntary medical device recall on 3M™ Prevena™ Plus 125 Therapy Units and Kits distributed between March 11, 2024, and December 18, 2024. This is the second notice of this voluntary recall.

Solventum has identified that devices distributed between this period were released with a leak alarm threshold lower than the design specification. This results in the therapy unit prematurely alarming although there is no leak within the system and potentially continuing to alarm during therapy although there is no visible leak.

This medical device recall is to inform customers who purchased the affected products of this issue to cease use and quarantine the impacted devices and to provide instructions on how to return devices for replacement.

Solventum has received 342 complaints for this issue, but zero (0) reports of injury have been reported. There are no risks to patients. The devices continue to provide negative pressure wound therapy, but due to the lower leak alarm threshold, nuisance alarms may be triggered and continue during therapy although there are no visible leaks.

Our records indicate that you have purchased one or more products affected by this voluntary field corrective action.

Actions Required by Your Facility:

1. Screen your facility or other storage locations for the impacted products and lot numbers listed in Attachment 1.
2. Immediately cease use and quarantine any identified product.
3. Complete the Field Notice Receipt and Customer Response Form acknowledging you have received a copy of this notice and select one of two options:
 - a. Acknowledge you do not have any product impacted by the field action.
 - b. Acknowledge you have product impacted by the field action and complete the catalog number, lot numbers and quantities you have identified and quarantined.
4. Return the form via email to kci3mfieldactionresponse@solventum.com.
5. Contact your sales representative or Solventum Customer Service at 1-800-275-4524 to help coordinate the return shipment and receive return labels. Products being returned to Solventum must be in an outer corrugated box with extra space filled to protect the individual cartons.
6. Products returned to Solventum will be assessed, reprogrammed, and returned or replaced.
 - a. Product SKUs and quantities shipped to Solventum will be returned or replaced with the same quantities.
 - b. Returned or replaced product will have a green dot visual identifier on the product label confirming the assessment and reprogramming of the product.
 - c. Returned or replaced product may have different lot numbers and expiry dates; Solventum plans to ship available product to help reduce wait times.



This Medical Device Recall must be distributed to all individuals who need to be aware within your organization and to any organization where the potentially affected devices have been transferred or sold. Please notify any facilities that you may have transferred devices to by using a letter and acknowledgement form.

Adverse reactions or quality problems experienced with the use of this product may be reported to the FDA's MedWatch Adverse Event Reporting program either online, by regular mail or by fax.

- Complete and submit the report Online.
- Regular Mail or Fax: Download form or call 1- 800-332-1088 to request a reporting form, then complete and return to the address on the pre-addressed form, or submit by fax to 1-800-FDA-0178

Solventum takes the quality of our products seriously, and we strive to meet or exceed our customers' expectations. We apologize for any inconvenience and greatly appreciate your understanding and cooperation as we take this voluntary action.

If you have any questions, please contact your local Solventum representative, visit Solventum.com or call 800-275-4524.

Sincerely,

A handwritten signature in cursive script that reads 'Harry Janiski'.

Harry Janiski
Vice President, Global Quality Assurance

Attachment 1: List of Products and Lots Affected – Page 1 of 6

CATALOG NUMBER	ITEM DESCRIPTION	LOT NUMBER
PRE3021USV005	Prevena Plus Duo Peel and Place 20cm System Kit, US	C16793V005
PRE3021USV005	Prevena Plus Duo Peel and Place 20cm System Kit, US	C17567V005
PRE3021USV005	Prevena Plus Duo Peel and Place 20cm System Kit, US	C18503V005
PRE3201USV006	Prevena Plus Peel and Place 35 cm system Kit, US	C16794V006
PRE3201USV006	Prevena Plus Peel and Place 35 cm system Kit, US	C17165V006
PRE3201USV006	Prevena Plus Peel and Place 35 cm system Kit, US	C17569V006
PRE3201USV006	Prevena Plus Peel and Place 35 cm system Kit, US	C17954V006
PRE3201USV006	Prevena Plus Peel and Place 35 cm system Kit, US	C18502V006
PRE3201USV006	Prevena Plus Peel and Place 35 cm system Kit, US	C18610V006
PRE3201USV006	Prevena Plus Peel and Place 35 cm system Kit, US	C19238V006
PRE3201USV007	Prevena Plus Peel and Place 35 cm system Kit, US	C18068V007
PRE3201USV007	Prevena Plus Peel and Place 35 cm system Kit, US	C18504V007
PRE3201V008	PREVENA PLUS System with PEEL & PLACE Dressing 35cm	C17273V008
PRE3201V008	PREVENA PLUS System with PEEL & PLACE Dressing 35cm	C17882V008
PRE3201V008	PREVENA PLUS System with PEEL & PLACE Dressing 35cm	C18391V008
PRE3201V008	PREVENA PLUS System with PEEL & PLACE Dressing 35cm	C19240V008
PRE3201V009	PREVENA PLUS System with PEEL & PLACE Dressing 35cm	C16796V009
PRE3201V009	PREVENA PLUS System with PEEL & PLACE Dressing 35cm	C17755V009
PRE3201V009	PREVENA PLUS System with PEEL & PLACE Dressing 35cm	C17955V009
PRE3201V010	PREVENA PLUS System with PEEL & PLACE Dressing 35cm	C16797V010
PRE3201V010	PREVENA PLUS System with PEEL & PLACE Dressing 35cm	C17956V010
PRE3201V010	PREVENA PLUS System with PEEL & PLACE Dressing 35cm	C18505V010
PRE3201V011	PREVENA PLUS System with PEEL & PLACE Dressing 35cm	C16798V011
PRE3201V011	PREVENA PLUS System with PEEL & PLACE Dressing 35cm	C17957V011
PRE3201V011	PREVENA PLUS System with PEEL & PLACE Dressing 35cm	C18506V011
PRE4000USV005	Prevena Plus 125 Therapy Unit, US	C17175V005
PRE4000USV005	Prevena Plus 125 Therapy Unit, US	C17884V005
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C16801V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C16874V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C16998V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C17030V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C17133V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C17183V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C17405V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C17541V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C17574V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C17757V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C17779V008

Attachment 1: List of Products and Lots Affected – Page 2 of 6

PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C17878V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C17951V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C18057V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C18149V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C18208V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C18392V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C18444V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C18719V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C18723V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C18810V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C18860V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C18956V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C19159V008
PRE4000USV008	PREVENA PLUS STANDALONE SYSTEM KIT, US	C19232V008
PRE4001AUV006	PREVENA PLUS SYSTEM KIT, AU	C16680V006
PRE4001AUV006	PREVENA PLUS SYSTEM KIT, AU	C17886V006
PRE4001AUV006	PREVENA PLUS SYSTEM KIT, AU	C18461V006
PRE4001AUV006	PREVENA PLUS SYSTEM KIT, AU	C18813V006
PRE4001CAV006	PREVENA PLUS SYSTEM KIT, CA	C17683V006
PRE4001CAV006	PREVENA PLUS SYSTEM KIT, CA	C17818V006
PRE4001CAV006	PREVENA PLUS SYSTEM KIT, CA	C19162V006
PRE4001UKV006	PREVENA PLUS SYSTEM KIT, UK	C17024V006
PRE4001UKV006	PREVENA PLUS SYSTEM KIT, UK	C17177V006
PRE4001UKV006	PREVENA PLUS SYSTEM KIT, UK	C17300V006
PRE4001UKV006	PREVENA PLUS SYSTEM KIT, UK	C17958V006
PRE4001USV005	PREVENA PLUS SYSTEM KIT, US	C18959V005
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C16753V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C16802V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C16875V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C16999V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C17032V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C17129V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C17184V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C17258V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C17303V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C17406V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C17565V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C17576V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C17780V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C17879V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C17952V006

Attachment 1: List of Products and Lots Affected – Page 3 of 6

PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C18058V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C18205V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C18227V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C18394V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C18466V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C18721V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C18725V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C18812V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C18862V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C18958V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C19060V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C19161V006
PRE4001USV006	PREVENA PLUS SYSTEM KIT, US	C19234V006
PRE4001V006	PREVENA PLUS SYSTEM KIT	C17178V006
PRE4001V006	PREVENA PLUS SYSTEM KIT	C18062V006
PRE4001V006	PREVENA PLUS SYSTEM KIT	C18814V006
PRE4001V008	PREVENA PLUS SYSTEM KIT	C18523V008
PRE4001ZAV006	PREVENA PLUS SYSTEM KIT, ZA	C18060V006
PRE4001ZAV006	PREVENA PLUS SYSTEM KIT, ZA	C18464V006
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C16715V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C16883V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17000V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17031V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17134V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17195V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17257V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17302V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17407V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17542V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17575V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17758V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17781V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17880V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17953V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18059V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18150V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18209V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18393V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18465V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18720V001

Attachment 1: List of Products and Lots Affected – Page 4 of 6

PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18724V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18811V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18861V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18957V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C19059V001
PRE4010V001	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C19160V001
PRE4010V003	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C16878V003
PRE4010V004	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C16879V004
PRE4010V004	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17684V004
PRE4010V005	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17007V005
PRE4010V005	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17534V005
PRE4010V005	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17959V005
PRE4010V007	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18815V007
PRE4010V016	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C16880V016
PRE4010V016	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17179V016
PRE4010V016	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17685V016
PRE4010V016	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17960V016
PRE4010V016	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18445V016
PRE4010V016	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18816V016
PRE4010V017	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C16881V017
PRE4010V017	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17180V017
PRE4010V017	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17304V017
PRE4010V017	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17686V017
PRE4010V017	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17961V017
PRE4010V017	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18446V017
PRE4010V017	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18817V017
PRE4010V018	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18077V018
PRE4010V018	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18863V018
PRE4010V019	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17181V019
PRE4010V019	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17687V019
PRE4010V019	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18447V019
PRE4010V020	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C16882V020
PRE4010V020	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17182V020
PRE4010V020	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17307V020
PRE4010V020	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17688V020
PRE4010V020	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18061V020
PRE4010V020	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18448V020
PRE4010V020	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18818V020
PRE4010V021	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C17305V021
PRE4010V021	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18819V021

Attachment 1: List of Products and Lots Affected – Page 5 of 6

PRE4010V023	PREVENA PLUS STANDALONE EXTENDED LIFE SYSTEM KIT	C18820V023
PRE5001V001	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 33X30 CM	C16928V001
PRE5001V001	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 33X30 CM	C17196V001
PRE5001V001	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 33X30 CM	C17759V001
PRE5001V001	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 33X30 CM	C18578V001
PRE5001V001	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 33X30 CM	C18625V001
PRE5001V001	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 33X30 CM	C19061V001
PRE5001V002	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 33X30 CM	C16748V002
PRE5001V002	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 33X30 CM	C17760V002
PRE5001V003	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 33X30 CM	C16749V003
PRE5001V003	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 33X30 CM	C16929V003
PRE5001V004	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 33X30 CM	C16750V004
PRE5001V004	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 33X30 CM	C18626V004
PRE5001V004	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 33X30 CM	C19062V004
PRE5001V005	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 33X30 CM	C16930V005
PRE5001V005	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 33X30 CM	C17761V005
PRE5001V005	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 33X30 CM	C18069V005
PRE5101V001	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 46X30 CM	C17762V001
PRE5101V002	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 46X30 CM	C17763V002
PRE5101V003	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 46X30 CM	C17219V003
PRE5101V004	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 46X30 CM	C17765V004
PRE5101V004	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 46X30 CM	C18070V004
PRE5101V005	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 46X30 CM	C18072V005
PRE5101V005	PREVENA RESTOR ARTHRO-FORM SYSTEM KIT 46X30 CM	C18627V005
PRE5221V001	PREVENA RESTOR BELLA-FORM SYSTEM KIT 21X19 CM	C18579V001
PRE5221V001	PREVENA RESTOR BELLA-FORM SYSTEM KIT 21X19 CM	C18628V001
PRE5221V001	PREVENA RESTOR BELLA-FORM SYSTEM KIT 21X19 CM	C19064V001
PRE5221V002	PREVENA RESTOR BELLA-FORM SYSTEM KIT 21X19 CM	C16931V002
PRE5221V002	PREVENA RESTOR BELLA-FORM SYSTEM KIT 21X19 CM	C17220V002
PRE5221V004	PREVENA RESTOR BELLA-FORM SYSTEM KIT 21X19 CM	C18580V004
PRE5221V004	PREVENA RESTOR BELLA-FORM SYSTEM KIT 21X19 CM	C18782V004
PRE5221V004	PREVENA RESTOR BELLA-FORM SYSTEM KIT 21X19 CM	C19066V004
PRE5221V005	PREVENA RESTOR BELLA-FORM SYSTEM KIT 21X19 CM	C17221V005
PRE5221V005	PREVENA RESTOR BELLA-FORM SYSTEM KIT 21X19 CM	C17767V005
PRE5221V005	PREVENA RESTOR BELLA-FORM SYSTEM KIT 21X19 CM	C18783V005
PRE5321V001	PREVENA RESTOR BELLA-FORM SYSTEM KIT 24X22 CM	C18612V001
PRE5321V001	PREVENA RESTOR BELLA-FORM SYSTEM KIT 24X22 CM	C18803V001
PRE5321V002	PREVENA RESTOR BELLA-FORM SYSTEM KIT 24X22 CM	C18613V002
PRE5321V002	PREVENA RESTOR BELLA-FORM SYSTEM KIT 24X22 CM	C18613V002
PRE5321V003	PREVENA RESTOR BELLA-FORM SYSTEM KIT 24X22 CM	C16751V003

Attachment 1: List of Products and Lots Affected – Page 6 of 6

PRE5321V003	PREVENA RESTOR BELLA-FORM SYSTEM KIT 24X22 CM	C18614V003
PRE5321V004	PREVENA RESTOR BELLA-FORM SYSTEM KIT 24X22 CM	C18615V004
PRE5321V004	PREVENA RESTOR BELLA-FORM SYSTEM KIT 24X22 CM	C19538V004
PRE5321V005	PREVENA RESTOR BELLA-FORM SYSTEM KIT 24X22 CM	C18616V005
PRE5421V001	PREVENA RESTOR BELLA-FORM SYSTEM KIT 29X27 CM	C16932V001
PRE5421V001	PREVENA RESTOR BELLA-FORM SYSTEM KIT 29X27 CM	C17222V001
PRE5421V001	PREVENA RESTOR BELLA-FORM SYSTEM KIT 29X27 CM	C17768V001
PRE5421V001	PREVENA RESTOR BELLA-FORM SYSTEM KIT 29X27 CM	C18073V001
PRE5421V001	PREVENA RESTOR BELLA-FORM SYSTEM KIT 29X27 CM	C18581V001
PRE5421V001	PREVENA RESTOR BELLA-FORM SYSTEM KIT 29X27 CM	C18784V001
PRE5421V001	PREVENA RESTOR BELLA-FORM SYSTEM KIT 29X27 CM	C19067V001
PRE5421V002	PREVENA RESTOR BELLA-FORM SYSTEM KIT 29X27 CM	C16933V002
PRE5421V002	PREVENA RESTOR BELLA-FORM SYSTEM KIT 29X27 CM	C17223V002
PRE5421V003	PREVENA RESTOR BELLA-FORM SYSTEM KIT 29X27 CM	C17224V003
PRE5421V003	PREVENA RESTOR BELLA-FORM SYSTEM KIT 29X27 CM	C17770V003
PRE5421V003	PREVENA RESTOR BELLA-FORM SYSTEM KIT 29X27 CM	C18074V003
PRE5421V003	PREVENA RESTOR BELLA-FORM SYSTEM KIT 29X27 CM	C18583V003
PRE5421V004	PREVENA RESTOR BELLA-FORM SYSTEM KIT 29X27 CM	C16752V004
PRE5421V004	PREVENA RESTOR BELLA-FORM SYSTEM KIT 29X27 CM	C16934V004
PRE5421V004	PREVENA RESTOR BELLA-FORM SYSTEM KIT 29X27 CM	C17225V004
PRE5421V004	PREVENA RESTOR BELLA-FORM SYSTEM KIT 29X27 CM	C17771V004
PRE5421V004	PREVENA RESTOR BELLA-FORM SYSTEM KIT 29X27 CM	C18584V004
PRE5421V005	PREVENA RESTOR BELLA-FORM SYSTEM KIT 29X27 CM	C18585V005
PRE5421V005	PREVENA RESTOR BELLA-FORM SYSTEM KIT 29X27 CM	C19068V005
PRE5501V001	PREVENA RESTOR AXIO-FORM SYSTEM KIT 29X28 CM	C16935V001
PRE5501V001	PREVENA RESTOR AXIO-FORM SYSTEM KIT 29X28 CM	C18076V001
PRE5501V001	PREVENA RESTOR AXIO-FORM SYSTEM KIT 29X28 CM	C18586V001



Field Notice Receipt and Customer Response Form Completion and Certification

Please provide account information.

Facility Name			
Account Number			
Address (no P.O. boxes)			
City, State, Postal Code			
Phone Number		Fax	
E-Mail Address			
Contact Name / Title			
Signature		Date	

Please provide shipping information if different from above.

Facility Name			
Account Number			
Address (no PO boxes)			
City, State, Postal Code			
Phone Number		Fax	
E-Mail Address			
Contact Name / Title			

PLEASE RETURN YOUR COMPLETED FORM TO: kci3mfieldactionresponse@solventum.com
with 3M™ Prevena™ Plus 125 Therapy Units and Kits Field Action in the subject line.