



URGENT –PRODUCT RECALL NOTICE
Option ELITE Vena Cava Filter System
1st Notice

September 12, 2025

To: S Nassau Comm Hosp
2445 Oceanside Rd
Oceanside, NY 11572
USA

Dear Argon Customer,

We are voluntarily recalling five affected lots of the Option ELITE Vena Cava Filter System products in response to reports of increased resistance when advancing the dilator within the introducer sheath.

The immediate health consequence risk is that if the dilator fails to advance through the introducer sheath, it may result in a procedural delay or kinking of the dilator, and may necessitate the opening of a new kit.

Any adverse events related to the use of this product should be reported to the FDA's MedWatch Adverse Event Reporting program either online at www.fda.gov/medwatch/report.htm, by regular mail, or by fax.

Argon records indicate that the following products were shipped to your facility:

Argon Part Number	Lot Number	Shipping Date to your facility	Number of units Shipped to your facility
352506070E	11633166	8/21/2025	2

Note: Appendix A contains the full listing of all part numbers and lot numbers affected.

Please take the following actions:

1. Immediately discontinue use of the affected devices listed above.
2. Segregate and quarantine the affected devices listed above.
3. Complete and return the response form included at the end of this notification via regular mail or email **by 9/26/2025**.
4. Return all affected devices to Argon Medical Devices using the instructions provided below.
5. Share this recall notice with the following parties to ensure your consignees receive a copy and respond accordingly:
 - a. Anyone within your organization who may use this device.
 - b. Any organization where the potentially affected products have been transferred or further distributed.
 - c. Anyone who needs to be aware within your organization.

This voluntary recall should be carried out across your organization and anywhere else these devices may have been transferred or distributed. If the affected product was transferred or distributed to another facility, please provide that facility with the attached blank form so that they can use it to properly account for the affected products, respond appropriately, and return the product.

Argon Medical Devices, Inc.

Please complete and return the response form as soon as possible so that we can account for and track the affected products in your possession. The product may be returned at Argon's expense to the address below with the returned goods authorization (RGA# 29074) clearly marked.

RGA# 29074
Argon Medical Devices, Inc.
1445 Flat Creek Road
Athens, TX 75751 USA
Attn: Arbee Cummings

For additional shipping information and details, please contact Arbee Cummings at 903-676-4220 or arbee.cummings@argonmedical.com.

For additional information on this letter, please email Clay Jones at Clayton.Jones@argonmedical.com or me at Scott.Bishop@argonmedical.com.

Our company is dedicated to delivering high-quality, effective medical devices. We take this seriously and understand that, on rare occasions, these actions are necessary to uphold that commitment. We apologize for any inconvenience this creates for your organization and appreciate your support in accounting for and returning the affected devices in your possession.

Sincerely,



Scott Bishop, MS
Vice President, Regulatory Affairs & Quality Assurance
Argon Medical Devices, Inc.

Cc: Clay Jones , Manager Quality & Compliance

Please proceed to next page to respond to inventory on hand

Option ELITE Vena Cava Filter System
Product Recall Response Form – 1st Notification
RGA# 29074

Customer Name & Address: S Nassau Comm Hosp
2445 Oceanside Rd
Oceanside, NY 11572
USA

Argon Part Number	Shipping Date to your facility	Lot Number	# of units Shipped to your facility	# Currently on hand at your facility	Number to be Returned to Argon
352506070E	8/21/2025	11633166	2		

How did you receive the affected products above? ☐ Directly From Argon ☐ From a Supplier/Distributor

Which Supplier or Distributor (if applicable)? _____

Were there any adverse events associated with recalled product? ☐ Yes ☐ No

If yes, please explain: _____

☐ I read and understand the instructions provided in the recall letter. I checked my stock and quarantined the items listed above.

Signature of Person Completing Form

Printed Name of Person Completing Form

Title of Person Completing Form

Date Signed

Phone Number

Email Address

Proposed Product Return Date:

Please return the complete form to:
Attn: Arbee Cummings
Argon Medical Devices, Inc.
1445 Flat Creek Road,
Athens, TX 75751 USA
Arbee.cummings@argonmedical.com

Appendix A: Complete List of Part and Lots Affected

Argon Medical Devices, Inc

Part Number	Option ELITE Vena Cava Filter System
352506070E	11633166
	11633125
	11629543
	11629273
	11628836

Option ELITE Vena Cava Filter System
Product Recall Response Form – 1st Notification
RGA# 29074

Customer Name & Address:

Argon Part Number	Shipping Date to your facility	Lot Number	# of units Shipped to your facility	# Currently on hand at your facility	Number to be Returned to Argon

How did you receive the affected products above? ☐ Directly From Argon ☐ From a Supplier/Distributor

Which Supplier or Distributor (if applicable)? _____

Were there any adverse events associated with recalled product? ☐ Yes ☐ No

If yes, please explain: _____

☐ I read and understand the instructions provided in the recall letter. I checked my stock and quarantined the items listed above.

Signature of Person Completing Form

Printed Name of Person Completing Form

Title of Person Completing Form

Date Signed

Phone Number

Email Address

Proposed Product Return Date:

Please return the complete form to:

Attn: Arbee Cummings

Argon Medical Devices, Inc.

1445 Flat Creek Road,

Athens, TX 75751 USA

Arbee.cummings@argonmedical.com