

Alert: Siemens Healthineers Issues Correction for 3 Tesla MRI Systems

This recall (</medical-devices/medical-device-recalls/what-medical-device-recall>) involves correcting certain devices and does not involve removing them from where they are used or sold. The FDA has identified this recall as the most serious type. This device may cause serious injury or death if you continue to use it without correction.

Affected Products

Product Name	Material Number	UDI - DI
MAGNETOM Verio	10276755	N/A
MAGNETOM Skyra (DE)	10432915	04056869006703
Biograph mMR	10433372	04056869006741
MAGNETOM Verio Dot	10684333	04056869006673
MAGNETOM Prisma	10849582	04056869006727
MAGNETOM Skyra fit	10849580	04056869006710
MAGNETOM Verio Dot Upgr.	10684334	04056869006680
MAGNETOM Vida (DE)	11060815	04056869039176
MAGNETOM Skyra (CN)	11410327	04056869227139
MAGNETOM Lumina (DE)	11344916	04056869230740
MAGNETOM Vida Fit	11410481	04056869245188
MAGNETOM Vida (CN)	11516152	04056869260143
MAGNETOM Lumina (CN)	11516153	04056869260136
MAGNETOM Skyra Fit BioMatrix	11516217	04056869299273
MAGNETOM Cima.X Upgr. (DE)	11689304	04056869975245
MAGNETOM Cima.X (DE)	11647158	04056869975221
MAGNETOM Connectom.X	11371480	N/A

Feedback

BIOGRAPH One (DE)	11689172	04068151020472
MAGNETOM Cima.X Fit	11647159	04056869975238
MAGNETOM Cima.X (CN)	11689616	04056869996509

What to Do

If the error message “*Magnet Supervision: The average heater power is above the alarm threshold. Please contact Siemens Healthineers Service*” appears on your host computer, no further scans can be performed. Follow additional safety actions immediately.

On August 27, 2025, Siemens Healthineers sent all affected customers an Urgent Medical Device Correction letter recommending the following actions:

- If the error message “*Magnet Supervision: The average heater power is above the alarm threshold. Please contact Siemens Healthineers Service*” appears on the host computer, take the following safety actions:
 - Do not reboot the system.
 - Remove the patient from the system.
 - Lock/secure all access to the scanning room with the system remaining energized at field.
 - Post a “DO NOT USE” notice on or near the device.
 - Call the facility’s local Siemens Healthineers service organization.
 - Make sure all relevant personnel are aware of this issue.
- Be aware that the Magnet Stop is only intended for emergencies. DO NOT initiate a manual quench of the magnet unless there is an emergency.
- Ensure all designated emergency exit routes remain unobstructed at all times and are properly marked. In particular, the exit from the examination room must be maintained at all times.
- Grant access to the system so an inspection can be carried out as directed by Siemens Healthineers service organization.
- Make sure all users of affected products and others who may need to be informed receive the relevant safety information provided within the notice and comply with its recommendations.

- If the equipment is no longer at the facility, forward the notice to the equipment's new owner.
- You may continue to use the system until an inspection has been performed.

Reason for Correction

Siemens Healthineers is correcting 3 Tesla Siemens Healthineers MRI systems due to the potential for an ice blockage to form or currently exist within the magnet venting system. In the event of a sudden loss of superconductivity, or magnet quench, helium gas may be unable to escape through designated venting paths, leading to pressure build-up within the helium containment system. This pressure build up could rupture the helium containment system, potentially resulting in a massive helium leak into the scanning room. Rupture of the helium containment system may also manifest as an explosion, with overpressure in the MRI room and debris expelled forcefully.

The use of affected product may cause serious adverse health consequences, including cold burns, asphyxiation, trauma, and death.

As of August 29, Siemens Healthineers has not reported any serious injuries or deaths related to this issue.

Device Use

Affected 3 Tesla Siemens Healthineers MRI systems include magnetic resonance diagnostic devices (MRDD) and mMR systems, which combine MRDD and Positron Emission Tomography (PET). These devices display the internal structure and/or function of the head, body, or extremities, and may be utilized by appropriately trained health care professionals to aid in the detection, localization, and diagnosis of diseases and disorders.

Contact Information

Customers in the U.S. with questions about this recall should contact the Siemens Healthineers service organization at 1-800-888-7436.

Additional FDA Resources:

- [FDA's Enforcement Report \(https://www.accessdata.fda.gov/scripts/ires/index.cfm?Event=97530\)](https://www.accessdata.fda.gov/scripts/ires/index.cfm?Event=97530).
- [Medical Device Recall Database \(https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfRES/res.cfm?start_search=1&event_id=97530\)](https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfRES/res.cfm?start_search=1&event_id=97530).

Unique Device Identifier (UDI)

The unique device identifier (UDI) helps identify individual medical devices sold in the United States from manufacturing through distribution to patient use. The UDI allows for more accurate reporting, reviewing, and analyzing of adverse event reports so that devices can be identified, and problems potentially corrected more quickly.

- [How do I recognize a UDI on a label? \(/medical-devices/unique-device-identification-system-udi-system/udi-basics\)](/medical-devices/unique-device-identification-system-udi-system/udi-basics)
- [AccessGUDID database - Identify Your Medical Device \(https://accessgudid.nlm.nih.gov/\)](https://accessgudid.nlm.nih.gov/)
- [Benefits of a UDI System \(/medical-devices/unique-device-identification-system-udi-system/benefits-udi-system\)](/medical-devices/unique-device-identification-system-udi-system/benefits-udi-system)

How do I report a problem?

Health care professionals and consumers may [report adverse reactions or quality problems \(https://www.accessdata.fda.gov/scripts/medwatch/index.cfm?action=reporting.home\)](https://www.accessdata.fda.gov/scripts/medwatch/index.cfm?action=reporting.home) they experienced using these devices to MedWatch: The FDA Safety Information and Adverse Event Reporting Program.

Was this page helpful? * (required)

Yes

No

Submit

 An official form of the United States government. Provided by [Touchpoints](#)