

URGENT MEDICAL DEVICE CORRECTION



April 25, 2025

GE HealthCare Ref. # 30116

To: Healthcare Administrator
Network Administrator
Director of IT Department

RE: Potential security vulnerability for MUSE 5 systems

Safety Issue

GE HealthCare has become aware of a potential security vulnerability for MUSE 5 systems. The MUSE Administrator Account active directory or local credential could potentially be used from your local network to log into a GE HealthCare MUSE 5 system and potentially manipulate recorded patient data.

There have been no injuries or unauthorized access to patient data reported to GE HealthCare as a result of this issue.

Actions to be taken by Customer/User

MUSE 5 systems have reached their End of Sales/Service Support as of November 1, 2014. Please see attached for the End of Sales/Service Support letter that was previously distributed. GE HealthCare recommends that customers replace or upgrade these systems.

If you wish to continue using your MUSE 5 system, you should ensure your IT Department provides ongoing prevention of unauthorized access to the system such as not using a well-known password.

If passwords are well-known, please update the credentials according to your organization's password policies or contact GE HealthCare Service or your local Service Representative if necessary.

Please ensure all potential users in your facility are made aware of this safety notification and the recommended actions.

Please retain this document for your records.

Please complete and return the attached Medical Device Notification Acknowledgement Response Form to DCAR.Recall@gehealthcare.com.

Risk Management Department
MAY 06 2025
Received

**Affected
Product Details**

MUSE 5 systems

Intended Use:

Large capacity client-server-based computer system that accesses stores and manages cardiovascular information. The information can consist of measurements, text, digitized waveforms and angiographic images. The AccuSketch Image Capture software in conjunction with the frame grabber board provides a computer-based method for digitizing angiographic images to allow for quantitative analysis. The system accepts video image inputs and performs manipulations on the digitized images to analyze vessel stenoses and ventricular wall motion.

The editable coronary tree is a tool used to electronically annotate and document the anatomy of the patient's vessels.

**Product
Correction**

Please refer to the instructions above under "Actions to be taken by Customer / User".

**Contact
Information**

If you have any questions or concerns regarding this notification, please contact GE HealthCare Service at 1-800-437-1171 or your local Service Representative.

Please be assured that maintaining a high level of safety and quality is our highest priority. If you have any questions, please contact us per the contact information above.

Sincerely,



Laila Gurney
Chief Quality & Regulatory Officer
GE HealthCare



Scott Kelley
Chief Medical & Safety Officer
GE HealthCare



GE HealthCare Ref.# 30116

**URGENT MEDICAL DEVICE CORRECTION ACKNOWLEDGEMENT
RESPONSE REQUIRED**

Please complete this form and return it to GE HealthCare promptly upon receipt and no later than 30 days from receipt. This will confirm receipt and understanding of the Urgent Medical Device Correction notification.

Facility Name: _____
Street Address: _____
City/State/ZIP/Country: _____
Customer Email Address: _____
Customer Phone Number: _____

By signing this form, we acknowledge receipt and understanding of the accompanying Urgent Medical Device Correction notification and that we have informed all potential users and have taken and will take appropriate actions in accordance with that Notification.

Please provide the name of the individual with responsibility who completed this form.

Signature: _____
Printed Name: _____
Position/Job Title: _____
Date (DD/MM/YYYY): _____

Please return completed form by scanning or taking a photo of the completed form and email to: DCAR.Recall@gehealthcare.com





GE Healthcare Lifecare Solutions

Bill Kohn
Global Product Manager

9900 Innovation Dr.
Wauwatosa, WI 53226
USA

T 414 721 2704
F 414 721 3863
William.kohn@ge.com

Recipient name
Recipient title
Recipient address
Recipient address

RE: CSS# or SSA #

31 October 2013
Subject: Product Service Discontinuance Notice

Dear Salutation Name

We know it is important to your long range planning to have advance notice of changes in the status of older medical equipment. GE Healthcare has established a policy of providing such advanced notice to our customers. The purpose of this letter is to provide you with information about our inability to write service contracts or provide service support for some of our older cardiology products.

Our records indicate that you own at least one or more pieces of the equipment identified in Table A. This includes the software, servers and client workstations supplied by GE on which the software was installed. As of the date of this letter, GE will no longer write new service contracts for the product(s) identified in Table A, nor will GE extend existing service contracts for these product(s) beyond the "Service Discontinued" date identified in Table A. Service will no longer support the addition of options or features including, but not limited to, HL7 modules, HL7 configuration changes or enhancements to MUSE v5E systems. If system options, features, or HL7 interface additions are desired the system must be upgraded to the current MUSE version.

As part of your current Software Support Agreement or Comprehensive with Software Support Agreement, this letter serves as 12 months prior written notice from GE of the removal of the products listed in Table A from your service agreement with GE, as of the "Service Discontinued" date. Your service agreement will continue for any products other than the end-of-support product(s) identified in the attached, with an appropriate adjustment of charges. As part of the agreement the MUSE system must be maintained no more than 2 software revisions behind the currently available version. All MUSE v5E systems are 2 revisions behind the current version.

The reasons for this decision are that parts and components are becoming increasingly difficult to procure, and therefore, our service department has fewer resources to support the older equipment, which impairs our ability to provide the kind of efficient, effective service that you expect.

GE Healthcare is committed to helping you provide state-of-the-art patient care. To discuss replacing your older equipment with a product we currently manufacture, please contact your GE Healthcare sales representative at 800-558-5120, Opt 4. If you have any questions regarding this notice, please contact our Technical Support Department at 1-800-558-7044, Option 2.

We appreciate your understanding in this matter.

Sincerely,



Bill Kohn
Global Product Manager

Table A

PRODUCTS DISCONTINUANCE LIST

GE Healthcare Product	Service Discontinued
MUSE Application Software Version 5.E XX	
HP DL360 G3 HIS Server ABK	
HP DL360 G4 HIS Server SCJ	
HP DL360 G5 HIS Server	
	1 November 2014
HP ML 370 G3 MUSE Server – Rack SAE	
HP ML 370 G3 MUSE Server – Tower SAC	
HP ML 370 G4 MUSE Server – Rack SBY	
HP ML 370 G4 MUSE Server – Tower SBX	
HP ML370 G5 MUSE Server – Rack	
HP ML370 G% MUSE Server – Tower	
HP DC 7100 MUSE Client – All	
HP DC 7600 MUSE Client – All	
HP RP 5700	