



July 23, 2025

  
SOUTH NASSAU COMM HOSPITAL  
1 HEALTHY WAY  
OCEANSIDE, NY 11572-1551



Customer Number: 86708

Re: Important Medical Device Information: **Flexicare-BritePro Solo Optic Laryngoscope Handles**

Dear Medical Inventory Controls/Supplies Coordinator/Buyer:

This letter pertains to an **Urgent Voluntary Medical Device Recall**. Please ensure this notice reaches the individual/department responsible for stocking and distributing supplies within your organization.

The product(s) impacted by the field safety notice, as well as the corresponding Item Number, are listed in the chart below. Our records indicate that you may have purchased these items.

| THS Item Number                 | Item Description                                                              | Lot Number               |
|---------------------------------|-------------------------------------------------------------------------------|--------------------------|
| 329-040-02-0110U                | BritePro Solo Single-Use Fiber Optic Mini Handle and Blade Mac 1              | 2007XXXXXX TO 2012XXXXXX |
| 329-040-02-0120U                | BritePro Solo Single-Use Fiber Optic Mini Handle and Blade Mac 2              | 2007XXXXXX TO 2012XXXXXX |
| 329-040-02-0130U                | BritePro Solo Single-Use Fiber Optic Mini Handle and Blade Mac 3              | 2007XXXXXX TO 2012XXXXXX |
| 329-040-02-0131U                | BritePro Solo Single-Use Fiber Optic Mini Handle and Blade Mac 3 Strong Curve | 2007XXXXXX TO 2012XXXXXX |
| 329-040-02-0140U                | BritePro Solo Single-Use Fiber Optic Mini Handle and Blade Mac 4              | 2007XXXXXX TO 2012XXXXXX |
| 329-040-02-0400U                | BritePro Solo Single-Use Fiber Optic Mini Handle and Blade Miller 00          | 2007XXXXXX TO 2012XXXXXX |
| 329-040-02-0410U                | BritePro Solo Single-Use Fiber Optic Mini Handle and Blade Miller 1           | 2007XXXXXX TO 2012XXXXXX |
| 329-040-02-0415U                | BritePro Solo Mini with Pro Miller Blade Size 1.5                             | 2007XXXXXX TO 2012XXXXXX |
| 329-040-02-0420U                | BritePro Solo Single-Use Fiber Optic Mini Handle and Blade Miller 2           | 2007XXXXXX TO 2012XXXXXX |
| 329-040-02-0430U                | BritePro Solo Single-Use Fiber Optic Mini Handle and Blade Miller 3           | 2007XXXXXX TO 2012XXXXXX |
| 329-040-03-0000U                | Britepro Solo Stubby Single-Use Fiber Optic Laryngoscope Handle               | 2007XXXXXX TO 2012XXXXXX |
| 329-040-309U,<br>329-040-309UEA | BritePro Solo Mini Single-Use Fiber Optic Laryngoscope Handle                 | 2007XXXXXX TO 2012XXXXXX |
| 329-040-310AU,<br>329-040-310U  | BritePro Solo Single Use Fiber Optic Laryngoscope Handle                      | 2007XXXXXX TO 2012XXXXXX |



Advancing Healthcare, Together

5000 Tuttle Crossing Blvd. Dublin, OH 43016

800.874.2646

www.tri-anim.com

|               |                                                                  |                          |
|---------------|------------------------------------------------------------------|--------------------------|
| 329-040-331U  | BritePro Solo Single-Use Fiber Optic Handle and Blade Mac 1      | 2007XXXXXX TO 2012XXXXXX |
| 329-040-332U  | BritePro Solo Single-Use Fiber Optic Handle and Blade Mac 2      | 2007XXXXXX TO 2012XXXXXX |
| 329-040-333U  | BritePro Solo Single-Use Fiber Optic Handle and Blade Mac 3      | 2007XXXXXX TO 2012XXXXXX |
| 329-040-334U  | BritePro Solo Single-Use Fiber Optic Handle and Blade Mac 4      | 2007XXXXXX TO 2012XXXXXX |
| 329-040-335U  | BritePro Solo Single-Use Fiber Optic Handle and Blade Mac 3.5    | 2007XXXXXX TO 2012XXXXXX |
| 329-040-341U  | BritePro Solo Single-Use Fiber Optic Handle and Blade Miller 1   | 2007XXXXXX TO 2012XXXXXX |
| 329-040-342U  | BritePro Solo Single-Use Fiber Optic Handle and Blade Miller 2   | 2007XXXXXX TO 2012XXXXXX |
| 329-040-343U  | BritePro Solo Single-Use Fiber Optic Handle and Blade Miller 3   | 2007XXXXXX TO 2012XXXXXX |
| 329-040-344U  | BritePro Solo Single-Use Fiber Optic Handle and Blade Miller 4   | 2007XXXXXX TO 2012XXXXXX |
| 329-040-361U  | BritePro Solo Single-Use Fiber Optic Handle and Blade Phillips 1 | 2007XXXXXX TO 2012XXXXXX |
| 329-040-3612U | BritePro Solo Single-Use Fiber Optic Handle and Blade Phillips 2 | 2007XXXXXX TO 2012XXXXXX |

According to the attached letter from Ambu, this **Urgent Voluntary Medical Device Recall** has been issued due to the BritePro Solo light not illuminating as intended.

Should you determine your inventory contains affected product, please quarantine the product. Pursuant to the attached recall notice, all affected devices must be manually checked to determine if the device fails to illuminate. If the device fails to illuminate contact Tri-anim Health Services Customer Care at 800-533-0523 to identify quarantine product, obtain an RMA for return and receive a credit for affected product.

We apologize for any inconvenience this action may cause. Thank you for allowing Tri-anim Health Services to be your source for emergency medical supplies and equipment.

Sincerely,

Michelle Samples

*Michelle Samples*

Regulatory Affairs  
Enclosure

Tri-anim Health Services, Inc. 5000 Tuttle Crossing Blvd.  
Dublin, OH 43016

614.760.5000  
www.tri-anim.com



# Flexicare Medical Limited

**Date:** July 11, 2025

## URGENT: VOLUNTARY MEDICAL DEVICE RECALL

### BritePro Solo Fiber Optic Laryngoscope Handles

**Dear Distributor / Customer,**

The purpose of this letter is to advise you that Flexicare Medical Limited (Flexicare) is voluntarily recalling specific lot numbers of BritePro Solo Fiber Optic Laryngoscope Handles manufactured in 2020. See table below for affected product codes and affected lots. Our records show that you have devices from the affected LOTS.

These devices are intended for direct vision laryngoscopy to obtain an illuminated view of the vocal cords and the glottis by inserting the laryngoscope blade into the oral cavity to manipulate the tongue or epiglottis.

#### **Attention**

All clinical users in relevant departments including, but not limited to, anesthesia, intensive care, emergency room, obstetrics, emergency crash carts, and pre-hospital emergency care.

#### **Description of the Problem**

Flexicare has identified an increase in global complaints relating to BritePro Solo product lines, reported as issues with the BritePro Solo light not illuminating as intended. We are aware of 203 product failures and 243 complaints associated with the problem. No adverse events (injuries/deaths) have been found to be caused by the device malfunction.

#### **Potential Risk**

Failure of the BritePro Solo light to illuminate could result in intubation being more challenging or not possible using the malfunctioning BritePro Solo. If the light malfunctions, it may prevent effective illumination of the patient airway, potentially causing delay in treatment or requiring intubation using a different device.

Page 1 of 9

---

#### **Flexicare Medical Limited**

Cynon Valley Business Park • Mountain Ash • CF45 4ER • UK  
+44 (0)1443 474647 • [enquiries@flexicare.com](mailto:enquiries@flexicare.com) • [www.flexicare.com](http://www.flexicare.com)

Registered Number: 2428573 • VAT Registration Number: GB 310 726 243





### **Actions to be taken by the Customer/User:**

1. Review Inventory for affected product and immediately quarantine.
2. All affected devices must be manually checked using the following procedure within the "additional information" section.
3. All devices which demonstrate a failure to illuminate shall be withdrawn from use, quarantined, and returned to the address below and as per the customer instructions.
4. If tested devices illuminate, they can continue to be used. Please ensure you subsequently continue to follow the Pre-Use Checks information as per the IFU.
5. Share this notice with any users of the product within your distribution network to ensure they are also aware of this Voluntary Medical Device Recall.
6. If your firm has further distributed the product, please provide them with the notification and inform them to notify their customers.
7. Adverse reactions or quality problems experienced with the use of this product may be reported to FDA's MedWatch Adverse Event Reporting program via:

**Web:** MedWatch website at [www.fda.gov/medwatch](http://www.fda.gov/medwatch)

**Phone:** 1-800-FDA-1088 (1-800-332-1088)

**Mail:** MedWatch, HF-2, FDA, 5600 Fisher's Lane, Rockville, MD 20852-9787

8. In the event of light failure, the caregiver should revert to the local standard of care, including oxygenation, sourcing a functional laryngoscope, or relying on another method to provide patient ventilation.
9. To minimize risk and identify an inoperative light prior to clinical use, pre-checks must be carried out according to standard hospital protocols and the Brite Pro Solo Instructions for Use (IFU).

### **Additional Information**

#### **Pre-Use Checks**

The light source on the BritePro Solo Fiber Optic Laryngoscope Handles should be checked within the packaging prior to use.

The method of detection is visual. The user can ascertain if the handle is illuminating without compromising the sterile packaging by following image 1 of the IFU below.

Spare BritePro Solo Fiber Optic Laryngoscope Handles should also be tested for their functionality and should be pre-tested as per standard hospital protocols and the IFU.

#### **BritePro Solo with blade**



#### **Pre-Use Check of handle Figure 1 of IFU**



**IFU Reference:** IP0000-593B-Airway Management BritePro Solo IFU (v. 052025)

Page 2 of 9



## **Customer Instructions**

Our records indicate that you have been supplied with devices from the affected lots, identified below. Products from the lot numbers should be checked immediately following the Pre-Use Checks (Refer to Additional Information section).

All devices which demonstrate a failure to illuminate shall be withdrawn from use, quarantined, and returned to the address below and as per the customer instructions.

Part numbers of the BritePro Solo Fiber Optic Laryngoscope products included in the Voluntary Medical Device Recall are listed in the table below.

Lot numbers affected start from 2007xxxxx to 2012xxxxx (July-2020 to Dec-2020).

| Part Code    | Product Description                                                           |
|--------------|-------------------------------------------------------------------------------|
| 040-02-0110U | BritePro Solo Single-Use Fiber Optic Mini Handle and Blade Mac 1              |
| 040-02-0120U | BritePro Solo Single-Use Fiber Optic Mini Handle and Blade Mac 2              |
| 040-02-0130U | BritePro Solo Single-Use Fiber Optic Mini Handle and Blade Mac 3              |
| 040-02-0131U | BritePro Solo Single-Use Fiber Optic Mini Handle and Blade Mac 3 Strong Curve |
| 040-02-0140U | BritePro Solo Single-Use Fiber Optic Mini Handle and Blade Mac 4              |
| 040-02-0400U | BritePro Solo Single-Use Fiber Optic Mini Handle and Blade Miller 00          |
| 040-02-0410U | BritePro Solo Single-Use Fiber Optic Mini Handle and Blade Miller 1           |
| 040-02-0415U | BritePro Solo Mini with Pro Miller Blade Size 1.5                             |
| 040-02-0420U | BritePro Solo Single-Use Fiber Optic Mini Handle and Blade Miller 2           |
| 040-02-0430U | BritePro Solo Single-Use Fiber Optic Mini Handle and Blade Miller 3           |
| 040-03-0000U | Britepro Solo Stubby Single-Use Fiber Optic Laryngoscope Handle               |
| 040-309U     | BritePro Solo Mini Single-Use Fiber Optic Laryngoscope Handle                 |
| 040-310AU    | BritePro Solo Single Use Fiber Optic Laryngoscope Handle                      |
| 040-310U     | BritePro Solo Single-Use Fiber Optic Laryngoscope Handle                      |
| 040-331U     | BritePro Solo Single-Use Fiber Optic Handle and Blade Mac 1                   |
| 040-332U     | BritePro Solo Single-Use Fiber Optic Handle and Blade Mac 2                   |
| 040-333U     | BritePro Solo Single-Use Fiber Optic Handle and Blade Mac 3                   |
| 040-334U     | BritePro Solo Single-Use Fiber Optic Handle and Blade Mac 4                   |
| 040-335U     | BritePro Solo Single-Use Fiber Optic Handle and Blade Mac 3.5                 |
| 040-341U     | BritePro Solo Single-Use Fiber Optic Handle and Blade Miller 1                |
| 040-342U     | BritePro Solo Single-Use Fiber Optic Handle and Blade Miller 2                |
| 040-343U     | BritePro Solo Single-Use Fiber Optic Handle and Blade Miller 3                |
| 040-344U     | BritePro Solo Single-Use Fiber Optic Handle and Blade Miller 4                |
| 040-361U     | BritePro Solo Single-Use Fiber Optic Handle and Blade Phillips 1              |
| 040-362U     | BritePro Solo Single-Use Fiber Optic Handle and Blade Phillips 2              |

### **Flexicare Medical Limited**

Cynon Valley Business Park • Mountain Ash • CF45 4ER • UK  
+44 (0)1443 474647 • enquiries@flexicare.com • www.flexicare.com

Registered Number: 2428573 • VAT Registration Number: GB 310 726 243





### Acknowledgement Instructions

Please complete the Acknowledgement Form by selecting the link or completing the form included with this communication.

Please complete the electronic link [Click on Link](#) or attached Acknowledgement Form and return to Flexicare by completing **the options** listed below:

**Flexicare Inc.**

Address: 15281 Barranca Pkwy Unit D, Irvine, CA 92618

Name: Regulatory and Quality Team

Phone: (949) 450-9999

Email: [safetynotice@flexicare.com](mailto:safetynotice@flexicare.com)

Fax: (949) 450-9992

If you have any questions, please call (949) 450-9999 and reference BritePro Solo Fiber Optic Laryngoscope Safety Notice and you will be directed to the appropriate person to assist you.

Once you have submitted the Acknowledgement Form, you will be contacted by a Flexicare representative to make necessary arrangements to comply with the Voluntary Medical Device Recall.

Flexicare has identified the cause of the light issues, and appropriate Corrective and Preventative Actions have been completed and implemented.

We apologize for any inconvenience this matter may cause you.

If you have any questions, please contact Flexicare:

**Contact Information:** (Office Hours: Monday through Friday, 8:00 AM – 5:00 PM PDT)

**Flexicare Inc.**

Address: 15281 Barranca Pkwy Unit D, Irvine, CA 92618

Name: Regulatory and Quality Team

Phone: (949) 450-9999

Email: [safetynotice@flexicare.com](mailto:safetynotice@flexicare.com)

Fax: (949) 450-9992

Authorized by:

Name: REBECCA FUNSTON

Signature:

  
Rebecca Concordia Funston | Jul 11, 2025 14:19 GMT+1

Title:

RA/QA DIRECTOR

Page 4 of 9

**Flexicare Medical Limited**

Cynon Valley Business Park • Mountain Ash • CF45 4ER • UK  
+44 (0)1443 474647 • enquiries@flexicare.com • www.flexicare.com

Registered Number: 2425573 • VAT Registration Number: GB 310 720 245





## Privacy Notice

Flexicare is the controller responsible for personal data submitted under this Voluntary Medical Device Recall. We will collect and process personal data when you provide the information requested on this form.

We will only collect the personal data requested on this Voluntary Medical Device Recall or that you otherwise provide to us as part of this Voluntary Medical Device Recall. This data may include identity data, contact data, transaction data, technical data, and profile data. We will process this data for the purpose of complying with our legal obligations as a manufacturer of regulated medical devices, contacting you regarding the Voluntary Medical Device Recall, and processing information regarding the devices identified in any responses submitted (the legal basis). By submitting the return form, you consent to the processing of personal data for the legal basis.

We may share personal data collected under this Voluntary Medical Device Recall within the Flexicare Group or with external third parties such as government regulators or service providers assisting with the Voluntary Medical Device Recall process. We may transfer, store and process your personal data in the United Kingdom where our headquarters are located (at the address below). Flexicare adheres to the UK GDPR, Data Protection Act 2018, and protects personal data with approved contractual clauses.

You have the right to access your personal information, request deletion of your personal information, and opt-out of the sale of your personal information (although Flexicare does not sell your personal information).

[Flexicare's Privacy Policy](#) contains more detailed information about our data processing.

If you have any questions about this privacy notice or about our data protection practices, or to exercise your rights, please contact our Data Protection Officer at the contact details below:

Mr. Rhodri Evans  
legal@flexicare.com  
+44 (0)1443 474647

---

### Flexicare Medical Limited

Cynon Valley Business Park • Mountain Ash • CF45 4ER • UK  
+44 (0)1443 474647 • enquiries@flexicare.com • www.flexicare.com

Registered Number: 2428573 • VAT Registration Number: GB 310 726 243

