



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

September 18, 2025

URGENT: MEDICAL DEVICE RECALL ZENITH ALPHA® 2 THORACIC ENDOVASCULAR GRAFT

ATTENTION:

Risk Management/Recall Administration

Our records indicate that you have received affected products.

Patient safety and customer satisfaction are Cook Medical's highest priorities. We appreciate your time and attention in reading this important notice.

Purpose of this Letter

The purpose of this letter is to inform you that Cook Medical is voluntarily recalling specific Reference Part Numbers (RPNs) and lot numbers pertaining to the ZENITH ALPHA® 2 THORACIC ENDOVASCULAR GRAFT, Proximal Component (tapered or non-tapered configurations), with diameters of 40mm, 42mm, 44mm, and 46mm, all lengths included.

This voluntary recall does **not** pertain to:

- ZENITH ALPHA® 2 THORACIC ENDOVASCULAR GRAFT Proximal Components with diameters of 24-38mm.
- ZENITH ALPHA® 2 THORACIC ENDOVASCULAR GRAFT Distal Components of any size
- ZENITH ALPHA® 2 THORACIC ENDOVASCULAR GRAFT Distal Extension Components of any size
- ZENITH ALPHA® THORACIC ENDOVASCULAR GRAFT (Previous Generation devices) of any size

Reason for Voluntary Recall

Cook Medical identified during manufacturing that specific ZENITH ALPHA® 2 THORACIC ENDOVASCULAR GRAFT devices may contain scrapings of polytetrafluoroethylene (PTFE) coating from the PTFE coated nitinol release-wire inside the stent graft which may be released during deployment. See Table 1 - "Affected Device Information" for list of affected devices.

Risk to Health

Potential adverse events that may occur if an affected device is used concern harms associated with:

- Device fragment (non-retrievable) resulting in no apparent sequelae and unlikely to cause future harm
- Device fragment embolization resulting in permanent impairment, life-threatening harm, and/or death (e.g., stroke).

To date, Cook Medical has received no customer complaints related to the adverse patient effects listed above for the affected devices.



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

Product Information

The ZENITH ALPHA® 2 THORACIC ENDOVASCULAR GRAFT consists of a stent graft preloaded into an introduction system.

The ZENITH ALPHA® 2 THORACIC ENDOVASCULAR GRAFT is a two-piece cylindrical endovascular system consisting of proximal (ZTA2-P(T)-/-) and distal (ZTA2-D-/-) stent graft components. The proximal component can be either tapered (ZTA2-PT-/-) or non-tapered (ZTA2-P-/-). An endovascular ancillary component is available as a distal extension (ZTA2-DE-/-). The stent grafts are constructed of woven polyester fabric sewn to self-expanding nitinol (a nickel-titanium based alloy) stents with braided polyester and monofilament polypropylene suture.

The proximal component has one proximal bare stent, one proximal internal sealing stent with fixation bars that protrude through the graft material, and one distal internal sealing stent. The distal component has three proximal internal sealing stents, one distal internal sealing stent, and one distal bare stent with bars. The distal extension has three proximal internal sealing stents and one distal internal sealing stent. Gold markers are placed on the graft at the proximal and distal aspects of the graft margins on all components.

Actions to be Taken by the Customer

1. Examine inventory immediately to determine if you have affected product(s) and quarantine any affected product that remains unused. Immediately cease all further distribution and use of the affected products.
2. Return the affected product(s) to Cook Medical with a copy of the Acknowledgement and Receipt Form to receive a product credit. Refer to the Acknowledgement and Receipt Form for return instructions.
3. Please complete the Acknowledgement and Receipt Form within **5 business** days of receiving this letter regardless of whether you have affected product(s) on hand. You may return the Form via fax (812.339.7316) or email (FieldActionsNA@CookMedical.com).
4. This notice must be shared with appropriate personnel, including down to the user level, within your organization or with any organization where the affected devices have been transferred.
5. Please immediately report adverse events to Cook Medical Customer Relations by phone at 800.457.4500 or 812.339.2235, Monday through Friday between 7:30am and 5:00pm (Eastern Time), or by email to: CustomerRelationsNA@CookMedical.com.

Action being Taken by Cook Medical

Cook Medical is removing potentially affected products from the market. Additionally, an investigation has been opened to determine the appropriate corrections to prevent the recurrence of a similar issue.



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

Other Information

This action is being taken with the knowledge of the Food and Drug Administration. Adverse events or quality problems experienced with the use of this product may also be reported to the FDA:

- Visit <http://www.fda.gov/Safety/MedWatch/HowToReport/default.htm> to obtain a form to fax or mail.
- Call the FDA at 800.FDA.1088.

We recognize that this situation is a potential disruption to your normal operations, and we sincerely apologize. Thank you for your immediate assistance in this matter. If you have any questions or concerns, please contact Cook Medical Customer Relations at 800.457.4500 or 812.339.2235. We look forward to your return of the Acknowledgement and Receipt Form.

A handwritten signature in black ink that reads "Aoife Doyle". The signature is fluid and cursive, with the first name "Aoife" being more prominent than the last name "Doyle".

Aoife Doyle
Director, Quality Assurance
Cook Medical



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

Table 1 – Affected Device Information

REFERENCE PART NUMBER (RPN)	ORDER NUMBER (GPN)	LOT NUMBER	BOX UNIQUE DEVICE IDENTIFIER (UDI)	MANUFACTURE DATE	EXPIRATION DATE	QUANTITY
ZTA2-P-40-117-W	G59915	E4626917	(01)10827002599151(17)271114(30)1(10)E4626917	14-nov-24	14-nov-27	1
ZTA2-P-40-117-W	G59915	E4627437	(01)10827002599151(17)271125(30)1(10)E4627437	25-nov-24	25-nov-27	1
ZTA2-P-40-117-W	G59915	E4627438	(01)10827002599151(17)271118(30)1(10)E4627438	18-nov-24	18-nov-27	1
ZTA2-P-40-117-W	G59915	E4627439	(01)10827002599151(17)271118(30)1(10)E4627439	18-nov-24	18-nov-27	1
ZTA2-P-40-117-W	G59915	E4627440	(01)10827002599151(17)271118(30)1(10)E4627440	18-nov-24	18-nov-27	1
ZTA2-P-40-167-W	G59921	E4627442	(01)10827002599212(17)271115(30)1(10)E4627442	15-nov-24	15-nov-27	1
ZTA2-P-40-167-W	G59921	E4627444	(01)10827002599212(17)271206(30)1(10)E4627444	06-dec-24	06-dec-27	1
ZTA2-P-40-167-W	G59921	E4627445	(01)10827002599212(17)271115(30)1(10)E4627445	15-nov-24	15-nov-27	1
ZTA2-P-40-167-W	G59921	E4627446	(01)10827002599212(17)271115(30)1(10)E4627446	15-nov-24	15-nov-27	1
ZTA2-P-40-167-W	G59921	E4627447	(01)10827002599212(17)271115(30)1(10)E4627447	15-nov-24	15-nov-27	1
ZTA2-P-40-167-W	G59921	E4627448	(01)10827002599212(17)271115(30)1(10)E4627448	15-nov-24	15-nov-27	1
ZTA2-P-40-167-W	G59921	E4627449	(01)10827002599212(17)271118(30)1(10)E4627449	18-nov-24	18-nov-27	1
ZTA2-P-40-167-W	G59921	E4627450	(01)10827002599212(17)271115(30)1(10)E4627450	15-nov-24	15-nov-27	1
ZTA2-P-40-167-W	G59921	E4627451	(01)10827002599212(17)271118(30)1(10)E4627451	18-nov-24	18-nov-27	1
ZTA2-P-40-167-W	G59921	E4627452	(01)10827002599212(17)271115(30)1(10)E4627452	15-nov-24	15-nov-27	1
ZTA2-PT-42-38-173-W	G60070	E4627453	(01)10827002600703(17)271118(30)1(10)E4627453	18-nov-24	18-nov-27	1
ZTA2-PT-42-38-173-W	G60070	E4627454	(01)10827002600703(17)271128(30)1(10)E4627454	28-nov-24	28-nov-27	1
ZTA2-PT-42-38-173-W	G60070	E4627455	(01)10827002600703(17)271118(30)1(10)E4627455	18-nov-24	18-nov-27	1
ZTA2-PT-42-38-173-W	G60070	E4627456	(01)10827002600703(17)271118(30)1(10)E4627456	18-nov-24	18-nov-27	1
ZTA2-PT-42-38-173-W	G60070	E4627457	(01)10827002600703(17)271118(30)1(10)E4627457	18-nov-24	18-nov-27	1
ZTA2-PT-42-38-173-W	G60070	E4627458	(01)10827002600703(17)271125(30)1(10)E4627458	25-nov-24	25-nov-27	1
ZTA2-PT-42-38-173-W	G60070	E4627459	(01)10827002600703(17)271118(30)1(10)E4627459	18-nov-24	18-nov-27	1
ZTA2-PT-42-38-173-W	G60070	E4627460	(01)10827002600703(17)271118(30)1(10)E4627460	18-nov-24	18-nov-27	1
ZTA2-P-42-173-W	G59932	E4628827	(01)10827002599328(17)271120(30)1(10)E4628827	20-nov-24	20-nov-27	1
ZTA2-P-42-173-W	G59932	E4628828	(01)10827002599328(17)271120(30)1(10)E4628828	20-nov-24	20-nov-27	1
ZTA2-P-42-173-W	G59932	E4628829	(01)10827002599328(17)271120(30)1(10)E4628829	20-nov-24	20-nov-27	1
ZTA2-P-42-173-W	G59932	E4628830	(01)10827002599328(17)271120(30)1(10)E4628830	20-nov-24	20-nov-27	1
ZTA2-P-42-173-W	G59932	E4628831	(01)10827002599328(17)271120(30)1(10)E4628831	20-nov-24	20-nov-27	1
ZTA2-P-42-173-W	G59932	E4628832	(01)10827002599328(17)271120(30)1(10)E4628832	20-nov-24	20-nov-27	1
ZTA2-P-42-173-W	G59932	E4628833	(01)10827002599328(17)271120(30)1(10)E4628833	20-nov-24	20-nov-27	1
ZTA2-P-42-173-W	G59932	E4628834	(01)10827002599328(17)271120(30)1(10)E4628834	20-nov-24	20-nov-27	1
ZTA2-P-42-173-W	G59932	E4628835	(01)10827002599328(17)271120(30)1(10)E4628835	20-nov-24	20-nov-27	1
ZTA2-P-42-173-W	G59932	E4628836	(01)10827002599328(17)271121(30)1(10)E4628836	21-nov-24	21-nov-27	1
ZTA2-P-46-125-W	G59950	E4628838	(01)10827002599502(17)271118(30)1(10)E4628838	18-nov-24	18-nov-27	1
ZTA2-P-46-125-W	G59950	E4628839	(01)10827002599502(17)271120(30)1(10)E4628839	20-nov-24	20-nov-27	1
ZTA2-P-46-125-W	G59950	E4628841	(01)10827002599502(17)271118(30)1(10)E4628841	18-nov-24	18-nov-27	1
ZTA2-P-46-125-W	G59950	E4628842	(01)10827002599502(17)271119(30)1(10)E4628842	19-nov-24	19-nov-27	1



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

REFERENCE PART NUMBER (RPN)	ORDER NUMBER (OPN)	LOT NUMBER	BOX UNIQUE DEVICE IDENTIFIER (UDI)	MANUFACTURE DATE	EXPIRATION DATE	QUANTITY
ZTA2-P-46-125-W	G59950	E4628843	(01)10827002599502(17)271119(30)1(10)E4628843	19-nov-24	19-nov-27	1
ZTA2-P-46-125-W	G59950	E4628844	(01)10827002599502(17)271120(30)1(10)E4628844	20-nov-24	20-nov-27	1
ZTA2-P-46-125-W	G59950	E4628845	(01)10827002599502(17)271119(30)1(10)E4628845	19-nov-24	19-nov-27	1
ZTA2-P-46-125-W	G59950	E4628846	(01)10827002599502(17)271119(30)1(10)E4628846	19-nov-24	19-nov-27	1
ZTA2-P-46-125-W	G59950	E4628848	(01)10827002599502(17)271118(30)1(10)E4628848	18-nov-24	18-nov-27	1
ZTA2-P-46-179-W	G59910	E4628849	(01)10827002599106(17)271203(30)1(10)E4628849	03-dec-24	03-dec-27	1
ZTA2-P-46-179-W	G59910	E4628850	(01)10827002599106(17)271205(30)1(10)E4628850	05-dec-24	05-dec-27	1
ZTA2-P-46-179-W	G59910	E4628851	(01)10827002599106(17)271121(30)1(10)E4628851	21-nov-24	21-nov-27	1
ZTA2-P-46-179-W	G59910	E4628852	(01)10827002599106(17)271126(30)1(10)E4628852	26-nov-24	26-nov-27	1
ZTA2-P-46-179-W	G59910	E4628853	(01)10827002599106(17)271128(30)1(10)E4628853	28-nov-24	28-nov-27	1
ZTA2-P-46-179-W	G59910	E4628854	(01)10827002599106(17)271120(30)1(10)E4628854	20-nov-24	20-nov-27	1
ZTA2-P-46-179-W	G59910	E4628855	(01)10827002599106(17)271121(30)1(10)E4628855	21-nov-24	21-nov-27	1
ZTA2-P-46-179-W	G59910	E4628856	(01)10827002599106(17)271128(30)1(10)E4628856	28-nov-24	28-nov-27	1
ZTA2-P-46-179-W	G59910	E4628857	(01)10827002599106(17)271120(30)1(10)E4628857	20-nov-24	20-nov-27	1
ZTA2-P-46-179-W	G59910	E4628858	(01)10827002599106(17)271204(30)1(10)E4628858	04-dec-24	04-dec-27	1
ZTA2-P-46-179-W	G59910	E4628859	(01)10827002599106(17)271126(30)1(10)E4628859	26-nov-24	26-nov-27	1
ZTA2-PT-40-36-167-W	G60064	E4628860	(01)10827002600642(17)271127(30)1(10)E4628860	27-nov-24	27-nov-27	1
ZTA2-PT-40-36-167-W	G60064	E4628861	(01)10827002600642(17)271127(30)1(10)E4628861	27-nov-24	27-nov-27	1
ZTA2-PT-40-36-167-W	G60064	E4628862	(01)10827002600642(17)271127(30)1(10)E4628862	27-nov-24	27-nov-27	1
ZTA2-PT-40-36-167-W	G60064	E4628863	(01)10827002600642(17)271127(30)1(10)E4628863	27-nov-24	27-nov-27	1
ZTA2-PT-40-36-167-W	G60064	E4628864	(01)10827002600642(17)271125(30)1(10)E4628864	25-nov-24	25-nov-27	1
ZTA2-PT-40-36-167-W	G60064	E4628865	(01)10827002600642(17)271125(30)1(10)E4628865	25-nov-24	25-nov-27	1
ZTA2-PT-40-36-167-W	G60064	E4628866	(01)10827002600642(17)271121(30)1(10)E4628866	21-nov-24	21-nov-27	1
ZTA2-PT-40-36-167-W	G60064	E4628867	(01)10827002600642(17)271205(30)1(10)E4628867	05-dec-24	05-dec-27	1
ZTA2-PT-40-36-167-W	G60064	E4628868	(01)10827002600642(17)271120(30)1(10)E4628868	20-nov-24	20-nov-27	1
ZTA2-PT-40-36-167-W	G60064	E4628869	(01)10827002600642(17)271120(30)1(10)E4628869	20-nov-24	20-nov-27	1
ZTA2-PT-40-36-167-W	G60064	E4628870	(01)10827002600642(17)271121(30)1(10)E4628870	21-nov-24	21-nov-27	1
ZTA2-PT-40-36-167-W	G60064	E4628871	(01)10827002600642(17)271125(30)1(10)E4628871	25-nov-24	25-nov-27	1
ZTA2-PT-40-36-167-W	G60064	E4628872	(01)10827002600642(17)271125(30)1(10)E4628872	25-nov-24	25-nov-27	1
ZTA2-PT-40-36-217-W	G60067	E4628873	(01)10827002600673(17)271119(30)1(10)E4628873	19-nov-24	19-nov-27	1
ZTA2-PT-40-36-217-W	G60067	E4628874	(01)10827002600673(17)271126(30)1(10)E4628874	26-nov-24	26-nov-27	1
ZTA2-PT-40-36-217-W	G60067	E4628875	(01)10827002600673(17)271119(30)1(10)E4628875	19-nov-24	19-nov-27	1
ZTA2-PT-40-36-217-W	G60067	E4628876	(01)10827002600673(17)271119(30)1(10)E4628876	19-nov-24	19-nov-27	1
ZTA2-PT-40-36-217-W	G60067	E4628877	(01)10827002600673(17)271119(30)1(10)E4628877	19-nov-24	19-nov-27	1
ZTA2-P-42-121-W	G59926	E4629315	(01)10827002599267(17)271126(30)1(10)E4629315	26-nov-24	26-nov-27	1
ZTA2-P-40-217-W	G59923	E4629331	(01)10827002599236(17)271126(30)1(10)E4629331	26-nov-24	26-nov-27	1
ZTA2-P-40-217-W	G59923	E4629332	(01)10827002599236(17)271126(30)1(10)E4629332	26-nov-24	26-nov-27	1
ZTA2-P-40-217-W	G59923	E4629335	(01)10827002599236(17)271126(30)1(10)E4629335	26-nov-24	26-nov-27	1



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

REFERENCE PART NUMBER (RPN)	ORDER NUMBER (GPN)	LOT NUMBER	BOX UNIQUE DEVICE IDENTIFIER (UDI)	MANUFACTURE DATE	EXPIRATION DATE	QUANTITY
ZTA2-P-40-217-W	G59923	E4629337	(01)10827002599236(17)271125(30)1(10)E4629337	25-nov-24	25-nov-27	1
ZTA2-P-40-217-W	G59923	E4629338	(01)10827002599236(17)271126(30)1(10)E4629338	26-nov-24	26-nov-27	1
ZTA2-P-40-217-W	G59923	E4629339	(01)10827002599236(17)271126(30)1(10)E4629339	26-nov-24	26-nov-27	1
ZTA2-P-42-121-W	G59926	E4629341	(01)10827002599267(17)271125(30)1(10)E4629341	25-nov-24	25-nov-27	1
ZTA2-P-42-121-W	G59926	E4629342	(01)10827002599267(17)271126(30)1(10)E4629342	26-nov-24	26-nov-27	1
ZTA2-P-42-121-W	G59926	E4629343	(01)10827002599267(17)271125(30)1(10)E4629343	25-nov-24	25-nov-27	1
ZTA2-P-42-121-W	G59926	E4629344	(01)10827002599267(17)271126(30)1(10)E4629344	26-nov-24	26-nov-27	1
ZTA2-P-42-121-W	G59926	E4629345	(01)10827002599267(17)271125(30)1(10)E4629345	25-nov-24	25-nov-27	1
ZTA2-P-42-121-W	G59926	E4629346	(01)10827002599267(17)271125(30)1(10)E4629346	25-nov-24	25-nov-27	1
ZTA2-P-42-121-W	G59926	E4629347	(01)10827002599267(17)271125(30)1(10)E4629347	25-nov-24	25-nov-27	1
ZTA2-P-42-121-W	G59926	E4629348	(01)10827002599267(17)271126(30)1(10)E4629348	26-nov-24	26-nov-27	1
ZTA2-P-42-121-W	G59926	E4629349	(01)10827002599267(17)271126(30)1(10)E4629349	26-nov-24	26-nov-27	1
ZTA2-PT-46-42-179-W	G60079	E4629937	(01)10827002600796(17)271128(30)1(10)E4629937	28-nov-24	28-nov-27	1
ZTA2-PT-46-42-179-W	G60079	E4629938	(01)10827002600796(17)271128(30)1(10)E4629938	28-nov-24	28-nov-27	1
ZTA2-PT-46-42-179-W	G60079	E4629939	(01)10827002600796(17)271204(30)1(10)E4629939	04-dec-24	04-dec-27	1
ZTA2-PT-46-42-179-W	G60079	E4629940	(01)10827002600796(17)271128(30)1(10)E4629940	28-nov-24	28-nov-27	1
ZTA2-PT-46-42-179-W	G60079	E4629941	(01)10827002600796(17)271128(30)1(10)E4629941	28-nov-24	28-nov-27	1
ZTA2-PT-46-42-179-W	G60079	E4629942	(01)10827002600796(17)271128(30)1(10)E4629942	28-nov-24	28-nov-27	1
ZTA2-PT-46-42-233-W	G60082	E4629943	(01)10827002600826(17)271126(30)1(10)E4629943	26-nov-24	26-nov-27	1
ZTA2-PT-46-42-233-W	G60082	E4629945	(01)10827002600826(17)271127(30)1(10)E4629945	27-nov-24	27-nov-27	1
ZTA2-PT-46-42-233-W	G60082	E4629947	(01)10827002600826(17)271127(30)1(10)E4629947	27-nov-24	27-nov-27	1
ZTA2-P-40-117-W	G59915	E4629948	(01)10827002599151(17)271127(30)1(10)E4629948	27-nov-24	27-nov-27	1
ZTA2-P-40-117-W	G59915	E4629949	(01)10827002599151(17)271127(30)1(10)E4629949	27-nov-24	27-nov-27	1
ZTA2-P-40-117-W	G59915	E4629950	(01)10827002599151(17)271127(30)1(10)E4629950	27-nov-24	27-nov-27	1
ZTA2-P-40-117-W	G59915	E4629951	(01)10827002599151(17)271127(30)1(10)E4629951	27-nov-24	27-nov-27	1
ZTA2-P-40-117-W	G59915	E4629952	(01)10827002599151(17)271127(30)1(10)E4629952	27-nov-24	27-nov-27	1
ZTA2-P-40-117-W	G59915	E4629953	(01)10827002599151(17)271127(30)1(10)E4629953	27-nov-24	27-nov-27	1
ZTA2-P-40-117-W	G59915	E4629954	(01)10827002599151(17)271126(30)1(10)E4629954	26-nov-24	26-nov-27	1
ZTA2-P-40-117-W	G59915	E4629955	(01)10827002599151(17)271127(30)1(10)E4629955	27-nov-24	27-nov-27	1
ZTA2-P-42-225-W	G59935	E4629956	(01)10827002599359(17)271127(30)1(10)E4629956	27-nov-24	27-nov-27	1
ZTA2-P-42-225-W	G59935	E4629957	(01)10827002599359(17)271127(30)1(10)E4629957	27-nov-24	27-nov-27	1
ZTA2-P-42-225-W	G59935	E4629958	(01)10827002599359(17)271128(30)1(10)E4629958	28-nov-24	28-nov-27	1
ZTA2-P-42-225-W	G59935	E4629959	(01)10827002599359(17)271127(30)1(10)E4629959	27-nov-24	27-nov-27	1
ZTA2-P-44-125-W	G59938	E4629960	(01)10827002599380(17)271127(30)1(10)E4629960	27-nov-24	27-nov-27	1
ZTA2-P-44-125-W	G59938	E4629961	(01)10827002599380(17)271128(30)1(10)E4629961	28-nov-24	28-nov-27	1
ZTA2-P-44-125-W	G59938	E4629962	(01)10827002599380(17)271205(30)1(10)E4629962	05-dec-24	05-dec-27	1
ZTA2-P-44-125-W	G59938	E4629963	(01)10827002599380(17)271127(30)1(10)E4629963	27-nov-24	27-nov-27	1
ZTA2-P-44-125-W	G59938	E4629964	(01)10827002599380(17)271127(30)1(10)E4629964	27-nov-24	27-nov-27	1



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

REFERENCE PART NUMBER (RPN)	ORDER NUMBER (GPN)	LOT NUMBER	BOX UNIQUE DEVICE IDENTIFIER (UDI)	MANUFACTURE DATE	EXPIRATION DATE	QUANTITY
ZTA2-P-44-125-W	G59938	E4629965	(01)10827002599380(17)271128(30)1(10)E4629965	28-nov-24	28-nov-27	1
ZTA2-P-46-233-W	G59912	E4629966	(01)10827002599120(17)271128(30)1(10)E4629966	28-nov-24	28-nov-27	1
ZTA2-P-46-233-W	G59912	E4629967	(01)10827002599120(17)271128(30)1(10)E4629967	28-nov-24	28-nov-27	1
ZTA2-P-46-233-W	G59912	E4629968	(01)10827002599120(17)271128(30)1(10)E4629968	28-nov-24	28-nov-27	1
ZTA2-PT-42-38-225-W	G60073	E4629969	(01)10827002600734(17)271127(30)1(10)E4629969	27-nov-24	27-nov-27	1
ZTA2-PT-42-38-225-W	G60073	E4629970	(01)10827002600734(17)271205(30)1(10)E4629970	05-dec-24	05-dec-27	1
ZTA2-PT-42-38-225-W	G60073	E4629971	(01)10827002600734(17)271127(30)1(10)E4629971	27-nov-24	27-nov-27	1
ZTA2-PT-42-38-225-W	G60073	E4629972	(01)10827002600734(17)271128(30)1(10)E4629972	28-nov-24	28-nov-27	1
ZTA2-PT-42-38-225-W	G60073	E4629973	(01)10827002600734(17)271128(30)1(10)E4629973	28-nov-24	28-nov-27	1
ZTA2-PT-46-42-179-W	G80079	E4629974	(01)10827002600796(17)271128(30)1(10)E4629974	28-nov-24	28-nov-27	1
ZTA2-P-40-167-W	G59921	E4631968	(01)10827002599212(17)271129(30)1(10)E4631968	29-nov-24	29-nov-27	1
ZTA2-P-40-167-W	G59921	E4631969	(01)10827002599212(17)271129(30)1(10)E4631969	29-nov-24	29-nov-27	1
ZTA2-P-40-167-W	G59921	E4631970	(01)10827002599212(17)271127(30)1(10)E4631970	27-nov-24	27-nov-27	1
ZTA2-P-44-179-W	G59944	E4631971	(01)10827002599441(17)271127(30)1(10)E4631971	27-nov-24	27-nov-27	1
ZTA2-P-44-179-W	G59944	E4631972	(01)10827002599441(17)271127(30)1(10)E4631972	27-nov-24	27-nov-27	1
ZTA2-P-44-179-W	G59944	E4631973	(01)10827002599441(17)271127(30)1(10)E4631973	27-nov-24	27-nov-27	1
ZTA2-P-44-179-W	G59944	E4631974	(01)10827002599441(17)271128(30)1(10)E4631974	28-nov-24	28-nov-27	1
ZTA2-P-44-179-W	G59944	E4631975	(01)10827002599441(17)271128(30)1(10)E4631975	28-nov-24	28-nov-27	1
ZTA2-P-44-179-W	G59944	E4631976	(01)10827002599441(17)271128(30)1(10)E4631976	28-nov-24	28-nov-27	1
ZTA2-P-44-233-W	G59947	E4631977	(01)10827002599472(17)271129(30)1(10)E4631977	29-nov-24	29-nov-27	1
ZTA2-P-44-233-W	G59947	E4631978	(01)10827002599472(17)271201(30)1(10)E4631978	01-dec-24	01-dec-27	1
ZTA2-P-44-233-W	G59947	E4631979	(01)10827002599472(17)271129(30)1(10)E4631979	29-nov-24	29-nov-27	1
ZTA2-P-44-233-W	G59947	E4631980	(01)10827002599472(17)271129(30)1(10)E4631980	29-nov-24	29-nov-27	1
ZTA2-PT-40-36-167-W	G60064	E4631981	(01)10827002600642(17)271129(30)1(10)E4631981	29-nov-24	29-nov-27	1
ZTA2-PT-40-36-167-W	G60064	E4631982	(01)10827002600642(17)271129(30)1(10)E4631982	29-nov-24	29-nov-27	1
ZTA2-PT-40-36-167-W	G60064	E4631983	(01)10827002600642(17)271202(30)1(10)E4631983	02-dec-24	02-dec-27	1
ZTA2-PT-40-36-167-W	G60064	E4631984	(01)10827002600642(17)271129(30)1(10)E4631984	29-nov-24	29-nov-27	1
ZTA2-PT-40-36-167-W	G60064	E4631985	(01)10827002600642(17)271203(30)1(10)E4631985	03-dec-24	03-dec-27	1
ZTA2-PT-40-36-167-W	G60064	E4631986	(01)10827002600642(17)271129(30)1(10)E4631986	29-nov-24	29-nov-27	1
ZTA2-PT-40-36-217-W	G60067	E4631987	(01)10827002600673(17)271129(30)1(10)E4631987	29-nov-24	29-nov-27	1
ZTA2-PT-40-36-217-W	G60067	E4631988	(01)10827002600673(17)271129(30)1(10)E4631988	29-nov-24	29-nov-27	1
ZTA2-PT-40-36-217-W	G60067	E4631989	(01)10827002600673(17)271129(30)1(10)E4631989	29-nov-24	29-nov-27	1
ZTA2-PT-40-36-217-W	G60067	E4631990	(01)10827002600673(17)271129(30)1(10)E4631990	29-nov-24	29-nov-27	1
ZTA2-PT-40-36-217-W	G60067	E4631991	(01)10827002600673(17)271129(30)1(10)E4631991	29-nov-24	29-nov-27	1
ZTA2-PT-40-36-217-W	G60067	E4631992	(01)10827002600673(17)271129(30)1(10)E4631992	29-nov-24	29-nov-27	1
ZTA2-PT-44-40-179-W	G60076	E4631993	(01)10827002600765(17)271128(30)1(10)E4631993	28-nov-24	28-nov-27	1
ZTA2-PT-44-40-179-W	G60076	E4631994	(01)10827002600765(17)271201(30)1(10)E4631994	01-dec-24	01-dec-27	1
ZTA2-PT-44-40-179-W	G60076	E4631995	(01)10827002600765(17)271202(30)1(10)E4631995	02-dec-24	02-dec-27	1



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

REFERENCE PART NUMBER (RPN)	ORDER NUMBER (GPN)	LOT NUMBER	BOX UNIQUE DEVICE IDENTIFIER (UDI)	MANUFACTURE DATE	EXPIRATION DATE	QUANTITY
ZTA2-PT-44-40-179-W	G60076	E4631996	(01)10827002600765(17)271129(30)1(10)E4631996	29-nov-24	29-nov-27	1
ZTA2-PT-44-40-179-W	G60076	E4631997	(01)10827002600765(17)271201(30)1(10)E4631997	01-dec-24	01-dec-27	1
ZTA2-PT-44-40-233-W	G60084	E4631998	(01)10827002600840(17)271129(30)1(10)E4631998	29-nov-24	29-nov-27	1
ZTA2-PT-44-40-233-W	G60084	E4631999	(01)10827002600840(17)271128(30)1(10)E4631999	28-nov-24	28-nov-27	1
ZTA2-PT-44-40-233-W	G60084	E4632000	(01)10827002600840(17)271129(30)1(10)E4632000	29-nov-24	29-nov-27	1
ZTA2-P-40-117-W	G59915	E4634911	(01)10827002599151(17)271204(30)1(10)E4634911	04-dec-24	04-dec-27	1
ZTA2-P-40-117-W	G59915	E4634913	(01)10827002599151(17)271204(30)1(10)E4634913	04-dec-24	04-dec-27	1
ZTA2-P-40-117-W	G59915	E4634916	(01)10827002599151(17)271204(30)1(10)E4634916	04-dec-24	04-dec-27	1
ZTA2-P-40-117-W	G59915	E4634917	(01)10827002599151(17)271204(30)1(10)E4634917	04-dec-24	04-dec-27	1
ZTA2-P-40-167-W	G59921	E4634918	(01)10827002599212(17)271203(30)1(10)E4634918	03-dec-24	03-dec-27	1
ZTA2-P-40-167-W	G59921	E4634919	(01)10827002599212(17)271205(30)1(10)E4634919	05-dec-24	05-dec-27	1
ZTA2-P-40-167-W	G59921	E4634920	(01)10827002599212(17)271203(30)1(10)E4634920	03-dec-24	03-dec-27	1
ZTA2-P-42-121-W	G59926	E4634922	(01)10827002599267(17)271204(30)1(10)E4634922	04-dec-24	04-dec-27	1
ZTA2-P-42-225-W	G59935	E4634923	(01)10827002599359(17)271204(30)1(10)E4634923	04-dec-24	04-dec-27	1
ZTA2-P-44-125-W	G59938	E4634924	(01)10827002599380(17)271204(30)1(10)E4634924	04-dec-24	04-dec-27	1
ZTA2-P-44-125-W	G59938	E4634925	(01)10827002599380(17)271204(30)1(10)E4634925	04-dec-24	04-dec-27	1
ZTA2-P-44-125-W	G59938	E4634926	(01)10827002599380(17)271205(30)1(10)E4634926	05-dec-24	05-dec-27	1
ZTA2-P-44-179-W	G59944	E4634927	(01)10827002599441(17)271205(30)1(10)E4634927	05-dec-24	05-dec-27	1
ZTA2-P-44-179-W	G59944	E4634928	(01)10827002599441(17)271204(30)1(10)E4634928	04-dec-24	04-dec-27	1
ZTA2-P-46-125-W	G59950	E4634929	(01)10827002599502(17)271203(30)1(10)E4634929	03-dec-24	03-dec-27	1
ZTA2-P-46-125-W	G59950	E4634930	(01)10827002599502(17)271203(30)1(10)E4634930	03-dec-24	03-dec-27	1
ZTA2-P-46-125-W	G59950	E4634931	(01)10827002599502(17)271203(30)1(10)E4634931	03-dec-24	03-dec-27	1
ZTA2-P-46-179-W	G59910	E4634932	(01)10827002599106(17)271203(30)1(10)E4634932	03-dec-24	03-dec-27	1
ZTA2-P-46-179-W	G59910	E4634933	(01)10827002599106(17)271203(30)1(10)E4634933	03-dec-24	03-dec-27	1
ZTA2-P-46-233-W	G59912	E4634936	(01)10827002599120(17)271204(30)1(10)E4634936	04-dec-24	04-dec-27	1
ZTA2-PT-42-38-173-W	G60070	E4634937	(01)10827002600703(17)271204(30)1(10)E4634937	04-dec-24	04-dec-27	1
ZTA2-PT-42-38-173-W	G60070	E4634938	(01)10827002600703(17)271204(30)1(10)E4634938	04-dec-24	04-dec-27	1
ZTA2-PT-42-38-173-W	G60070	E4634939	(01)10827002600703(17)271204(30)1(10)E4634939	04-dec-24	04-dec-27	1
ZTA2-PT-42-38-173-W	G60070	E4634940	(01)10827002600703(17)271204(30)1(10)E4634940	04-dec-24	04-dec-27	1
ZTA2-PT-46-42-179-W	G60079	E4634941	(01)10827002600796(17)271203(30)1(10)E4634941	03-dec-24	03-dec-27	1
ZTA2-PT-46-42-179-W	G60079	E4634942	(01)10827002600796(17)271203(30)1(10)E4634942	03-dec-24	03-dec-27	1
ZTA2-PT-46-42-179-W	G60079	E4634943	(01)10827002600796(17)271205(30)1(10)E4634943	05-dec-24	05-dec-27	1
ZTA2-PT-46-42-179-W	G60079	E4634944	(01)10827002600796(17)271203(30)1(10)E4634944	03-dec-24	03-dec-27	1
ZTA2-PT-46-42-233-W	G60082	E4634947	(01)10827002600826(17)271204(30)1(10)E4634947	04-dec-24	04-dec-27	1
ZTA2-P-40-167-W	G59921	E4635174	(01)10827002599212(17)271204(30)1(10)E4635174	04-dec-24	04-dec-27	1
ZTA2-P-42-225-W	G59935	E4635175	(01)10827002599359(17)271204(30)1(10)E4635175	04-dec-24	04-dec-27	1
ZTA2-P-46-125-W	G59950	E4635177	(01)10827002599502(17)271204(30)1(10)E4635176	04-dec-24	04-dec-27	1
ZTA2-PT-40-36-217-W	G60067	E4635194	(01)10827002600673(17)271206(30)1(10)E4635194	06-dec-24	06-dec-27	1



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

REFERENCE PART NUMBER (RPN)	ORDER NUMBER (GPN)	LOT NUMBER	BOX UNIQUE DEVICE IDENTIFIER (UDI)	MANUFACTURE DATE	EXPIRATION DATE	QUANTITY
ZTA2-PT-40-36-217-W	G60067	E4635195	(01)10827002600673(17)271205(30)1(10)E4635195	05-dec-24	05-dec-27	1
ZTA2-PT-40-36-217-W	G60067	E4635196	(01)10827002600673(17)271204(30)1(10)E4635196	04-dec-24	04-dec-27	1
ZTA2-PT-40-36-217-W	G60067	E4635197	(01)10827002600673(17)271204(30)1(10)E4635197	04-dec-24	04-dec-27	1
ZTA2-PT-42-38-173-W	G60070	E4635198	(01)10827002600703(17)271204(30)1(10)E4635198	04-dec-24	04-dec-27	1
ZTA2-PT-42-38-225-W	G60073	E4635199	(01)10827002600734(17)271204(30)1(10)E4635199	04-dec-24	04-dec-27	1
ZTA2-PT-42-38-225-W	G60073	E4635200	(01)10827002600734(17)271204(30)1(10)E4635200	04-dec-24	04-dec-27	1
ZTA2-PT-42-38-225-W	G60073	E4635201	(01)10827002600734(17)271204(30)1(10)E4635201	04-dec-24	04-dec-27	1
ZTA2-PT-44-40-179-W	G60076	E4635202	(01)10827002600765(17)271204(30)1(10)E4635202	04-dec-24	04-dec-27	1
ZTA2-PT-44-40-179-W	G60076	E4635203	(01)10827002600765(17)271204(30)1(10)E4635203	04-dec-24	04-dec-27	1
ZTA2-PT-44-40-179-W	G60076	E4635204	(01)10827002600765(17)271204(30)1(10)E4635204	04-dec-24	04-dec-27	1
ZTA2-PT-44-40-233-W	G60084	E4635205	(01)10827002600840(17)271204(30)1(10)E4635205	04-dec-24	04-dec-27	1
ZTA2-P-40-117-W	G59915	E4646219	(01)10827002599151(17)280110(30)1(10)E4646219	10-jan-25	10-jan-28	1
ZTA2-P-40-167-W	G59921	E4646220	(01)10827002599212(17)280110(30)1(10)E4646220	10-jan-25	10-jan-28	1
ZTA2-P-40-217-W	G59923	E4646221	(01)10827002599236(17)280114(30)1(10)E4646221	14-jan-25	14-jan-28	1
ZTA2-P-42-173-W	G59932	E4646222	(01)10827002599328(17)280122(30)1(10)E4646222	22-jan-25	22-jan-28	1
ZTA2-P-42-173-W	G59932	E4646223	(01)10827002599328(17)280110(30)1(10)E4646223	10-jan-25	10-jan-28	1
ZTA2-P-42-173-W	G59932	E4646224	(01)10827002599328(17)280110(30)1(10)E4646224	10-jan-25	10-jan-28	1
ZTA2-P-42-173-W	G59932	E4646225	(01)10827002599328(17)280220(30)1(10)E4646225	20-feb-25	20-feb-28	1
ZTA2-P-42-173-W	G59932	E4646226	(01)10827002599328(17)280110(30)1(10)E4646226	10-jan-25	10-jan-28	1
ZTA2-P-42-225-W	G59935	E4646228	(01)10827002599359(17)280110(30)1(10)E4646228	10-jan-25	10-jan-28	1
ZTA2-P-42-225-W	G59935	E4646229	(01)10827002599359(17)280110(30)1(10)E4646229	10-jan-25	10-jan-28	1
ZTA2-P-42-225-W	G59935	E4646230	(01)10827002599359(17)280110(30)1(10)E4646230	10-jan-25	10-jan-28	1
ZTA2-P-44-179-W	G59944	E4646231	(01)10827002599441(17)280122(30)1(10)E4646231	22-jan-25	22-jan-28	1
ZTA2-P-44-233-W	G59947	E4646232	(01)10827002599472(17)280112(30)1(10)E4646232	12-jan-25	12-jan-28	1
ZTA2-P-44-233-W	G59947	E4646233	(01)10827002599472(17)280110(30)1(10)E4646233	10-jan-25	10-jan-28	1
ZTA2-P-44-233-W	G59947	E4646234	(01)10827002599472(17)280115(30)1(10)E4646234	15-jan-25	15-jan-28	1
ZTA2-P-44-233-W	G59947	E4646235	(01)10827002599472(17)280113(30)1(10)E4646235	13-jan-25	13-jan-28	1
ZTA2-P-44-233-W	G59947	E4646236	(01)10827002599472(17)280110(30)1(10)E4646236	10-jan-25	10-jan-28	1
ZTA2-P-46-125-W	G59950	E4646237	(01)10827002599502(17)280114(30)1(10)E4646237	14-jan-25	14-jan-28	1
ZTA2-P-46-125-W	G59950	E4646238	(01)10827002599502(17)280113(30)1(10)E4646238	13-jan-25	13-jan-28	1
ZTA2-P-46-179-W	G59910	E4646239	(01)10827002599106(17)280112(30)1(10)E4646239	12-jan-25	12-jan-28	1
ZTA2-P-46-179-W	G59910	E4646240	(01)10827002599106(17)280112(30)1(10)E4646240	12-jan-25	12-jan-28	1
ZTA2-P-46-179-W	G59910	E4646241	(01)10827002599106(17)280112(30)1(10)E4646241	12-jan-25	12-jan-28	1
ZTA2-P-46-179-W	G59910	E4646242	(01)10827002599106(17)280112(30)1(10)E4646242	12-jan-25	12-jan-28	1
ZTA2-P-46-233-W	G59912	E4646243	(01)10827002599120(17)280121(30)1(10)E4646243	21-jan-25	21-jan-28	1
ZTA2-P-46-233-W	G59912	E4646244	(01)10827002599120(17)280114(30)1(10)E4646244	14-jan-25	14-jan-28	1
ZTA2-P-46-233-W	G59912	E4646245	(01)10827002599120(17)280113(30)1(10)E4646245	13-jan-25	13-jan-28	1
ZTA2-PT-40-36-167-W	G60064	E4646246	(01)10827002600642(17)280113(30)1(10)E4646246	13-jan-25	13-jan-28	1



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

REFERENCE PART NUMBER (RPN)	ORDER NUMBER (GPN)	LOT NUMBER	BOX UNIQUE DEVICE IDENTIFIER (UDI)	MANUFACTURE DATE	EXPIRATION DATE	QUANTITY
ZTA2-PT-40-36-167-W	G60064	E4646247	(01)10827002600642(17)280113(30)1(10)E4646247	13-jan-25	13-jan-28	1
ZTA2-PT-42-38-173-W	G60070	E4646248	(01)10827002600703(17)280113(30)1(10)E4646248	13-jan-25	13-jan-28	1
ZTA2-PT-42-38-173-W	G60070	E4646250	(01)10827002600703(17)280113(30)1(10)E4646250	13-jan-25	13-jan-28	1
ZTA2-PT-42-38-225-W	G60073	E4646251	(01)10827002600734(17)280113(30)1(10)E4646251	13-jan-25	13-jan-28	1
ZTA2-PT-44-40-179-W	G60076	E4646252	(01)10827002600765(17)280123(30)1(10)E4646252	23-jan-25	23-jan-28	1
ZTA2-PT-44-40-179-W	G60076	E4646253	(01)10827002600765(17)280113(30)1(10)E4646253	13-jan-25	13-jan-28	1
ZTA2-PT-44-40-233-W	G60084	E4646254	(01)10827002600840(17)280114(30)1(10)E4646254	14-jan-25	14-jan-28	1
ZTA2-PT-44-40-233-W	G60084	E4646255	(01)10827002600840(17)280114(30)1(10)E4646255	14-jan-25	14-jan-28	1
ZTA2-PT-44-40-233-W	G60084	E4646257	(01)10827002600840(17)280114(30)1(10)E4646257	14-jan-25	14-jan-28	1
ZTA2-PT-44-40-233-W	G60084	E4646258	(01)10827002600840(17)280114(30)1(10)E4646258	14-jan-25	14-jan-28	1
ZTA2-PT-46-42-179-W	G60079	E4646259	(01)10827002600796(17)280110(30)1(10)E4646259	10-jan-25	10-jan-28	1
ZTA2-PT-46-42-233-W	G60082	E4646261	(01)10827002600826(17)280110(30)1(10)E4646261	10-jan-25	10-jan-28	1
ZTA2-P-40-117-W	G59915	E4656519	(01)10827002599151(17)280212(30)1(10)E4656519	12-feb-25	12-feb-28	1
ZTA2-P-40-117-W	G59915	E4656520	(01)10827002599151(17)280210(30)1(10)E4656520	10-feb-25	10-feb-28	1
ZTA2-P-40-117-W	G59915	E4656521	(01)10827002599151(17)280218(30)1(10)E4656521	18-feb-25	18-feb-28	1
ZTA2-P-40-167-W	G59921	E4656525	(01)10827002599212(17)280210(30)1(10)E4656525	10-feb-25	10-feb-28	1
ZTA2-P-40-167-W	G59921	E4656526	(01)10827002599212(17)280210(30)1(10)E4656526	10-feb-25	10-feb-28	1
ZTA2-P-40-167-W	G59921	E4656527	(01)10827002599212(17)280210(30)1(10)E4656527	10-feb-25	10-feb-28	1
ZTA2-P-40-167-W	G59921	E4656528	(01)10827002599212(17)280210(30)1(10)E4656528	10-feb-25	10-feb-28	1
ZTA2-P-40-167-W	G59921	E4656529	(01)10827002599212(17)280210(30)1(10)E4656529	10-feb-25	10-feb-28	1
ZTA2-P-42-121-W	G59926	E4656530	(01)10827002599267(17)280211(30)1(10)E4656530	11-feb-25	11-feb-28	1
ZTA2-P-42-121-W	G59926	E4656531	(01)10827002599267(17)280210(30)1(10)E4656531	10-feb-25	10-feb-28	1
ZTA2-P-42-121-W	G59926	E4656532	(01)10827002599267(17)280210(30)1(10)E4656532	10-feb-25	10-feb-28	1
ZTA2-P-46-125-W	G59950	E4656533	(01)10827002599502(17)280211(30)1(10)E4656533	11-feb-25	11-feb-28	1
ZTA2-P-46-125-W	G59950	E4656534	(01)10827002599502(17)280210(30)1(10)E4656534	10-feb-25	10-feb-28	1
ZTA2-P-46-125-W	G59950	E4656535	(01)10827002599502(17)280210(30)1(10)E4656535	10-feb-25	10-feb-28	1
ZTA2-P-46-179-W	G59910	E4656536	(01)10827002599106(17)280217(30)1(10)E4656536	17-feb-25	17-feb-28	1
ZTA2-P-46-179-W	G59910	E4656537	(01)10827002599106(17)280210(30)1(10)E4656537	10-feb-25	10-feb-28	1
ZTA2-PT-40-36-167-W	G60064	E4656538	(01)10827002600642(17)280210(30)1(10)E4656538	10-feb-25	10-feb-28	1
ZTA2-PT-40-36-217-W	G60067	E4656539	(01)10827002600673(17)280210(30)1(10)E4656539	10-feb-25	10-feb-28	1
ZTA2-PT-40-36-217-W	G60067	E4656540	(01)10827002600673(17)280213(30)1(10)E4656540	13-feb-25	13-feb-28	1
ZTA2-PT-42-38-173-W	G60070	E4656541	(01)10827002600703(17)280211(30)1(10)E4656541	11-feb-25	11-feb-28	1
ZTA2-PT-42-38-173-W	G60070	E4656542	(01)10827002600703(17)280210(30)1(10)E4656542	10-feb-25	10-feb-28	1
ZTA2-PT-42-38-173-W	G60070	E4656543	(01)10827002600703(17)280210(30)1(10)E4656543	10-feb-25	10-feb-28	1
ZTA2-PT-42-38-173-W	G60070	E4656544	(01)10827002600703(17)280210(30)1(10)E4656544	10-feb-25	10-feb-28	1
ZTA2-PT-46-42-179-W	G60079	E4656545	(01)10827002600796(17)280210(30)1(10)E4656545	10-feb-25	10-feb-28	1
ZTA2-P-40-117-W	G59915	E4659972	(01)10827002599151(17)280211(30)1(10)E4659972	11-feb-25	11-feb-28	1
ZTA2-P-40-117-W	G59915	E4659973	(01)10827002599151(17)280211(30)1(10)E4659973	11-feb-25	11-feb-28	1



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

REFERENCE PART NUMBER (RPN)	ORDER NUMBER (GPN)	LOT NUMBER	BOX UNIQUE DEVICE IDENTIFIER (UDI)	MANUFACTURE DATE	EXPIRATION DATE	QUANTITY
ZTA2-P-40-117-W	G59915	E4659974	(01)10827002599151(17)280212(30)1(10)E4659974	12-feb-25	12-feb-28	1
ZTA2-P-40-167-W	G59921	E4659975	(01)10827002599212(17)280211(30)1(10)E4659975	11-feb-25	11-feb-28	1
ZTA2-P-42-121-W	G59926	E4659976	(01)10827002599267(17)280227(30)1(10)E4659976	27-feb-25	27-feb-28	1
ZTA2-P-42-121-W	G59926	E4659977	(01)10827002599267(17)280211(30)1(10)E4659977	11-feb-25	11-feb-28	1
ZTA2-P-42-173-W	G59932	E4659979	(01)10827002599328(17)280217(30)1(10)E4659979	17-feb-25	17-feb-28	1
ZTA2-P-42-173-W	G59932	E4659980	(01)10827002599328(17)280212(30)1(10)E4659980	12-feb-25	12-feb-28	1
ZTA2-P-42-173-W	G59932	E4659981	(01)10827002599328(17)280211(30)1(10)E4659981	11-feb-25	11-feb-28	1
ZTA2-P-44-125-W	G59938	E4659982	(01)10827002599380(17)280211(30)1(10)E4659982	11-feb-25	11-feb-28	1
ZTA2-P-44-179-W	G59944	E4659983	(01)10827002599441(17)280212(30)1(10)E4659983	12-feb-25	12-feb-28	1
ZTA2-P-44-233-W	G59947	E4659984	(01)10827002599472(17)280214(30)1(10)E4659984	14-feb-25	14-feb-28	1
ZTA2-P-46-125-W	G59950	E4659985	(01)10827002599502(17)280226(30)1(10)E4659985	26-feb-25	26-feb-28	1
ZTA2-P-46-179-W	G59910	E4659987	(01)10827002599106(17)280226(30)1(10)E4659987	26-feb-25	26-feb-28	1
ZTA2-P-46-179-W	G59910	E4659988	(01)10827002599106(17)280212(30)1(10)E4659988	12-feb-25	12-feb-28	1
ZTA2-PT-40-36-167-W	G60064	E4659990	(01)10827002600642(17)280212(30)1(10)E4659990	12-feb-25	12-feb-28	1
ZTA2-PT-40-36-167-W	G60064	E4659991	(01)10827002600642(17)280212(30)1(10)E4659991	12-feb-25	12-feb-28	1
ZTA2-PT-40-36-167-W	G60064	E4659992	(01)10827002600642(17)280212(30)1(10)E4659992	12-feb-25	12-feb-28	1
ZTA2-PT-40-36-217-W	G60067	E4659993	(01)10827002600673(17)280217(30)1(10)E4659993	17-feb-25	17-feb-28	1
ZTA2-PT-42-38-173-W	G60070	E4659994	(01)10827002600703(17)280214(30)1(10)E4659994	14-feb-25	14-feb-28	1
ZTA2-PT-42-38-225-W	G60073	E4659995	(01)10827002600734(17)280212(30)1(10)E4659995	12-feb-25	12-feb-28	1
ZTA2-PT-42-38-225-W	G60073	E4659996	(01)10827002600734(17)280217(30)1(10)E4659996	17-feb-25	17-feb-28	1
ZTA2-PT-46-42-179-W	G60079	E4659997	(01)10827002600796(17)280218(30)1(10)E4659997	18-feb-25	18-feb-28	1
ZTA2-PT-46-42-179-W	G60079	E4659998	(01)10827002600796(17)280211(30)1(10)E4659998	11-feb-25	11-feb-28	1
ZTA2-PT-46-42-179-W	G60079	E4659999	(01)10827002600796(17)280212(30)1(10)E4659999	12-feb-25	12-feb-28	1
ZTA2-P-40-217-W	G59923	E4662229	(01)10827002599236(17)280217(30)1(10)E4662229	17-feb-25	17-feb-28	1
ZTA2-P-40-217-W	G59923	E4662230	(01)10827002599236(17)280214(30)1(10)E4662230	14-feb-25	14-feb-28	1
ZTA2-P-40-217-W	G59923	E4662231	(01)10827002599236(17)280214(30)1(10)E4662231	14-feb-25	14-feb-28	1
ZTA2-P-44-125-W	G59938	E4662232	(01)10827002599380(17)280220(30)1(10)E4662232	20-feb-25	20-feb-28	1
ZTA2-P-44-125-W	G59938	E4662233	(01)10827002599380(17)280214(30)1(10)E4662233	14-feb-25	14-feb-28	1
ZTA2-P-44-179-W	G59944	E4662234	(01)10827002599441(17)280214(30)1(10)E4662234	14-feb-25	14-feb-28	1
ZTA2-PT-40-36-167-W	G60064	E4662235	(01)10827002600642(17)280214(30)1(10)E4662235	14-feb-25	14-feb-28	1
ZTA2-PT-44-40-179-W	G60076	E4662236	(01)10827002600765(17)280214(30)1(10)E4662236	14-feb-25	14-feb-28	1
ZTA2-PT-46-42-233-W	G60082	E4662238	(01)10827002600826(17)280214(30)1(10)E4662238	14-feb-25	14-feb-28	1
ZTA2-P-40-117-W	G59915	E4665450	(01)10827002599151(17)280227(30)1(10)E4665450	27-feb-25	27-feb-28	1
ZTA2-P-40-167-W	G59921	E4665451	(01)10827002599212(17)280224(30)1(10)E4665451	24-feb-25	24-feb-28	1
ZTA2-P-40-167-W	G59921	E4665452	(01)10827002599212(17)280306(30)1(10)E4665452	06-mar-25	06-mar-28	1
ZTA2-P-40-167-W	G59921	E4665453	(01)10827002599212(17)280220(30)1(10)E4665453	20-feb-25	20-feb-28	1
ZTA2-P-40-217-W	G59923	E4665454	(01)10827002599236(17)280220(30)1(10)E4665454	20-feb-25	20-feb-28	1
ZTA2-P-40-217-W	G59923	E4665455	(01)10827002599236(17)280220(30)1(10)E4665455	20-feb-25	20-feb-28	1



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

REFERENCE PART NUMBER (RPN)	ORDER NUMBER (GPN)	LOT NUMBER	BOX UNIQUE DEVICE IDENTIFIER (UDI)	MANUFACTURE DATE	EXPIRATION DATE	QUANTITY
ZTA2-P-40-217-W	G59923	E4665456	(01)10827002599236(17)280220(30)1(10)E4665456	20-feb-25	20-feb-28	1
ZTA2-P-44-125-W	G59938	E4665458	(01)10827002599380(17)280224(30)1(10)E4665458	24-feb-25	24-feb-28	1
ZTA2-P-44-179-W	G59944	E4665459	(01)10827002599441(17)280220(30)1(10)E4665459	20-feb-25	20-feb-28	1
ZTA2-P-44-179-W	G59944	E4665460	(01)10827002599441(17)280226(30)1(10)E4665460	26-feb-25	26-feb-28	1
ZTA2-PT-40-36-167-W	G60064	E4665461	(01)10827002600642(17)280220(30)1(10)E4665461	20-feb-25	20-feb-28	1
ZTA2-PT-40-36-167-W	G60064	E4665462	(01)10827002600642(17)280220(30)1(10)E4665462	20-feb-25	20-feb-28	1
ZTA2-PT-40-36-167-W	G60064	E4665463	(01)10827002600642(17)280220(30)1(10)E4665463	20-feb-25	20-feb-28	1
ZTA2-PT-42-38-173-W	G60070	E4665464	(01)10827002600703(17)280224(30)1(10)E4665464	24-feb-25	24-feb-28	1
ZTA2-PT-42-38-173-W	G60070	E4665465	(01)10827002600703(17)280220(30)1(10)E4665465	20-feb-25	20-feb-28	1
ZTA2-PT-46-42-179-W	G60079	E4665466	(01)10827002600796(17)280317(30)1(10)E4665466	17-mar-25	17-mar-28	1
ZTA2-P-40-167-W	G59921	E4670244	(01)10827002599212(17)280303(30)1(10)E4670244	03-mar-25	03-mar-28	1
ZTA2-P-40-167-W	G59921	E4670245	(01)10827002599212(17)280303(30)1(10)E4670245	03-mar-25	03-mar-28	1
ZTA2-P-40-217-W	G59923	E4670246	(01)10827002599236(17)280228(30)1(10)E4670246	28-feb-25	28-feb-28	1
ZTA2-P-40-217-W	G59923	E4670247	(01)10827002599236(17)280303(30)1(10)E4670247	03-mar-25	03-mar-28	1
ZTA2-P-42-225-W	G59935	E4670248	(01)10827002599359(17)280228(30)1(10)E4670248	28-feb-25	28-feb-28	1
ZTA2-P-44-125-W	G59938	E4670249	(01)10827002599380(17)280303(30)1(10)E4670249	03-mar-25	03-mar-28	1
ZTA2-P-44-125-W	G59938	E4670250	(01)10827002599380(17)280303(30)1(10)E4670250	03-mar-25	03-mar-28	1
ZTA2-P-46-125-W	G59950	E4670251	(01)10827002599502(17)280303(30)1(10)E4670251	03-mar-25	03-mar-28	1
ZTA2-PT-40-36-167-W	G60064	E4670252	(01)10827002600642(17)280303(30)1(10)E4670252	03-mar-25	03-mar-28	1
ZTA2-PT-40-36-167-W	G60064	E4670253	(01)10827002600642(17)280306(30)1(10)E4670253	06-mar-25	06-mar-28	1
ZTA2-PT-40-36-167-W	G60064	E4670254	(01)10827002600642(17)280303(30)1(10)E4670254	03-mar-25	03-mar-28	1
ZTA2-PT-40-36-217-W	G60067	E4670255	(01)10827002600673(17)280303(30)1(10)E4670255	03-mar-25	03-mar-28	1
ZTA2-P-40-167-W	G59921	E4672425	(01)10827002599212(17)280306(30)1(10)E4672425	06-mar-25	06-mar-28	1
ZTA2-P-40-167-W	G59921	E4672426	(01)10827002599212(17)280305(30)1(10)E4672426	05-mar-25	05-mar-28	1
ZTA2-P-40-167-W	G59921	E4672427	(01)10827002599212(17)280306(30)1(10)E4672427	06-mar-25	06-mar-28	1
ZTA2-P-40-217-W	G59923	E4672431	(01)10827002599236(17)280305(30)1(10)E4672431	05-mar-25	05-mar-28	1
ZTA2-P-40-217-W	G59923	E4672432	(01)10827002599236(17)280305(30)1(10)E4672432	05-mar-25	05-mar-28	1
ZTA2-P-42-121-W	G59926	E4672433	(01)10827002599267(17)280306(30)1(10)E4672433	06-mar-25	06-mar-28	1
ZTA2-P-42-121-W	G59926	E4672434	(01)10827002599267(17)280306(30)1(10)E4672434	06-mar-25	06-mar-28	1
ZTA2-P-42-121-W	G59926	E4672435	(01)10827002599267(17)280306(30)1(10)E4672435	06-mar-25	06-mar-28	1
ZTA2-P-42-121-W	G59926	E4672436	(01)10827002599267(17)280306(30)1(10)E4672436	06-mar-25	06-mar-28	1
ZTA2-P-42-121-W	G59926	E4672437	(01)10827002599267(17)280313(30)1(10)E4672437	13-mar-25	13-mar-28	1
ZTA2-P-42-121-W	G59926	E4672438	(01)10827002599267(17)280306(30)1(10)E4672438	06-mar-25	06-mar-28	1
ZTA2-P-42-121-W	G59926	E4672439	(01)10827002599267(17)280306(30)1(10)E4672439	06-mar-25	06-mar-28	1
ZTA2-P-42-173-W	G59932	E4672440	(01)10827002599328(17)280311(30)1(10)E4672440	11-mar-25	11-mar-28	1
ZTA2-P-42-173-W	G59932	E4672441	(01)10827002599328(17)280305(30)1(10)E4672441	05-mar-25	05-mar-28	1
ZTA2-P-42-173-W	G59932	E4672442	(01)10827002599328(17)280306(30)1(10)E4672442	06-mar-25	06-mar-28	1
ZTA2-P-42-173-W	G59932	E4672443	(01)10827002599328(17)280305(30)1(10)E4672443	05-mar-25	05-mar-28	1



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

REFERENCE PART NUMBER (RPN)	ORDER NUMBER (GPN)	LOT NUMBER	BOX UNIQUE DEVICE IDENTIFIER (UDI)	MANUFACTURE DATE	EXPIRATION DATE	QUANTITY
ZTA2-P-42-225-W	G59935	E4672444	(01)10827002599359(17)280305(30)1(10)E4672444	05-mar-25	05-mar-28	1
ZTA2-P-44-125-W	G59938	E4672445	(01)10827002599380(17)280306(30)1(10)E4672445	06-mar-25	06-mar-28	1
ZTA2-P-44-125-W	G59938	E4672446	(01)10827002599380(17)280306(30)1(10)E4672446	06-mar-25	06-mar-28	1
ZTA2-P-44-179-W	G59944	E4672447	(01)10827002599441(17)280305(30)1(10)E4672447	05-mar-25	05-mar-28	1
ZTA2-P-46-125-W	G59950	E4672448	(01)10827002599502(17)280306(30)1(10)E4672448	06-mar-25	06-mar-28	1
ZTA2-P-46-125-W	G59950	E4672449	(01)10827002599502(17)280306(30)1(10)E4672449	06-mar-25	06-mar-28	1
ZTA2-P-46-125-W	G59950	E4672450	(01)10827002599502(17)280306(30)1(10)E4672450	06-mar-25	06-mar-28	1
ZTA2-P-46-179-W	G59910	E4672451	(01)10827002599106(17)280306(30)1(10)E4672451	06-mar-25	06-mar-28	1
ZTA2-P-46-179-W	G59910	E4672452	(01)10827002599106(17)280306(30)1(10)E4672452	06-mar-25	06-mar-28	1
ZTA2-P-46-179-W	G59910	E4672453	(01)10827002599106(17)280307(30)1(10)E4672453	07-mar-25	07-mar-28	1
ZTA2-P-46-179-W	G59910	E4672454	(01)10827002599106(17)280307(30)1(10)E4672454	07-mar-25	07-mar-28	1
ZTA2-P-46-179-W	G59910	E4672455	(01)10827002599106(17)280307(30)1(10)E4672455	07-mar-25	07-mar-28	1
ZTA2-PT-40-36-217-W	G60067	E4672463	(01)10827002600673(17)280307(30)1(10)E4672463	07-mar-25	07-mar-28	1
ZTA2-PT-42-38-173-W	G60070	E4672465	(01)10827002600703(17)280305(30)1(10)E4672465	05-mar-25	05-mar-28	1
ZTA2-PT-42-38-173-W	G60070	E4672466	(01)10827002600703(17)280312(30)1(10)E4672466	12-mar-25	12-mar-28	1
ZTA2-PT-42-38-173-W	G60070	E4672467	(01)10827002600703(17)280305(30)1(10)E4672467	05-mar-25	05-mar-28	1
ZTA2-PT-42-38-173-W	G60070	E4672468	(01)10827002600703(17)280307(30)1(10)E4672468	07-mar-25	07-mar-28	1
ZTA2-PT-42-38-173-W	G60070	E4672469	(01)10827002600703(17)280307(30)1(10)E4672469	07-mar-25	07-mar-28	1
ZTA2-PT-42-38-225-W	G60073	E4672470	(01)10827002600734(17)280306(30)1(10)E4672470	06-mar-25	06-mar-28	1
ZTA2-PT-44-40-179-W	G60076	E4672471	(01)10827002600765(17)280306(30)1(10)E4672471	06-mar-25	06-mar-28	1
ZTA2-PT-44-40-179-W	G60076	E4672472	(01)10827002600765(17)280312(30)1(10)E4672472	12-mar-25	12-mar-28	1
ZTA2-P-40-167-W	G59921	E4675854	(01)10827002599212(17)280320(30)1(10)E4675854	20-mar-25	20-mar-28	1
ZTA2-P-40-217-W	G59923	E4675858	(01)10827002599236(17)280319(30)1(10)E4675858	19-mar-25	19-mar-28	1
ZTA2-P-42-121-W	G59926	E4675859	(01)10827002599267(17)280313(30)1(10)E4675859	13-mar-25	13-mar-28	1
ZTA2-P-42-173-W	G59932	E4675860	(01)10827002599328(17)280313(30)1(10)E4675860	13-mar-25	13-mar-28	1
ZTA2-P-42-173-W	G59932	E4675861	(01)10827002599328(17)280313(30)1(10)E4675861	13-mar-25	13-mar-28	1
ZTA2-P-42-225-W	G59935	E4675862	(01)10827002599359(17)280314(30)1(10)E4675862	14-mar-25	14-mar-28	1
ZTA2-P-42-225-W	G59935	E4675863	(01)10827002599359(17)280320(30)1(10)E4675863	20-mar-25	20-mar-28	1
ZTA2-P-44-125-W	G59938	E4675864	(01)10827002599380(17)280313(30)1(10)E4675864	13-mar-25	13-mar-28	1
ZTA2-P-44-125-W	G59938	E4675865	(01)10827002599380(17)280314(30)1(10)E4675865	14-mar-25	14-mar-28	1
ZTA2-P-46-179-W	G59910	E4675868	(01)10827002599106(17)280314(30)1(10)E4675868	14-mar-25	14-mar-28	1
ZTA2-PT-40-36-167-W	G60064	E4675871	(01)10827002600642(17)280314(30)1(10)E4675871	14-mar-25	14-mar-28	1
ZTA2-PT-40-36-167-W	G60064	E4675872	(01)10827002600642(17)280313(30)1(10)E4675872	13-mar-25	13-mar-28	1
ZTA2-PT-42-38-173-W	G60070	E4675876	(01)10827002600703(17)280314(30)1(10)E4675876	14-mar-25	14-mar-28	1
ZTA2-PT-42-38-173-W	G60070	E4675877	(01)10827002600703(17)280313(30)1(10)E4675877	13-mar-25	13-mar-28	1
ZTA2-PT-42-38-225-W	G60073	E4675878	(01)10827002600734(17)280313(30)1(10)E4675878	13-mar-25	13-mar-28	1
ZTA2-PT-46-42-179-W	G60079	E4675879	(01)10827002600796(17)280321(30)1(10)E4675879	21-mar-25	21-mar-28	1
ZTA2-P-40-167-W	G59921	E4677105	(01)10827002599212(17)280318(30)1(10)E4677105	18-mar-25	18-mar-28	1



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

REFERENCE PART NUMBER (RPN)	ORDER NUMBER (GPN)	LOT NUMBER	BOX UNIQUE DEVICE IDENTIFIER (UDI)	MANUFACTURE DATE	EXPIRATION DATE	QUANTITY
ZTA2-P-42-121-W	G59926	E4677112	(01)10827002599267(17)280319(30)1(10)E4677112	19-mar-25	19-mar-28	1
ZTA2-P-42-121-W	G59926	E4677114	(01)10827002599267(17)280319(30)1(10)E4677114	19-mar-25	19-mar-28	1
ZTA2-P-42-121-W	G59926	E4677115	(01)10827002599267(17)280319(30)1(10)E4677115	19-mar-25	19-mar-28	1
ZTA2-P-42-121-W	G59926	E4677117	(01)10827002599267(17)280319(30)1(10)E4677117	19-mar-25	19-mar-28	1
ZTA2-P-44-125-W	G59938	E4677118	(01)10827002599380(17)280317(30)1(10)E4677118	17-mar-25	17-mar-28	1
ZTA2-P-44-125-W	G59938	E4677119	(01)10827002599380(17)280319(30)1(10)E4677119	19-mar-25	19-mar-28	1
ZTA2-P-44-125-W	G59938	E4677120	(01)10827002599380(17)280319(30)1(10)E4677120	19-mar-25	19-mar-28	1
ZTA2-P-44-125-W	G59938	E4677121	(01)10827002599380(17)280317(30)1(10)E4677121	17-mar-25	17-mar-28	1
ZTA2-P-44-179-W	G59944	E4677122	(01)10827002599441(17)280318(30)1(10)E4677122	18-mar-25	18-mar-28	1
ZTA2-P-44-179-W	G59944	E4677123	(01)10827002599441(17)280318(30)1(10)E4677123	18-mar-25	18-mar-28	1
ZTA2-P-44-179-W	G59944	E4677125	(01)10827002599441(17)280318(30)1(10)E4677125	18-mar-25	18-mar-28	1
ZTA2-P-44-179-W	G59944	E4677126	(01)10827002599441(17)280318(30)1(10)E4677126	18-mar-25	18-mar-28	1
ZTA2-P-44-179-W	G59944	E4677127	(01)10827002599441(17)280319(30)1(10)E4677127	19-mar-25	19-mar-28	1
ZTA2-P-46-125-W	G59950	E4677129	(01)10827002599502(17)280319(30)1(10)E4677129	19-mar-25	19-mar-28	1
ZTA2-PT-40-36-167-W	G60064	E4677137	(01)10827002600642(17)280317(30)1(10)E4677137	17-mar-25	17-mar-28	1
ZTA2-PT-40-36-167-W	G60064	E4677138	(01)10827002600642(17)280317(30)1(10)E4677138	17-mar-25	17-mar-28	1
ZTA2-PT-40-36-167-W	G60064	E4677139	(01)10827002600642(17)280319(30)1(10)E4677139	19-mar-25	19-mar-28	1
ZTA2-PT-46-42-179-W	G60079	E4677141	(01)10827002600796(17)280401(30)1(10)E4677141	01-apr-25	01-apr-28	1
ZTA2-PT-46-42-179-W	G60079	E4677142	(01)10827002600796(17)280319(30)1(10)E4677142	19-mar-25	19-mar-28	1
ZTA2-PT-46-42-179-W	G60079	E4677143	(01)10827002600796(17)280318(30)1(10)E4677143	18-mar-25	18-mar-28	1
ZTA2-P-40-117-W	G59915	E4679401	(01)10827002599151(17)280320(30)1(10)E4679401	20-mar-25	20-mar-28	1
ZTA2-P-40-117-W	G59915	E4679406	(01)10827002599151(17)280319(30)1(10)E4679406	19-mar-25	19-mar-28	1
ZTA2-P-40-167-W	G59921	E4679407	(01)10827002599212(17)280324(30)1(10)E4679407	24-mar-25	24-mar-28	1
ZTA2-P-40-167-W	G59921	E4679408	(01)10827002599212(17)280320(30)1(10)E4679408	20-mar-25	20-mar-28	1
ZTA2-P-40-167-W	G59921	E4679409	(01)10827002599212(17)280321(30)1(10)E4679409	21-mar-25	21-mar-28	1
ZTA2-P-40-167-W	G59921	E4679410	(01)10827002599212(17)280321(30)1(10)E4679410	21-mar-25	21-mar-28	1
ZTA2-P-40-167-W	G59921	E4679411	(01)10827002599212(17)280321(30)1(10)E4679411	21-mar-25	21-mar-28	1
ZTA2-P-40-167-W	G59921	E4679412	(01)10827002599212(17)280324(30)1(10)E4679412	24-mar-25	24-mar-28	1
ZTA2-P-40-167-W	G59921	E4679413	(01)10827002599212(17)280325(30)1(10)E4679413	25-mar-25	25-mar-28	1
ZTA2-P-40-167-W	G59921	E4679414	(01)10827002599212(17)280324(30)1(10)E4679414	24-mar-25	24-mar-28	1
ZTA2-P-40-167-W	G59921	E4679415	(01)10827002599212(17)280324(30)1(10)E4679415	24-mar-25	24-mar-28	1
ZTA2-P-42-173-W	G59932	E4679420	(01)10827002599328(17)280321(30)1(10)E4679420	21-mar-25	21-mar-28	1
ZTA2-P-42-225-W	G59935	E4679424	(01)10827002599359(17)280320(30)1(10)E4679424	20-mar-25	20-mar-28	1
ZTA2-P-44-125-W	G59938	E4679427	(01)10827002599380(17)280320(30)1(10)E4679427	20-mar-25	20-mar-28	1
ZTA2-P-44-233-W	G59947	E4679429	(01)10827002599472(17)280320(30)1(10)E4679429	20-mar-25	20-mar-28	1
ZTA2-P-44-233-W	G59947	E4679430	(01)10827002599472(17)280324(30)1(10)E4679430	24-mar-25	24-mar-28	1
ZTA2-PT-40-36-167-W	G60064	E4679432	(01)10827002600642(17)280324(30)1(10)E4679432	24-mar-25	24-mar-28	1
ZTA2-PT-40-36-167-W	G60064	E4679433	(01)10827002600642(17)280324(30)1(10)E4679433	24-mar-25	24-mar-28	1



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

REFERENCE PART NUMBER (RPN)	ORDER NUMBER (GPN)	LOT NUMBER	BOX UNIQUE DEVICE IDENTIFIER (UDI)	MANUFACTURE DATE	EXPIRATION DATE	QUANTITY
ZTA2-PT-40-36-167-W	G60064	E4679434	(01)10827002600642(17)280321(30)1(10)E4679434	21-mar-25	21-mar-28	1
ZTA2-PT-40-36-167-W	G60064	E4679435	(01)10827002600642(17)280321(30)1(10)E4679435	21-mar-25	21-mar-28	1
ZTA2-PT-40-36-167-W	G60064	E4679436	(01)10827002600642(17)280324(30)1(10)E4679436	24-mar-25	24-mar-28	1
ZTA2-PT-42-38-173-W	G60070	E4679439	(01)10827002600703(17)280324(30)1(10)E4679439	24-mar-25	24-mar-28	1
ZTA2-PT-42-38-173-W	G60070	E4679440	(01)10827002600703(17)280324(30)1(10)E4679440	24-mar-25	24-mar-28	1
ZTA2-PT-42-38-173-W	G60070	E4679441	(01)10827002600703(17)280324(30)1(10)E4679441	24-mar-25	24-mar-28	1
ZTA2-PT-42-38-173-W	G60070	E4679442	(01)10827002600703(17)280321(30)1(10)E4679442	21-mar-25	21-mar-28	1
ZTA2-PT-44-40-179-W	G60076	E4679443	(01)10827002600765(17)280324(30)1(10)E4679443	24-mar-25	24-mar-28	1
ZTA2-PT-44-40-233-W	G60084	E4679445	(01)10827002600840(17)280327(30)1(10)E4679445	27-mar-25	27-mar-28	1
ZTA2-PT-46-42-179-W	G60079	E4679446	(01)10827002600796(17)280324(30)1(10)E4679446	24-mar-25	24-mar-28	1
ZTA2-P-40-117-W	G59915	E4683233	(01)10827002599151(17)280403(30)1(10)E4683233	03-apr-25	03-apr-28	1
ZTA2-P-40-117-W	G59915	E4683236	(01)10827002599151(17)280401(30)1(10)E4683236	01-apr-25	01-apr-28	1
ZTA2-P-40-117-W	G59915	E4683237	(01)10827002599151(17)280401(30)1(10)E4683237	01-apr-25	01-apr-28	1
ZTA2-P-40-117-W	G59915	E4683238	(01)10827002599151(17)280401(30)1(10)E4683238	01-apr-25	01-apr-28	1
ZTA2-P-40-167-W	G59921	E4683243	(01)10827002599212(17)280402(30)1(10)E4683243	02-apr-25	02-apr-28	1
ZTA2-P-40-217-W	G59923	E4683246	(01)10827002599236(17)280402(30)1(10)E4683246	02-apr-25	02-apr-28	1
ZTA2-P-40-217-W	G59923	E4683247	(01)10827002599236(17)280424(30)1(10)E4683247	24-apr-25	24-apr-28	1
ZTA2-P-42-121-W	G59926	E4683249	(01)10827002599267(17)280404(30)1(10)E4683249	04-apr-25	04-apr-28	1
ZTA2-P-42-121-W	G59926	E4683251	(01)10827002599267(17)280402(30)1(10)E4683251	02-apr-25	02-apr-28	1
ZTA2-P-42-173-W	G59932	E4683254	(01)10827002599328(17)280401(30)1(10)E4683254	01-apr-25	01-apr-28	1
ZTA2-P-42-173-W	G59932	E4683255	(01)10827002599328(17)280401(30)1(10)E4683255	01-apr-25	01-apr-28	1
ZTA2-P-44-233-W	G59947	E4683261	(01)10827002599472(17)280331(30)1(10)E4683261	31-mar-25	31-mar-28	1
ZTA2-P-46-125-W	G59950	E4683262	(01)10827002599502(17)280424(30)1(10)E4683262	24-apr-25	24-apr-28	1
ZTA2-P-46-233-W	G59912	E4683270	(01)10827002599120(17)280401(30)1(10)E4683270	01-apr-25	01-apr-28	1
ZTA2-P-46-233-W	G59912	E4683272	(01)10827002599120(17)280401(30)1(10)E4683272	01-apr-25	01-apr-28	1
ZTA2-P-46-233-W	G59912	E4683273	(01)10827002599120(17)280403(30)1(10)E4683273	03-apr-25	03-apr-28	1
ZTA2-PT-40-36-167-W	G60064	E4683276	(01)10827002600642(17)280328(30)1(10)E4683276	28-mar-25	28-mar-28	1
ZTA2-PT-40-36-167-W	G60064	E4683278	(01)10827002600642(17)280401(30)1(10)E4683278	01-apr-25	01-apr-28	1
ZTA2-PT-40-36-167-W	G60064	E4683279	(01)10827002600642(17)280402(30)1(10)E4683279	02-apr-25	02-apr-28	1
ZTA2-PT-42-38-173-W	G60070	E4683282	(01)10827002600703(17)280401(30)1(10)E4683282	01-apr-25	01-apr-28	1
ZTA2-PT-42-38-173-W	G60070	E4683283	(01)10827002600703(17)280403(30)1(10)E4683283	03-apr-25	03-apr-28	1
ZTA2-PT-42-38-173-W	G60070	E4683284	(01)10827002600703(17)280401(30)1(10)E4683284	01-apr-25	01-apr-28	1
ZTA2-PT-42-38-173-W	G60070	E4683285	(01)10827002600703(17)280401(30)1(10)E4683285	01-apr-25	01-apr-28	1
ZTA2-PT-44-40-179-W	G60076	E4683286	(01)10827002600765(17)280401(30)1(10)E4683286	01-apr-25	01-apr-28	1
ZTA2-PT-46-42-179-W	G60079	E4683288	(01)10827002600796(17)280403(30)1(10)E4683288	03-apr-25	03-apr-28	1
ZTA2-PT-46-42-233-W	G60082	E4683289	(01)10827002600826(17)280402(30)1(10)E4683289	02-apr-25	02-apr-28	1
ZTA2-PT-46-42-233-W	G60082	E4683290	(01)10827002600826(17)280402(30)1(10)E4683290	02-apr-25	02-apr-28	1
ZTA2-P-40-117-W	G59915	E4684592	(01)10827002599151(17)280403(30)1(10)E4684592	03-apr-25	03-apr-28	1



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

REFERENCE PART NUMBER (RPN)	ORDER NUMBER (GPN)	LOT NUMBER	BOX UNIQUE DEVICE IDENTIFIER (UDI)	MANUFACTURE DATE	EXPIRATION DATE	QUANTITY
ZTA2-P-40-117-W	G59915	E4684593	(01)10827002599151(17)280403(30)1(10)E4684593	03-apr-25	03-apr-28	1
ZTA2-P-40-167-W	G59921	E4684597	(01)10827002599212(17)280403(30)1(10)E4684597	03-apr-25	03-apr-28	1
ZTA2-P-40-167-W	G59921	E4684598	(01)10827002599212(17)280403(30)1(10)E4684598	03-apr-25	03-apr-28	1
ZTA2-P-40-167-W	G59921	E4684601	(01)10827002599212(17)280404(30)1(10)E4684601	04-apr-25	04-apr-28	1
ZTA2-P-42-121-W	G59926	E4684603	(01)10827002599267(17)280422(30)1(10)E4684603	22-apr-25	22-apr-28	1
ZTA2-P-42-121-W	G59926	E4684604	(01)10827002599267(17)280403(30)1(10)E4684604	03-apr-25	03-apr-28	1
ZTA2-P-42-121-W	G59926	E4684606	(01)10827002599267(17)280403(30)1(10)E4684606	03-apr-25	03-apr-28	1
ZTA2-P-42-173-W	G59932	E4684608	(01)10827002599328(17)280403(30)1(10)E4684608	03-apr-25	03-apr-28	1
ZTA2-P-42-173-W	G59932	E4684609	(01)10827002599328(17)280403(30)1(10)E4684609	03-apr-25	03-apr-28	1
ZTA2-P-44-125-W	G59938	E4684610	(01)10827002599380(17)280403(30)1(10)E4684610	03-apr-25	03-apr-28	1
ZTA2-P-44-125-W	G59938	E4684612	(01)10827002599380(17)280403(30)1(10)E4684612	03-apr-25	03-apr-28	1
ZTA2-P-44-179-W	G59944	E4684614	(01)10827002599441(17)280404(30)1(10)E4684614	04-apr-25	04-apr-28	1
ZTA2-P-44-179-W	G59944	E4684616	(01)10827002599441(17)280404(30)1(10)E4684616	04-apr-25	04-apr-28	1
ZTA2-P-44-179-W	G59944	E4684618	(01)10827002599441(17)280404(30)1(10)E4684618	04-apr-25	04-apr-28	1
ZTA2-P-46-125-W	G59950	E4684619	(01)10827002599502(17)280403(30)1(10)E4684619	03-apr-25	03-apr-28	1
ZTA2-P-46-125-W	G59950	E4684620	(01)10827002599502(17)280403(30)1(10)E4684620	03-apr-25	03-apr-28	1
ZTA2-P-46-125-W	G59950	E4684621	(01)10827002599502(17)280402(30)1(10)E4684621	02-apr-25	02-apr-28	1
ZTA2-P-46-125-W	G59950	E4684624	(01)10827002599502(17)280402(30)1(10)E4684624	02-apr-25	02-apr-28	1
ZTA2-P-46-179-W	G59910	E4684627	(01)10827002599106(17)280403(30)1(10)E4684627	03-apr-25	03-apr-28	1
ZTA2-P-46-233-W	G59912	E4684633	(01)10827002599120(17)280403(30)1(10)E4684633	03-apr-25	03-apr-28	1
ZTA2-P-46-233-W	G59912	E4684635	(01)10827002599120(17)280402(30)1(10)E4684635	02-apr-25	02-apr-28	1
ZTA2-PT-40-36-167-W	G60064	E4684636	(01)10827002600642(17)280403(30)1(10)E4684636	03-apr-25	03-apr-28	1
ZTA2-PT-40-36-217-W	G60067	E4684638	(01)10827002600673(17)280422(30)1(10)E4684638	22-apr-25	22-apr-28	1
ZTA2-PT-42-38-173-W	G60070	E4684641	(01)10827002600703(17)280402(30)1(10)E4684641	02-apr-25	02-apr-28	1
ZTA2-PT-42-38-173-W	G60070	E4684642	(01)10827002600703(17)280402(30)1(10)E4684642	02-apr-25	02-apr-28	1
ZTA2-PT-42-38-173-W	G60070	E4684643	(01)10827002600703(17)280402(30)1(10)E4684643	02-apr-25	02-apr-28	1
ZTA2-PT-42-38-225-W	G60073	E4684644	(01)10827002600734(17)280402(30)1(10)E4684644	02-apr-25	02-apr-28	1
ZTA2-PT-42-38-225-W	G60073	E4684645	(01)10827002600734(17)280402(30)1(10)E4684645	02-apr-25	02-apr-28	1
ZTA2-PT-42-38-225-W	G60073	E4684646	(01)10827002600734(17)280402(30)1(10)E4684646	02-apr-25	02-apr-28	1
ZTA2-PT-44-40-179-W	G60076	E4684647	(01)10827002600765(17)280402(30)1(10)E4684647	02-apr-25	02-apr-28	1
ZTA2-PT-44-40-179-W	G60076	E4684648	(01)10827002600765(17)280402(30)1(10)E4684648	02-apr-25	02-apr-28	1
ZTA2-PT-46-42-179-W	G60079	E4684650	(01)10827002600796(17)280403(30)1(10)E4684650	03-apr-25	03-apr-28	1
ZTA2-PT-46-42-179-W	G60079	E4684651	(01)10827002600796(17)280416(30)1(10)E4684651	16-apr-25	16-apr-28	1
ZTA2-P-40-167-W	G59921	E4690174	(01)10827002599212(17)280429(30)1(10)E4690174	29-apr-25	29-apr-28	1
ZTA2-P-40-167-W	G59921	E4690176	(01)10827002599212(17)280428(30)1(10)E4690176	28-apr-25	28-apr-28	1
ZTA2-P-40-167-W	G59921	E4690178	(01)10827002599212(17)280415(30)1(10)E4690178	15-apr-25	15-apr-28	1
ZTA2-P-40-167-W	G59921	E4690180	(01)10827002599212(17)280411(30)1(10)E4690180	11-apr-25	11-apr-28	1
ZTA2-P-42-173-W	G59932	E4690182	(01)10827002599328(17)280414(30)1(10)E4690182	14-apr-25	14-apr-28	1



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

REFERENCE PART NUMBER (RPN)	ORDER NUMBER (GPN)	LOT NUMBER	BOX UNIQUE DEVICE IDENTIFIER (UDI)	MANUFACTURE DATE	EXPIRATION DATE	QUANTITY
ZTA2-P-42-173-W	G59932	E4690183	(01)10827002599328(17)280422(30)1(10)E4690183	22-apr-25	22-apr-28	1
ZTA2-P-42-173-W	G59932	E4690184	(01)10827002599328(17)280422(30)1(10)E4690184	22-apr-25	22-apr-28	1
ZTA2-P-42-225-W	G59935	E4690185	(01)10827002599359(17)280514(30)1(10)E4690185	14-maj-25	14-maj-28	1
ZTA2-P-42-225-W	G59935	E4690186	(01)10827002599359(17)280415(30)1(10)E4690186	15-apr-25	15-apr-28	1
ZTA2-P-46-125-W	G59950	E4690188	(01)10827002599502(17)280423(30)1(10)E4690188	23-apr-25	23-apr-28	1
ZTA2-P-46-125-W	G59950	E4690189	(01)10827002599502(17)280410(30)1(10)E4690189	10-apr-25	10-apr-28	1
ZTA2-P-46-179-W	G59910	E4690194	(01)10827002599106(17)280415(30)1(10)E4690194	15-apr-25	15-apr-28	1
ZTA2-PT-40-36-167-W	G60064	E4690195	(01)10827002600642(17)280422(30)1(10)E4690195	22-apr-25	22-apr-28	1
ZTA2-PT-40-36-167-W	G60064	E4690196	(01)10827002600642(17)280422(30)1(10)E4690196	22-apr-25	22-apr-28	1
ZTA2-PT-40-36-167-W	G60064	E4690197	(01)10827002600642(17)280422(30)1(10)E4690197	22-apr-25	22-apr-28	1
ZTA2-PT-40-36-167-W	G60064	E4690198	(01)10827002600642(17)280423(30)1(10)E4690198	23-apr-25	23-apr-28	1
ZTA2-PT-40-36-167-W	G60064	E4690200	(01)10827002600642(17)280414(30)1(10)E4690200	14-apr-25	14-apr-28	1
ZTA2-PT-40-36-167-W	G60064	E4690201	(01)10827002600642(17)280413(30)1(10)E4690201	13-apr-25	13-apr-28	1
ZTA2-PT-40-36-217-W	G60067	E4690202	(01)10827002600673(17)280422(30)1(10)E4690202	22-apr-25	22-apr-28	1
ZTA2-PT-42-38-225-W	G60073	E4690203	(01)10827002600734(17)280430(30)1(10)E4690203	30-apr-25	30-apr-28	1
ZTA2-PT-42-38-225-W	G60073	E4690204	(01)10827002600734(17)280430(30)1(10)E4690204	30-apr-25	30-apr-28	1
ZTA2-PT-42-38-225-W	G60073	E4690205	(01)10827002600734(17)280428(30)1(10)E4690205	28-apr-25	28-apr-28	1
ZTA2-PT-42-38-225-W	G60073	E4690206	(01)10827002600734(17)280414(30)1(10)E4690206	14-apr-25	14-apr-28	1
ZTA2-PT-44-40-179-W	G60076	E4690209	(01)10827002600765(17)280422(30)1(10)E4690209	22-apr-25	22-apr-28	1
ZTA2-PT-44-40-179-W	G60076	E4690211	(01)10827002600765(17)280414(30)1(10)E4690211	14-apr-25	14-apr-28	1
ZTA2-PT-44-40-179-W	G60076	E4690212	(01)10827002600765(17)280413(30)1(10)E4690212	13-apr-25	13-apr-28	1
ZTA2-PT-46-42-179-W	G60079	E4690213	(01)10827002600796(17)280422(30)1(10)E4690213	22-apr-25	22-apr-28	1
ZTA2-PT-46-42-233-W	G60082	E4690214	(01)10827002600826(17)280422(30)1(10)E4690214	22-apr-25	22-apr-28	1
ZTA2-P-42-121-W	G59926	E4693147	(01)10827002599267(17)280422(30)1(10)E4693147	22-apr-25	22-apr-28	1
ZTA2-P-42-121-W	G59926	E4693148	(01)10827002599267(17)280512(30)1(10)E4693148	12-maj-25	12-maj-28	1
ZTA2-P-42-173-W	G59932	E4693150	(01)10827002599328(17)280512(30)1(10)E4693150	12-maj-25	12-maj-28	1
ZTA2-P-42-173-W	G59932	E4693153	(01)10827002599328(17)280428(30)1(10)E4693153	28-apr-25	28-apr-28	1
ZTA2-PT-40-36-167-W	G60064	E4693157	(01)10827002600642(17)280507(30)1(10)E4693157	07-maj-25	07-maj-28	1
ZTA2-PT-44-40-179-W	G60076	E4693161	(01)10827002600765(17)280428(30)1(10)E4693161	28-apr-25	28-apr-28	1
ZTA2-PT-44-40-179-W	G60076	E4693162	(01)10827002600765(17)280514(30)1(10)E4693162	14-maj-25	14-maj-28	1
ZTA2-PT-46-42-179-W	G60079	E4693164	(01)10827002600796(17)280423(30)1(10)E4693164	23-apr-25	23-apr-28	1
ZTA2-PT-46-42-179-W	G60079	E4693165	(01)10827002600796(17)280423(30)1(10)E4693165	23-apr-25	23-apr-28	1
ZTA2-PT-46-42-179-W	G60079	E4693166	(01)10827002600796(17)280428(30)1(10)E4693166	28-apr-25	28-apr-28	1
ZTA2-PT-46-42-233-W	G60082	E4693170	(01)10827002600826(17)280429(30)1(10)E4693170	29-apr-25	29-apr-28	1
ZTA2-PT-46-42-233-W	G60082	E4693171	(01)10827002600826(17)280423(30)1(10)E4693171	23-apr-25	23-apr-28	1
ZTA2-P-40-217-W	G59923	E4695558	(01)10827002599236(17)280428(30)1(10)E4695558	28-apr-25	28-apr-28	1
ZTA2-P-46-125-W	G59950	E4695559	(01)10827002599502(17)280428(30)1(10)E4695559	28-apr-25	28-apr-28	1
ZTA2-P-46-125-W	G59950	E4695560	(01)10827002599502(17)280428(30)1(10)E4695560	28-apr-25	28-apr-28	1



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

REFERENCE PART NUMBER (RPN)	ORDER NUMBER (GPN)	LOT NUMBER	BOX UNIQUE DEVICE IDENTIFIER (UDI)	MANUFACTURE DATE	EXPIRATION DATE	QUANTITY
ZTA2-P-46-125-W	G59950	E4695561	(01)10827002599502(17)280428(30)1(10)E4695561	28-apr-25	28-apr-28	1
ZTA2-PT-46-42-179-W	G60079	E4695567	(01)10827002600796(17)280428(30)1(10)E4695567	28-apr-25	28-apr-28	1
ZTA2-P-40-117-W	G59915	E4701373	(01)10827002599151(17)280514(30)1(10)E4701373	14-maj-25	14-maj-28	1
ZTA2-P-40-117-W	G59915	E4701375	(01)10827002599151(17)280514(30)1(10)E4701375	14-maj-25	14-maj-28	1
ZTA2-P-40-117-W	G59915	E4701376	(01)10827002599151(17)280514(30)1(10)E4701376	14-maj-25	14-maj-28	1
ZTA2-P-40-117-W	G59915	E4701377	(01)10827002599151(17)280512(30)1(10)E4701377	12-maj-25	12-maj-28	1
ZTA2-P-40-117-W	G59915	E4701378	(01)10827002599151(17)280511(30)1(10)E4701378	11-maj-25	11-maj-28	1
ZTA2-P-40-117-W	G59915	E4701379	(01)10827002599151(17)280511(30)1(10)E4701379	11-maj-25	11-maj-28	1
ZTA2-P-40-117-W	G59915	E4701380	(01)10827002599151(17)280514(30)1(10)E4701380	14-maj-25	14-maj-28	1
ZTA2-P-40-167-W	G59921	E4701381	(01)10827002599212(17)280513(30)1(10)E4701381	13-maj-25	13-maj-28	1
ZTA2-P-40-167-W	G59921	E4701382	(01)10827002599212(17)280522(30)1(10)E4701382	22-maj-25	22-maj-28	1
ZTA2-P-40-167-W	G59921	E4701383	(01)10827002599212(17)280514(30)1(10)E4701383	14-maj-25	14-maj-28	1
ZTA2-P-40-167-W	G59921	E4701384	(01)10827002599212(17)280513(30)1(10)E4701384	13-maj-25	13-maj-28	1
ZTA2-P-40-167-W	G59921	E4701385	(01)10827002599212(17)280513(30)1(10)E4701385	13-maj-25	13-maj-28	1
ZTA2-P-40-167-W	G59921	E4701386	(01)10827002599212(17)280514(30)1(10)E4701386	14-maj-25	14-maj-28	1
ZTA2-P-42-121-W	G59926	E4701387	(01)10827002599267(17)280514(30)1(10)E4701387	14-maj-25	14-maj-28	1
ZTA2-P-42-121-W	G59926	E4701389	(01)10827002599267(17)280512(30)1(10)E4701389	12-maj-25	12-maj-28	1
ZTA2-P-42-173-W	G59932	E4701390	(01)10827002599328(17)280512(30)1(10)E4701390	12-maj-25	12-maj-28	1
ZTA2-P-44-125-W	G59938	E4701391	(01)10827002599380(17)280512(30)1(10)E4701391	12-maj-25	12-maj-28	1
ZTA2-P-44-125-W	G59938	E4701392	(01)10827002599380(17)280512(30)1(10)E4701392	12-maj-25	12-maj-28	1
ZTA2-P-44-125-W	G59938	E4701393	(01)10827002599380(17)280512(30)1(10)E4701393	12-maj-25	12-maj-28	1
ZTA2-P-44-125-W	G59938	E4701394	(01)10827002599380(17)280514(30)1(10)E4701394	14-maj-25	14-maj-28	1
ZTA2-P-44-179-W	G59944	E4701395	(01)10827002599441(17)280511(30)1(10)E4701395	11-maj-25	11-maj-28	1
ZTA2-P-44-179-W	G59944	E4701396	(01)10827002599441(17)280511(30)1(10)E4701396	11-maj-25	11-maj-28	1
ZTA2-P-44-179-W	G59944	E4701397	(01)10827002599441(17)280512(30)1(10)E4701397	12-maj-25	12-maj-28	1
ZTA2-P-44-179-W	G59944	E4701398	(01)10827002599441(17)280511(30)1(10)E4701398	11-maj-25	11-maj-28	1
ZTA2-PT-40-36-217-W	G60067	E4701399	(01)10827002600673(17)280522(30)1(10)E4701399	22-maj-25	22-maj-28	1
ZTA2-PT-40-36-217-W	G60067	E4701400	(01)10827002600673(17)280513(30)1(10)E4701400	13-maj-25	13-maj-28	1
ZTA2-PT-40-36-217-W	G60067	E4701401	(01)10827002600673(17)280514(30)1(10)E4701401	14-maj-25	14-maj-28	1
ZTA2-P-40-117-W	G59915	E4701402	(01)10827002599151(17)280512(30)1(10)E4701402	12-maj-25	12-maj-28	1
ZTA2-P-40-117-W	G59915	E4704235	(01)10827002599151(17)280519(30)1(10)E4704235	19-maj-25	19-maj-28	1
ZTA2-P-40-167-W	G59921	E4704236	(01)10827002599212(17)280515(30)1(10)E4704236	15-maj-25	15-maj-28	1
ZTA2-P-40-217-W	G59923	E4704238	(01)10827002599236(17)280515(30)1(10)E4704238	15-maj-25	15-maj-28	1
ZTA2-P-40-217-W	G59923	E4704239	(01)10827002599236(17)280519(30)1(10)E4704239	19-maj-25	19-maj-28	1
ZTA2-P-42-121-W	G59926	E4704240	(01)10827002599267(17)280515(30)1(10)E4704240	15-maj-25	15-maj-28	1
ZTA2-P-42-121-W	G59926	E4704241	(01)10827002599267(17)280519(30)1(10)E4704241	19-maj-25	19-maj-28	1
ZTA2-P-42-173-W	G59932	E4704242	(01)10827002599328(17)280515(30)1(10)E4704242	15-maj-25	15-maj-28	1
ZTA2-P-42-173-W	G59932	E4704244	(01)10827002599328(17)280515(30)1(10)E4704244	15-maj-25	15-maj-28	1



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

REFERENCE PART NUMBER (RPN)	ORDER NUMBER (GPN)	LOT NUMBER	BOX UNIQUE DEVICE IDENTIFIER (UDI)	MANUFACTURE DATE	EXPIRATION DATE	QUANTITY
ZTA2-P-44-233-W	G59947	E4704248	(01)10827002599472(17)280522(30)1(10)E4704248	22-maj-25	22-maj-28	1
ZTA2-PT-40-36-167-W	G60064	E4704249	(01)10827002600642(17)280515(30)1(10)E4704249	15-maj-25	15-maj-28	1
ZTA2-PT-40-36-167-W	G60064	E4704250	(01)10827002600642(17)280515(30)1(10)E4704250	15-maj-25	15-maj-28	1
ZTA2-PT-40-36-167-W	G60064	E4704251	(01)10827002600642(17)280519(30)1(10)E4704251	19-maj-25	19-maj-28	1
ZTA2-PT-40-36-217-W	G60067	E4704252	(01)10827002600673(17)280522(30)1(10)E4704252	22-maj-25	22-maj-28	1
ZTA2-PT-40-36-217-W	G60067	E4704253	(01)10827002600673(17)280515(30)1(10)E4704253	15-maj-25	15-maj-28	1
ZTA2-PT-42-38-173-W	G60070	E4704255	(01)10827002600703(17)280515(30)1(10)E4704255	15-maj-25	15-maj-28	1
ZTA2-PT-42-38-173-W	G60070	E4704256	(01)10827002600703(17)280515(30)1(10)E4704256	15-maj-25	15-maj-28	1
ZTA2-PT-42-38-173-W	G60070	E4704257	(01)10827002600703(17)280515(30)1(10)E4704257	15-maj-25	15-maj-28	1
ZTA2-PT-42-38-173-W	G60070	E4704258	(01)10827002600703(17)280520(30)1(10)E4704258	20-maj-25	20-maj-28	1
ZTA2-PT-42-38-173-W	G60070	E4704259	(01)10827002600703(17)280519(30)1(10)E4704259	19-maj-25	19-maj-28	1
ZTA2-P-42-121-W	G59926	E4707614	(01)10827002599267(17)280523(30)1(10)E4707614	23-maj-25	23-maj-28	1
ZTA2-P-42-225-W	G59935	E4707620	(01)10827002599359(17)280526(30)1(10)E4707620	26-maj-25	26-maj-28	1
ZTA2-P-44-233-W	G59947	E4707621	(01)10827002599472(17)280526(30)1(10)E4707621	26-maj-25	26-maj-28	1
ZTA2-P-46-179-W	G59910	E4707623	(01)10827002599106(17)280522(30)1(10)E4707623	22-maj-25	22-maj-28	1
ZTA2-P-46-179-W	G59910	E4707625	(01)10827002599106(17)280522(30)1(10)E4707625	22-maj-25	22-maj-28	1
ZTA2-PT-40-36-167-W	G60064	E4707626	(01)10827002600642(17)280601(30)1(10)E4707626	01-jun-25	01-jun-28	1
ZTA2-PT-40-36-167-W	G60064	E4707627	(01)10827002600642(17)280522(30)1(10)E4707627	22-maj-25	22-maj-28	1
ZTA2-PT-40-36-167-W	G60064	E4707628	(01)10827002600642(17)280522(30)1(10)E4707628	22-maj-25	22-maj-28	1
ZTA2-PT-40-36-167-W	G60064	E4707629	(01)10827002600642(17)280522(30)1(10)E4707629	22-maj-25	22-maj-28	1
ZTA2-PT-40-36-217-W	G60067	E4707631	(01)10827002600673(17)280601(30)1(10)E4707631	01-jun-25	01-jun-28	1
ZTA2-PT-42-38-173-W	G60070	E4707632	(01)10827002600703(17)280522(30)1(10)E4707632	22-maj-25	22-maj-28	1
ZTA2-PT-42-38-225-W	G60073	E4707636	(01)10827002600734(17)280522(30)1(10)E4707636	22-maj-25	22-maj-28	1
ZTA2-PT-42-38-225-W	G60073	E4707637	(01)10827002600734(17)280522(30)1(10)E4707637	22-maj-25	22-maj-28	1
ZTA2-PT-44-40-179-W	G60076	E4707638	(01)10827002600765(17)280523(30)1(10)E4707638	23-maj-25	23-maj-28	1
ZTA2-PT-44-40-179-W	G60076	E4707639	(01)10827002600765(17)280522(30)1(10)E4707639	22-maj-25	22-maj-28	1
ZTA2-PT-44-40-179-W	G60076	E4707640	(01)10827002600765(17)280522(30)1(10)E4707640	22-maj-25	22-maj-28	1
ZTA2-P-42-173-W	G59932	E4710863	(01)10827002599328(17)280602(30)1(10)E4710863	02-jun-25	02-jun-28	1
ZTA2-PT-40-36-167-W	G60064	E4710866	(01)10827002600642(17)280601(30)1(10)E4710866	01-jun-25	01-jun-28	1
ZTA2-PT-40-36-217-W	G60067	E4710867	(01)10827002600673(17)280602(30)1(10)E4710867	02-jun-25	02-jun-28	1
ZTA2-PT-42-38-173-W	G60070	E4710868	(01)10827002600703(17)280602(30)1(10)E4710868	02-jun-25	02-jun-28	1
ZTA2-PT-42-38-225-W	G60073	E4710869	(01)10827002600734(17)280602(30)1(10)E4710869	02-jun-25	02-jun-28	1
ZTA2-PT-42-38-225-W	G60073	E4710870	(01)10827002600734(17)280602(30)1(10)E4710870	02-jun-25	02-jun-28	1
ZTA2-PT-42-38-225-W	G60073	E4710871	(01)10827002600734(17)280602(30)1(10)E4710871	02-jun-25	02-jun-28	1
ZTA2-PT-46-42-179-W	G60079	E4710875	(01)10827002600796(17)280602(30)1(10)E4710875	02-jun-25	02-jun-28	1
ZTA2-PT-46-42-179-W	G60079	E4710876	(01)10827002600796(17)280601(30)1(10)E4710876	01-jun-25	01-jun-28	1
ZTA2-PT-46-42-179-W	G60079	E4710877	(01)10827002600796(17)280602(30)1(10)E4710877	02-jun-25	02-jun-28	1
ZTA2-PT-46-42-179-W	G60079	E4710878	(01)10827002600796(17)280602(30)1(10)E4710878	02-jun-25	02-jun-28	1



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

REFERENCE PART NUMBER (RPN)	ORDER NUMBER (GPN)	LOT NUMBER	BOX UNIQUE DEVICE IDENTIFIER (UDI)	MANUFACTURE DATE	EXPIRATION DATE	QUANTITY
ZTA2-PT-46-42-179-W	G60079	E4710879	(01)10827002600796(17)280601(30)1(10)E4710879	01-jun-25	01-jun-28	1
ZTA2-PT-46-42-233-W	G60082	E4710880	(01)10827002600826(17)280601(30)1(10)E4710880	01-jun-25	01-jun-28	1
ZTA2-PT-40-36-217-W	G60067	E4712861	(01)10827002600673(17)280610(30)1(10)E4712861	10-jun-25	10-jun-28	1
ZTA2-PT-42-38-173-W	G60070	E4712862	(01)10827002600703(17)280611(30)1(10)E4712862	11-jun-25	11-jun-28	1
ZTA2-PT-42-38-173-W	G60070	E4712863	(01)10827002600703(17)280606(30)1(10)E4712863	06-jun-25	06-jun-28	1
ZTA2-PT-42-38-173-W	G60070	E4712864	(01)10827002600703(17)280604(30)1(10)E4712864	04-jun-25	04-jun-28	1
ZTA2-PT-46-42-179-W	G60079	E4712865	(01)10827002600796(17)280604(30)1(10)E4712865	04-jun-25	04-jun-28	1
ZTA2-PT-46-42-179-W	G60079	E4712866	(01)10827002600796(17)280610(30)1(10)E4712866	10-jun-25	10-jun-28	1
ZTA2-PT-46-42-179-W	G60079	E4712867	(01)10827002600796(17)280610(30)1(10)E4712867	10-jun-25	10-jun-28	1
ZTA2-P-40-167	G59920	E4715752	(01)10827002599205(17)280622(30)1(10)E4715752	22-jun-25	22-jun-28	1
ZTA2-P-42-173	G59931	E4715753	(01)10827002599311(17)280616(30)1(10)E4715753	16-jun-25	16-jun-28	1
ZTA2-P-42-173	G59931	E4715754	(01)10827002599311(17)280613(30)1(10)E4715754	13-jun-25	13-jun-28	1
ZTA2-P-42-173	G59931	E4715755	(01)10827002599311(17)280616(30)1(10)E4715755	16-jun-25	16-jun-28	1
ZTA2-P-46-125	G59949	E4715756	(01)10827002599496(17)280613(30)1(10)E4715756	13-jun-25	13-jun-28	1
ZTA2-P-46-179	G59907	E4715757	(01)10827002599076(17)280616(30)1(10)E4715757	16-jun-25	16-jun-28	1
ZTA2-P-40-167	G59920	E4718032	(01)10827002599205(17)280619(30)1(10)E4718032	19-jun-25	19-jun-28	1
ZTA2-P-40-167	G59920	E4718033	(01)10827002599205(17)280624(30)1(10)E4718033	24-jun-25	24-jun-28	1
ZTA2-P-40-167	G59920	E4718034	(01)10827002599205(17)280619(30)1(10)E4718034	19-jun-25	19-jun-28	1
ZTA2-P-40-167	G59920	E4718035	(01)10827002599205(17)280619(30)1(10)E4718035	19-jun-25	19-jun-28	1
ZTA2-P-42-121	G59925	E4718044	(01)10827002599250(17)280619(30)1(10)E4718044	19-jun-25	19-jun-28	1
ZTA2-P-44-233	G59946	E4718047	(01)10827002599465(17)280703(30)1(10)E4718047	03-jul-25	03-jul-28	1
ZTA2-P-44-233	G59946	E4718048	(01)10827002599465(17)280624(30)1(10)E4718048	24-jun-25	24-jun-28	1
ZTA2-P-44-233	G59946	E4718049	(01)10827002599465(17)280624(30)1(10)E4718049	24-jun-25	24-jun-28	1
ZTA2-P-46-125	G59949	E4718052	(01)10827002599496(17)280619(30)1(10)E4718052	19-jun-25	19-jun-28	1
ZTA2-P-46-179	G59907	E4718053	(01)10827002599076(17)280619(30)1(10)E4718053	19-jun-25	19-jun-28	1
ZTA2-P-46-179	G59907	E4718055	(01)10827002599076(17)280619(30)1(10)E4718055	19-jun-25	19-jun-28	1
ZTA2-P-46-179	G59907	E4718056	(01)10827002599076(17)280624(30)1(10)E4718056	24-jun-25	24-jun-28	1
ZTA2-P-46-179	G59907	E4718057	(01)10827002599076(17)280619(30)1(10)E4718057	19-jun-25	19-jun-28	1
ZTA2-PT-40-36-167	G60063	E4718064	(01)10827002600635(17)280619(30)1(10)E4718064	19-jun-25	19-jun-28	1
ZTA2-PT-40-36-167	G60063	E4718065	(01)10827002600635(17)280619(30)1(10)E4718065	19-jun-25	19-jun-28	1
ZTA2-PT-42-38-173	G60069	E4718066	(01)10827002600697(17)280619(30)1(10)E4718066	19-jun-25	19-jun-28	1
ZTA2-P-42-173	G59931	E4721401	(01)10827002599311(17)280625(30)1(10)E4721401	25-jun-25	25-jun-28	1
ZTA2-P-42-173	G59931	E4721402	(01)10827002599311(17)280625(30)1(10)E4721402	25-jun-25	25-jun-28	1
ZTA2-P-42-173	G59931	E4721403	(01)10827002599311(17)280626(30)1(10)E4721403	26-jun-25	26-jun-28	1
ZTA2-PT-40-36-217	G60066	E4721425	(01)10827002600666(17)280626(30)1(10)E4721425	26-jun-25	26-jun-28	1
ZTA2-PT-40-36-217	G60066	E4721427	(01)10827002600666(17)280702(30)1(10)E4721427	02-jul-25	02-jul-28	1
ZTA2-PT-40-36-217	G60066	E4721428	(01)10827002600666(17)280626(30)1(10)E4721428	26-jun-25	26-jun-28	1
ZTA2-PT-42-38-173	G60069	E4721429	(01)10827002600697(17)280624(30)1(10)E4721429	24-jun-25	24-jun-28	1



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

REFERENCE PART NUMBER (RPN)	ORDER NUMBER (GPN)	LOT NUMBER	BOX UNIQUE DEVICE IDENTIFIER (UDI)	MANUFACTURE DATE	EXPIRATION DATE	QUANTITY
ZTA2-PT-42-38-173	G60069	E4721432	(01)10827002600697(17)280625(30)1(10)E4721432	25-jun-25	25-jun-28	1
ZTA2-PT-42-38-225	G60072	E4721435	(01)10827002600727(17)280624(30)1(10)E4721435	24-jun-25	24-jun-28	1
ZTA2-PT-42-38-225	G60072	E4721436	(01)10827002600727(17)280624(30)1(10)E4721436	24-jun-25	24-jun-28	1
ZTA2-PT-42-38-225	G60072	E4721438	(01)10827002600727(17)280624(30)1(10)E4721438	24-jun-25	24-jun-28	1
ZTA2-PT-42-38-225	G60072	E4721439	(01)10827002600727(17)280624(30)1(10)E4721439	24-jun-25	24-jun-28	1
ZTA2-PT-46-42-179	G60078	E4721440	(01)10827002600789(17)280702(30)1(10)E4721440	02-jul-25	02-jul-28	1
ZTA2-PT-46-42-179	G60078	E4721441	(01)10827002600789(17)280702(30)1(10)E4721441	02-jul-25	02-jul-28	1
ZTA2-PT-46-42-179	G60078	E4721443	(01)10827002600789(17)280703(30)1(10)E4721443	03-jul-25	03-jul-28	1
ZTA2-PT-46-42-179	G60078	E4721445	(01)10827002600789(17)280702(30)1(10)E4721445	02-jul-25	02-jul-28	1
ZTA2-PT-46-42-179	G60078	E4721446	(01)10827002600789(17)280703(30)1(10)E4721446	03-jul-25	03-jul-28	1
ZTA2-PT-46-42-179	G60078	E4721447	(01)10827002600789(17)280625(30)1(10)E4721447	25-jun-25	25-jun-28	1
ZTA2-PT-46-42-179	G60078	E4721448	(01)10827002600789(17)280702(30)1(10)E4721448	02-jul-25	02-jul-28	1
ZTA2-PT-46-42-233	G60081	E4721450	(01)10827002600819(17)280626(30)1(10)E4721450	26-jun-25	26-jun-28	1
ZTA2-PT-46-42-233	G60081	E4721452	(01)10827002600819(17)280624(30)1(10)E4721452	24-jun-25	24-jun-28	1
ZTA2-PT-46-42-233	G60081	E4721453	(01)10827002600819(17)280624(30)1(10)E4721453	24-jun-25	24-jun-28	1
ZTA2-PT-46-42-233	G60081	E4721454	(01)10827002600819(17)280624(30)1(10)E4721454	24-jun-25	24-jun-28	1
ZTA2-PT-46-42-233	G60081	E4721455	(01)10827002600819(17)280624(30)1(10)E4721455	24-jun-25	24-jun-28	1
ZTA2-PT-46-42-233	G60081	E4721456	(01)10827002600819(17)280624(30)1(10)E4721456	24-jun-25	24-jun-28	1
ZTA2-PT-46-42-233	G60081	E4721457	(01)10827002600819(17)280624(30)1(10)E4721457	24-jun-25	24-jun-28	1
ZTA2-PT-46-42-233	G60081	E4721458	(01)10827002600819(17)280624(30)1(10)E4721458	24-jun-25	24-jun-28	1
ZTA2-P-42-225	G59934	E4724678	(01)10827002599342(17)280708(30)1(10)E4724678	08-jul-25	08-jul-28	1
ZTA2-P-42-225	G59934	E4724679	(01)10827002599342(17)280708(30)1(10)E4724679	08-jul-25	08-jul-28	1
ZTA2-P-42-225	G59934	E4724680	(01)10827002599342(17)280708(30)1(10)E4724680	08-jul-25	08-jul-28	1
ZTA2-P-42-225	G59934	E4724681	(01)10827002599342(17)280708(30)1(10)E4724681	08-jul-25	08-jul-28	1
ZTA2-PT-44-40-179	G60075	E4724723	(01)10827002600758(17)280708(30)1(10)E4724723	08-jul-25	08-jul-28	1
ZTA2-PT-44-40-179	G60075	E4724724	(01)10827002600758(17)280708(30)1(10)E4724724	08-jul-25	08-jul-28	1
ZTA2-PT-44-40-179	G60075	E4724725	(01)10827002600758(17)280708(30)1(10)E4724725	08-jul-25	08-jul-28	1
ZTA2-PT-44-40-179	G60075	E4724726	(01)10827002600758(17)280708(30)1(10)E4724726	08-jul-25	08-jul-28	1
ZTA2-P-44-125	G59937	E4727703	(01)10827002599373(17)280711(30)1(10)E4727703	11-jul-25	11-jul-28	1
ZTA2-P-44-125	G59937	E4727704	(01)10827002599373(17)280711(30)1(10)E4727704	11-jul-25	11-jul-28	1
ZTA2-P-44-179	G59943	E4727711	(01)10827002599434(17)280717(30)1(10)E4727711	17-jul-25	17-jul-28	1
ZTA2-P-44-179	G59943	E4727712	(01)10827002599434(17)280717(30)1(10)E4727712	17-jul-25	17-jul-28	1
ZTA2-P-40-167-W	G59921	E4627443	(01)10827002599212(17)271115(30)1(10)E4627443	15-nov-24	15-nov-27	1
ZTA2-P-42-173-W	G59932	E4628837	(01)10827002599328(17)271121(30)1(10)E4628837	21-nov-24	21-nov-27	1
ZTA2-PT-46-42-179-W	G60079	E4634945	(01)10827002600796(17)271203(30)1(10)E4634945	03-dec-24	03-dec-27	1
ZTA2-PT-42-38-173-W	G60070	E4675875	(01)10827002600703(17)280314(30)1(10)E4675875	14-mar-25	14-mar-28	1
ZTA2-PT-44-40-179-W	G60076	E4662237	(01)10827002600765(17)280214(30)1(10)E4662237	14-feb-25	14-feb-28	1
ZTA2-P-42-173-W	G59932	E4646227	(01)10827002599328(17)280110(30)1(10)E4646227	10-jan-25	10-jan-28	1



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

REFERENCE PART NUMBER (RPN)	ORDER NUMBER (GPN)	LOT NUMBER	BOX UNIQUE DEVICE IDENTIFIER (UDI)	MANUFACTURE DATE	EXPIRATION DATE	QUANTITY
ZTA2-P-40-167-W	G59921	E4635173	(01)10827002599212(17)271204(30)1(10)E4635173	04-dec-24	04-dec-27	1
ZTA2-P-42-121-W	G59926	E4634921	(01)10827002599267(17)271204(30)1(10)E4634921	04-dec-24	04-dec-27	1
ZTA2-P-40-117-W	G59915	E4634915	(01)10827002599151(17)271204(30)1(10)E4634915	04-dec-24	04-dec-27	1
ZTA2-P-46-125-W	G59950	E4628847	(01)10827002599502(17)271119(30)1(10)E4628847	19-nov-24	19-nov-27	1
ZTA2-P-46-125-W	G59950	E4635176	(01)10827002599502(17)271204(30)1(10)E4635176	04-dec-24	04-dec-27	1
ZTA2-PT-46-42-179-W	G60079	E4656546	(01)10827002600796(17)280210(30)1(10)E4656546	10-feb-25	10-feb-28	1
ZTA2-P-40-217-W	G59923	E4672428	(01)10827002599236(17)280306(30)1(10)E4672428	06-mar-25	06-mar-28	1
ZTA2-PT-40-36-217-W	G60067	E4672464	(01)10827002600673(17)280307(30)1(10)E4672464	07-mar-25	07-mar-28	1
ZTA2-P-46-233-W	G59912	E4683269	(01)10827002599120(17)280401(30)1(10)E4683269	01-apr-25	01-apr-28	1
ZTA2-P-46-179-W	G59910	E4684625	(01)10827002599106(17)280403(30)1(10)E4684625	03-apr-25	03-apr-28	1
ZTA2-PT-40-36-217-W	G60067	E4704254	(01)10827002600673(17)280515(30)1(10)E4704254	15-maj-25	15-maj-28	1
ZTA2-PT-40-36-217-W	G60067	E4712860	(01)10827002600673(17)280612(30)1(10)E4712860	12-jun-25	12-jun-28	1



Risk Management Department

SEP 24 2025

Received

COOK MEDICAL LLC
 400 DANIELS WAY PO BOX 4195
 BLOOMINGTON, IN 47402-4195
 PHONE: 812.339.2235 TOLL FREE: 800.457.4500 FAX: 800.554.8335
 WWW.COOKMEDICAL.COM

URGENT: PROMPT RESPONSE REQUIRED

Acknowledgement and Receipt Form

You are required to fill out and return this form.

Please return it within 5 business days of receipt, even if you do not have any of the affected product(s).

Zenith Alpha® 2 Thoracic Endovascular Graft

Your Customer Information

C15623-0
 THE MOUNT SINAI HOSP
 Attn: Risk Management
 101 ST/BETWEEN 5TH & MADISON
 RECEIVING DEPT
 NEW YORK, NY 10029

Actions You Have Taken

1. I have read and understand the recall instructions provided in the letter dated September 18th, 2025. ☒ Yes ☐ No

Additional Actions for Distributors

Share the "URGENT: MEDICAL DEVICE RECALL" letter with appropriate personnel within your organization or with any organization where the affected product(s) have been transferred.

Please confirm the following:

1. I have notified the appropriate personnel of the recall. ☒ Yes ☐ No ☐ N/A

Affected Product Information

Please list the quantity of affected product at your facility to be returned or write "none" if the affected product is no longer in your inventory or is no longer with your customers.

QUANTITY	PART NUMBER	LOT NUMBER
0	0	0

Form Completed By (Printed Name and Department):

Chris Paradise, Supply Chain

Date:

9. 24. 25

Email Address (Please print clearly for product return information):

christopher.paradise@mountsinai.org



COOK MEDICAL LLC
400 DANIELS WAY PO BOX 4195
BLOOMINGTON, IN 47402-4195
PHONE: 812.339.2235 TOLL FREE: 800-457-4500 FAX: 800.554.8335
WWW.COOKMEDICAL.COM

Returning this form and affected product(s)

- Upon completion of this form, please return it to the Cook Medical Customer Relations department via fax or email.
Fax: 812.339.7316 Email: FieldActionsNA@CookMedical.com
- Should you have any questions regarding the recall, please contact the Field Actions Department at Cook Medical Customer Relations. Please use FieldActionsNA@CookMedical.com, or call us toll free – 800.457.4500, press 5, then enter extension 151090. Monday through Friday, 8:00AM to 4:30PM, Eastern Time.
- Before returning affected product to Cook Medical, please obtain an RGA#.
 - Please ensure your email address is clear and complete. The RGA#, return address, and instructions will be emailed to the address provided.
 - You may also contact 800.457.4500 press option 5 and enter extension 151090 to obtain an RGA# and return address.
- To report adverse events, contact Cook Medical Customer Relations at 800.457.4500 extension 152146 or 812.339.2235 extension 152146

