

iVAS Elite Kits

Attn: Risk Manager, Materials Manager, OR Director


April 16, 2025

This notification is to inform you that Stryker Instruments is voluntarily recalling specific products of the iVAS Elite Kits.

Catalog Number	Product Description	GTIN	Lot Number
0808-010-000	10MM (10G) iVAS Elite Kit	07613327142686	See Appendix A
0808-010-800	10MM (10G) iVAS Elite Bipod Kit	07613327358308	See Appendix A
0808-015-000	15MM (10G) iVAS Elite Kit	07613327142785	See Appendix A
0808-015-800	15MM (10G) iVAS Elite Bipod Kit	07613327358278	See Appendix A
0808-020-000	20MM (10G) iVAS Elite Kit	07613327142693	See Appendix A
0808-020-800	iVAS 10G 20MM Multi-Kit Multi	07613327358292	See Appendix A
0808-110-000	10MM (11G) iVAS Elite Kit	07613327142709	See Appendix B
0808-110-800	10MM (11G) iVAS Elite Bipod Kit	07613327358315	See Appendix B
0808-115-000	15MM (11G) iVAS Elite Kit	07613327142716	See Appendix B
0808-115-800	15MM (11G) iVAS Elite Bipod Kit	07613327358339	See Appendix B

Product Description

The iVAS Elite kits are designed for use in vertebral augmentation procedures to provide controlled and contained cement delivery to the surgical site.

**Product Issue**

There is potential for the device to exhibit two failures: out of specification and/or leakage. For the out of specification devices, the narrow inner diameter of the cement tube prevents the flow of low-viscosity cement through the instrument. Devices impacted by the leak defect may introduce cement to unintended locations. Refer to above image for location of product issue.

Potential Risks

There is a potential risk of user annoyance or dissatisfaction due to the products failure to deliver cement. Both failure modes would be recognized intraoperatively. There may be a short delay if the product needs to be replaced. The product is optional and not required for completion of the procedure.

Actions Needed

1. Immediately review your inventory to locate and quarantine any affected products at your facility.
2. Return the enclosed Business Reply Form (BRF), even if the affected product is no longer in inventory, to the email Instruments.Recalls@stryker.com.
3. Upon receipt of the completed BRF, Stryker will provide a shipping label to return recalled product on-hand.
4. A replacement will be provided upon receipt of the recalled product.
5. Maintain awareness of this communication internally and inform Stryker if any of the subject devices have been distributed to other organizations. If so, provide contact details so Stryker can inform the recipients accordingly.

On behalf of Stryker, we thank you sincerely for your help and support in completing this action and regret any inconvenience that may be caused. We would like to reassure you that Stryker is committed to ensuring that only conforming devices, meeting our high internal quality standards and your expectations, remain on the market.

Sincerely,

Brooke Thompson | Recall Coordinator | (269) 290-0524 | Instruments.Recalls@stryker.com

For reporting any serious adverse events or product quality problems contact Stryker at 1-800-253-3210, ext. 5555 8:00 A.M. to 7:00 P.M. (Eastern Time), Monday – Friday or by email at inst.stryker.cs@stryker.com. In addition to contacting Stryker, healthcare professionals and consumers may report any potential quality problems or adverse reactions, or events associated with the use of this product may be reported to the FDA's MedWatch Adverse Event Reporting program either online, or by fax or phone. Online: www.fda.gov/Safety/MedWatch/HowToReport/default.htm Fax: (800) FDA-0178 Phone: (800) FDA-1088.

Stryker Corporation or its affiliates own, use or have applied for the following trademarks or service marks: Stryker All other trademarks are trademarks of their respective owners or holders.

Appendix A – 10G Elite Kits

A total of 173 lot numbers across six catalog numbers are listed below. Distribution dates in scope range from 12/11/2020 – 10/24/2024.

10MM (10G) iVAS Elite Kit (0808-010-000)			20MM (10G) iVAS Elite Kit (0808-020-000)	10MM (10G) iVAS Elite Bipod Kit (0808-010-800)		
20325012	20335012	20363012		20352012	21032012	21049012
21035012	21056012	21078012		21106012	21323012	22090012
21124012	21316012	21343012		22109012	22109022	22111012
22076022	22094012	22096012		22251012	22264012	22327012
22208012	22269012	22276012	22178012	22348012	23039012	23062012
22290012	22356012	23086012	22270012	23074012	23115012	23157012
23087012	23191012	23201012	23040012	23202012	23317012	23321012
23307012	23312012		23320012			

15MM (10G) iVAS Elite Kit (0808-015-000)					
20349012	20350012	20357012	20364012	21033012	21034012
21036012	21040012	21053012	21054012	21062012	21082012
21112012	21123012	21145012	21162012	21215012	21231012
21238012	21239012	21253012	21265012	21273012	21287012
21364012	22004012	22024012	22049012	22052012	22069012
22102012	22110012	22175012	22244012	22244022	22257012
22257022	22280012	22297012	22323012	22325012	22332012
22335012	23031012	23032012	23041012	23051012	23192012
23192022	23199012	23214012	23263012	23265012	23291012

15MM (10G) iVAS Elite Bipod Kit (0808-015-800)							
20330012	20344012	20345012	21027012	21048012	21109012	21110012	21111012
21118012	21148012	21167012	21207012	21228012	21229012	21230012	21252012
21264012	22005012	22010012	22013012	22018012	22091012	22137012	22138012
22145012	22173012	22173022	22204012	22204022	22207012	22207022	22252012
22252022	22266012	22281012	22291012	22291022	22301012	22307012	22307022
23010012	23024012	23073012	23093012	23095012	23110012	23135012	23150012
23174012	23186012	23212012	23268012	23277012	23313012	23317022	

iVAS 10G 20MM Multi-Kit Multi (0808-020-800)			
20325022	20351012	21028012	21299012
22139012	22196012	22312022	22349012
23030012	23094012	23108012	23151012
23151022	23198012	23269012	23314012

For reporting any serious adverse events or product quality problems contact Stryker at 1-800-253-3210, ext. 5555 8:00 A.M. to 7:00 P.M. (Eastern Time), Monday – Friday or by email at inst.stryker.cs@stryker.com. In addition to contacting Stryker, healthcare professionals and consumers may report any potential quality problems or adverse reactions, or events associated with the use of this product may be reported to the FDA's MedWatch Adverse Event Reporting program either online, or by fax or phone. Online: www.fda.gov/Safety/MedWatch/HowToReport/default.htm Fax: (800) FDA-0178 Phone: (800) FDA-1088.

Stryker Corporation or its affiliates own, use or have applied for the following trademarks or service marks: Stryker All other trademarks are trademarks of their respective owners or holders.

Appendix B – 11G Elite Kits

A total of 84 lot numbers across four catalog numbers are listed below. Distribution dates in scope range from 12/11/2020 – 10/24/2024.

10MM (11G) iVAS Elite Kit (0808-110-000)				
21055012	21120012	21120022	21212012	21216012
21257012	21291012	22144012	22144022	22245012
22283012	22306012	22321012	22334012	23185012

10MM (11G) iVAS Elite Biped Kit (0808-110-800)				
21029012	21042012	21173012	21237012	21300012
21342012	22035012	22092022	22140012	22195012
22273012	22364012	23089012	23153012	23181012

15MM (11G) iVAS Elite Kit (0808-115-000)				
20346012	21015012	21022012	21025012	21025022
21026012	21133012	21138012	21139012	21144012
21176012	21292012	21319012	21321012	21321022
22025012	22040012	22206012	22206022	22209012
22209022	22255012	22256012	22256022	22274012
22292012	23047012	23058012	23114012	23156012
23178012	23200012	23213012	23278012	23324012

15MM (11G) iVAS Elite Biped Kit (0808-115-800)				
20337012	20338012	21140012	21141012	21244012
21326012	22108012	22140022	22181012	22265012
22312012	22313012	22326012	22333012	22333022
22355012	23171012	23255012	23315012	

For reporting any serious adverse events or product quality problems contact Stryker at 1-800-253-3210, ext. 5555 8:00 A.M. to 7:00 P.M. (Eastern Time), Monday – Friday or by email at inst.stryker.cs@stryker.com. In addition to contacting Stryker, healthcare professionals and consumers may report any potential quality problems or adverse reactions, or events associated with the use of this product may be reported to the FDA's MedWatch Adverse Event Reporting program either online, or by fax or phone. Online: www.fda.gov/Safety/MedWatch/HowToReport/default.htm Fax: (800) FDA-0178 Phone: (800) FDA-1088.

Stryker Corporation or its affiliates own, use or have applied for the following trademarks or service marks: Stryker All other trademarks are trademarks of their respective owners or holders.

Business Reply Form

iVAS Elite Kits

April 16, 2025

Please select from the options below and complete this form. Email the completed form to Instruments.Recalls@stryker.com.

RESPONSE IS REQUIRED.

- ☐ No remaining affected products on-hand.
- ☐ I, the customer, choose to return the following product(s) for replacement

Catalog Number	Product Description	Lot number	Quantity on hand (EA)*
0808-010-000	10MM (10G) iVAS Elite Kit		
0808-010-800	10MM (10G) iVAS Elite Biped Kit		
0808-015-000	15MM (10G) iVAS Elite Kit		
0808-015-800	15MM (10G) iVAS Elite Biped Kit		
0808-020-000	20MM (10G) iVAS Elite Kit		
0808-020-800	iVAS 10G 20MM Multi-Kit Multi		
0808-110-000	10MM (11G) iVAS Elite Kit		
0808-110-800	10MM (11G) iVAS Elite Biped Kit		
0808-115-000	15MM (11G) iVAS Elite Kit		
0808-115-800	15MM (11G) iVAS Elite Biped Kit		

*If all devices have been used and none remain for return, please indicate 0 (zero) for quantity on hand. Please respond even if there is no record of receiving affected inventory.

Form completed by:

Facility Name			
Facility Address			
Printed Name		Title	
Email		Phone	
Signature		Date	

Note: Your signature indicates that you have received and understand the enclosed notification and that you have performed all actions requested.

If you have further distributed any affected product, please indicate to whom:

Facility Name		Contact Person	
Full Address			

For reporting any serious adverse events or product quality problems contact Stryker at 1-800-253-3210, ext. 5555 8:00 A.M. to 7:00 P.M. (Eastern Time), Monday – Friday or by email at inst.stryker.cs@stryker.com. In addition to contacting Stryker, healthcare professionals and consumers may report any potential quality problems or adverse reactions, or events associated with the use of this product may be reported to the FDA's MedWatch Adverse Event Reporting program either online, or by fax or phone. Online: www.fda.gov/Safety/MedWatch/HowToReport/default.htm Fax: (800) FDA-0178 Phone: (800) FDA-1088.

Stryker Corporation or its affiliates own, use or have applied for the following trademarks or service marks: Stryker All other trademarks are trademarks of their respective owners or holders.