


URGENT: MEDICAL DEVICE RECALL	
Description	Specific lots of Alcon Constellation Pak® with Damaged Tyvek®
Product Reference	Alcon Constellation Pak®
Market Action Identifier	2025.007

May 1, 2025

MOUNT SINAI HOSPITAL
3 E 101ST ST
NEW YORK, NY 10029-6528
Account # 100055760



Dear MOUNT SINAI HOSPITAL,

The purpose of this letter is to notify you that Alcon has initiated a Medical Device Recall for specific lots of its Alcon Constellation Pak®. Alcon is conducting this Medical Device Recall as there is potential for some trays within impacted lots to have damage in the Tyvek® cover.

The following affected product has been shipped to your facility:

Product Number	Product Description	Affected Lot(s)
8065753106	25+TTLPL VPK 20000CPM BEV VAL	1776DC

This event was identified internally, and to date Alcon has not received any reports of customer complaints or adverse events related to this issue.

Description of the Issue

There is the potential that the Tyvek® cover of some units within specific lots of Alcon Constellation Pak® were damaged during the manufacturing process. Please see photo below for an example. Due to the risk of Constellation Pak® sterility being compromised, Alcon is recalling potentially affected lots. The use of non-sterile surgical products may increase the risk of post-operative infection, which may require additional medical and/or surgical intervention.

Risk Management Department

MAY 06 2025

Received

Exhibit 1: Image of Constellation Pak® with Damaged Tyvek®



Actions to be taken by the Customer / User

We are asking that you locate and dispose of any affected lots of Alcon Constellation Pak® remaining in your inventory. To comply with this Medical Device Recall and request the replacement of any unused product, please take the following steps:

1. Review your inventory to determine if you have any unused affected product within your facility. See table on page 1 for affected Constellation Pak® lots shipped to your location.
2. Segregate and **dispose** of any unused affected product from your inventory.
3. Call Alcon Customer Service to arrange for replacement of your affected inventory of Alcon Constellation Pak®.
4. Post this notification letter near where affected products are stored to ensure awareness of the recall and continued shipment of corrected Constellation Pak®.
5. Respond to Alcon indicating your understanding of these instructions **even if you have zero (0) units remaining in inventory** by completing the online response form at <https://qrco.de/constellation-response> or by completing and returning this "Response Form" and returning to Alcon via email or fax.
6. Please forward this notification to all departments within your organization who may be in possession of this affected product; and any other organization to which this product may have been transferred.

Alcon Corrective Action

The Constellation Pak® is essential to support posterior segment (i.e., vitreoretinal) procedures using the Constellation® Vision System. To prevent immediate market shortage and the risk of cancelling emergency posterior segment surgeries, Alcon will reinspect the remaining inventory of affected Constellation Pak® lots. In the future, you may receive shipments of identified lots. The reinspected packs will be marked with a green indicator sticker to visually distinguish them from packs distributed prior to the recall initiation. The green sticker indicates the Constellation Pak® has been inspected and confirmed to be undamaged and is safe for use. An example of a Constellation Pak® with a green indicator sticker is provided in Exhibit 2 below.

Exhibit 2: Image of example green indicator sticker on Constellation Pak®



Contact for Further Questions about this Medical Device Recall

In the event you have experienced adverse events or product quality issues related to this communication, please contact Alcon at notifeye.alcon.com or by phone (1-800-757-9780, Option 5).

Adverse events or quality problems experienced with the use of this product may also be reported to the FDA MedWatch Adverse Event Reporting program. For more information, please visit fda.gov/Safety/MedWatch/HowToReport

Should you have any questions or concerns about this matter, please feel free to call Alcon Customer Service at 1-800-862-5266 or contact your Alcon Sales Representative.

Sincerely,

Heather Attra
SVP, Chief Quality & Regulatory Affairs Officer



Respond Electronically using QR Code or link below

<https://qrco.de/constellation-response>

As an alternative, respond via email/fax per the steps below

RESPONSE FORM

**Specific lots of Alcon Constellation Pak®
with Damaged Tyvek®
MA# 2025.007**

MOUNT SINAI HOSPITAL
3 E 101ST ST
NEW YORK, NY 10029-6528
Account# 100055760

To comply with this Medical Device Recall, please take the following steps:

1. Review your inventory to determine if you have any unused affected product within your facility.

Product Number	Product Description	Affected Lot(s)	Units Disposed
8065753106	25+TTLPL VPK 20000CPM BEV VAL	1776DC	

2. Segregate and **dispose** of any unused affected product from your inventory.
3. Call Alcon Customer Service to arrange for replacement of your affected inventory of Alcon Constellation Pak®.
4. Post this notification letter near where affected products are stored to ensure awareness of the recall and continued shipment of corrected Constellation Pak®.
5. Respond to Alcon indicating your understanding of these instructions **even if you have zero (0) units remaining in inventory** by completing the online response form at <https://qrco.de/constellation-response> or by completing and returning this "Response Form" and returning to Alcon via email or fax.
6. Please forward this notification to all departments within your organization who may be in possession of this affected product; and any other organization to which this product may have been transferred.

Please return this Response Form via fax or email to Alcon:

Fax: 817-302-4337

Email: Market.Actions@Alcon.com

Your signature below attests that you have read and understood this notification.

Signature:

Date:

Printed Name:

Title: