

MAY 06 2025

Received



GE HealthCare

**IMPORTANT PRODUCT
INFORMATION**

April 28, 2025

GE HealthCare Ref. # 30115

To: Healthcare Administrator
Network Administrator
Director of IT Department

RE: Potential security vulnerability for MUSE and DICOM Gateway Pro systems

Dear Customer,

GE HealthCare has provided security instructions with MUSE v7, MUSE v8, MUSE v9, MUSE NX R1, MUSE NX R2, DICOM Gateway Pro v9, DICOM Gateway Pro NX R1, and DICOM Gateway Pro NX R2 systems (see table below). **If these security instructions are not followed, it may result in unauthorized access from your local network to these devices.**

NOTE: This notification is **NOT** related to a product safety issue.

Product	Device Product Code
MUSE v7	Product Code = SBP
MUSE v8	Product Code = SHF (MUSE v8 core sw) Product Code = SHG (MUSE v7-to-v8 upgrade) Product Code = SHJ (MUSE v8 System Application DVD) Product Code = SKW (MUSE v8 core sw w/CCG)
MUSE v9	Product Code = SNK (MUSE v9 ATO) Product Code = SN6 (MUSE v9 kits)
DICOM Gateway Pro v9	Product Code = SN7
MUSE NX R1	Product Code = SSR (production) Product Code = SSS (test)
DICOM Gateway Pro NX R1	Product Code = STJ
MUSE NX R2	Product Code = SWJ (production) Product Code = SWK (test)
DICOM Gateway Pro NX R2	Product Code = SWN

Please follow the Security Compliance checklist in the Privacy & Security Guide that accompanied these systems and **implement good security practices** by ensuring that your systems are **not using a well-known password**. If passwords are well-known, update the credentials according to your organization's password policies. Please feel free to contact GE HealthCare Service or your local Service Representative if you need assistance.

Please ensure all potential users in your facility are made aware of this Important Product notification and the recommended actions. Please retain a copy of this letter for your records.

Please be assured that maintaining a high level of safety and quality is our highest priority. If you have any questions or concerns regarding this notification, please contact GE HealthCare Service at 1-800-437-1171 or your local Service Representative.

Sincerely,

A handwritten signature in black ink that reads "Amit Tandon". The signature is written in a cursive style and is underlined with a single horizontal stroke.

Amit Tandon
General Manager, Diagnostic Cardiology
GE HealthCare