

MAY 06 2025

Received

March 25, 2025

Instructions for use (IFU) for da Vinci X/Xi/SP Systems: Handheld Camera

Dear Intuitive Customer,

Intuitive is updating the Instructions for Use (IFUs) to support new products and features for da Vinci X, Xi, and SP systems. The Handheld Camera is a 2D camera head that enables the use of third-party laparoscopes with an eyepiece cap compliant with ISO/TS 18339 with the da Vinci X, Xi, and SP Vision Cart. It is compatible with 5 to 10-mm rigid laparoscopes and enables optical first entry with the da Vinci 8-mm optical trocar.

To provide continuous product improvements to our customers, Intuitive has updated the da Vinci Handheld Camera to improve reliability, improve performance with select third-party laparoscopes, and enable compatibility with the da Vinci SP system.

The IFUs listed in the table below are now available on the Intuitive IFU website (<https://manuals.intuitivesurgical.com/>) to provide instantaneous access and reduce environmental impact.

There are two options to obtain updated IFUs:

- 1- Access electronic copies of the IFUs online:
<https://manuals.intuitivesurgical.com/>
- 2- Call customer service to request a mailed copy

Please take the following actions to maintain your IFUs:

da Vinci X/Xi IFUs

Part Number	Description	Actions to perform
554074-02	da Vinci Handheld Camera User Manual Addendum	Add to the user manual binder
554075-02	da Vinci Handheld Camera Reprocessing Instructions	Add to the reprocessing instructions appendices

For any questions regarding this announcement, please contact your local Intuitive sales representative or Intuitive Customer Service at: +1.800.876.1310 or email customerservice@intusurg.com

Sincerely,
Intuitive