

NOV 05 2025

Received

stryker

URGENT

Medical Device Recall

NICO Myriad™ Handpiece | NICO Myriad™ Illumination Pack | NICO BrainPath

Attn: OR Director, Materials Manager, Risk Manager | October 22, 2025

Products Affected

Catalog Number	Product Description	GTIN	Lot Numbers
Refer to page 2 - 8			

Product Description

The **NICO Myriad System** consists of a sterile, single patient use handpiece designed to be utilized only with an approved NICO console and foot pedal. The **NICO Myriad Illumination Pack** consists of Handpiece Sleeve(s) and an Illumination Fiber. The Pack is provided sterile and designed to be utilized only with the NICO Myriad Handpiece and a NICO light source (LX or SPECTRA). The **NICO BrainPath** consists of multiple-sized reusable and re-sterilizable obturators with coordinating single patient use disposable sheaths. The obturator and sheath are assembled in the operating room immediately prior to use. After placement, the obturator is removed leaving behind the sheath which provides a 13.5 mm or 11 mm surgical corridor.

Product Issue

It has been identified that the tape used to secure NICO components within the packaging contains latex, despite the products being labelled as latex-free. There are 0 complaints or adverse events reported for this issue.

Potential Risks

Individuals without known latex allergies are not at risk. However latex exposure may pose a health risk to individuals with latex allergies. If a latex sensitized individual comes into direct and sustained contact with the tape, there is a potential for an allergic reaction which may include localized skin irritation. More serious allergic reactions such as hives, swelling, respiratory distress, or anaphylaxis could occur. These reactions require immediate medical intervention and are considered rare in this context due to limited exposure and the use of gloves in surgical settings.

Given the nature of the exposure and the typical use environment (e.g., operating room), the overall risk of severe harm is considered low, but not negligible. Users with known latex allergies should avoid direct contact with the affected product.

Actions Needed

1. Immediately review this notification and distribute this letter internally to all relevant users of the NICO product families.
2. To ensure internal awareness, please mark the product with "CONTAINS LATEX" and place a copy of this notice at the product's storage location before use.
3. If you have received this notice directly from Stryker, please complete and return the enclosed Business Reply Form (BRF) to Instruments.recalls@stryker.com, even if the affected product is no longer in your inventory. Submitting the BRF prevents duplicate notices from being sent.
4. Product returns are not required. However, if you would like to receive a replacement, you may choose one of the following options:
 - o **Option 1:** Destroy the affected products onsite with your Stryker representative and send a completed Certificate of Destruction (COD, pg. 9) to Instruments.recalls@stryker.com. Replacements will be issued by Stryker once both the BRF and COD are received.
 - o **Option 2:** Return the products to Stryker using a provided shipping label. A shipping label will be provided once you submit a completed BRF and replacements will be issued once the product is received.
5. If you have distributed any of the affected devices to other organizations, please inform Stryker and provide their contact details so we can notify them directly.

On behalf of Stryker, we thank you sincerely for your help and support in completing this action and regret any inconvenience that may be caused. We would like to reassure you that Stryker is committed to ensuring that only conforming devices, meeting our high internal quality standards and your expectations, remain on the market. If you have any questions about this matter, contact Stryker Instruments by email at Instruments.recalls@stryker.com. For reporting any serious adverse events or product quality problems contact Stryker at 1-800-253-3210, ext. 5555 8:00 A.M. to 7:00 P.M. (Eastern Time), Monday - Friday or by email at inst.stryker.cs@stryker.com.

In addition to contacting Stryker, healthcare professionals and consumers may report any potential quality problems or adverse reactions, or events associated with the use of this product may be reported to the FDA's MedWatch Adverse Event Reporting program either online, or by fax or phone. Online: www.fda.gov/Safety/MedWatch/HowToReport/default.htm Fax: (800) FDA-0178 Phone: (800) FDA-1088.

Sincerely,

Brooke Thompson | Recall Coordinator | (269) 290-0524 | Instruments.Recalls@stryker.com

Stryker Corporation or its affiliates own, use or have applied for the following trademarks or service marks: Stryker, Myriad. All other trademarks are trademarks of their respective owners or holders.

Products Affected

Material #	Product Description						GTIN			
NN-1000	11g Illumination Sleeve/Fiber						00816744021009			
	Lot numbers									
	37AC5163	37AE9811	37AJ1438	37BJ0513	37CB3410	37CB3416	37CD6003	37CE3030	37EB5720	
	37FA4260	37GD2216	37GJ1603	37IC7023	37IJ0208	37JJ0122	37KI1300	37LA0158	37LI0320	37LJ1522
Material #	Product Description						GTIN			
NN-1001	13g Illumination Sleeve/Fiber						00816744021016			
	Lot numbers									
	37AC5177	37AE9942	37AE9955	37AJ1439	37BJ0516	37CB3380	37CB3419	37CD6005	37CJ0204	
	37DE4864	37EB5721	37FA4261	37GD2219	37GJ1607	37IC7027	37IJ0209	37JD5753	37JJ0123	
	37KI1301	37LA0159	37LI0322	37LJ1525	37EJ2020					
Material #	Product Description						GTIN			
NN-2000	11x50 BP, 13x13 MHP/PB - MIPS Pack						00816744022006			
	Lot numbers									
	37AI1069	37AJ2002	37CH0527	37CI0836	37CI1929	37CJ2202	37DC0092	37ED9009	37EI2165	
	37EJ2015	37FA3139	37FD0972	37FH0128	37FH2524	37GH2417	37GJ1614	37HA6437	37IH0510	
	37JH0925	37JI2169	37KD7602	37KH1319	37KI1901	37KJ0700	37LB4003	37LC1480		
Material #	Product Description						GTIN			
NN-2001	11x60 BP, 13x13 MHP/PB - MIPS Pack						00816744022013			
	Lot numbers									
	37AB1573	37AC5618	37AJ2003	37CH0530	37CI0837	37CI1927	37DA9502	37DC0093	37EI2166	
	37EJ2017	37FA3529	37FD0981	37FH1528	37GH2513	37HB8964	37IH0511	37JH0932	37JI2171	
	37KH1567	37KI1902	37KJ0701	37LB4002						
Material #	Product Description						GTIN			
NN-2002	11x75 BP, 13x13 MHP/PB - MIPS Pack						00816744022020			
	Lot numbers									
	37AC5619	37AJ2004	37CB4046	37CH0532	37CI0838	37CI1926	37CJ2203	37DC0094	37EI2167	
	37EJ2018	37FH1118	37GH3035	37IA7267	37IH0512	37JH0938	37JI2172	37KH1582	37KI1903	
	37KJ0702	37LB3889	37LC1487							
Material #	Product Description						GTIN			
NN-2003	13.5x50 BP, 13x13 MHP/PB - MIPS Pack						00816744022037			
	Lot numbers									
	37AA1900	37AB1578	37AC5620	37AI1072	37AJ2006	37BI1612	37CA0223	37CB4060	37CE2876	
	37CH0529	37CI0839	37CI1925	37CJ2204	37DA9503	37DC0095	37ED9117	37EH2123	37EI2169	
	37EJ2020	37FA3632	37FD0995	37FH2634	37FH2929	37GC4807	37GH0201	37GJ1616	37HA6438	
	37HB8966	37HH2932	37HI1649	37IA7272	37ID4854	37IH1334	37IH2709	37IJ0314	37JH0424	
	37JH0942	37JI2173	37KD7603	37KH1566	37KI1904	37KJ0703	37LB4004	37LC1489	37LD8383	37LH0606
Material #	Product Description						GTIN			
NN-2004	13.5x60 BP, 13x13 MHP/PB - MIPS Pack						00816744022044			
	Lot numbers									
	37AA1901	37AB1580	37AC5623	37AI1525	37AI2810	37AJ2007	37BI1614	37CA0224	37CB4065	
	37CH0533	37CI0840	37CI1923	37CJ2205	37DA9504	37DC0112	37DJ2212	37ED9118	37EH2208	
	37EI2182	37FD0996	37FH2636	37FH2920	37GC4810	37GH0200	37GJ1617	37HA6439	37HB8968	
	37HI1651	37IA7273	37ID4856	37IH0400	37IH1349	37IH2715	37JH0426	37JH1034	37JI2174	
	37KD7608	37KH1570	37KI1905	37KJ0704	37LB4005	37LC1493	37LD8389	37LH0621		
Material #	Product Description						GTIN			
NN-2005	13.5x75 BP, 13x13 MHP/PB - MIPS Pack						00816744022051			
	Lot numbers									
	37AA1902	37AB1581	37AC5624	37AE0073	37AE1259	37AI1615	37AI2908	37AJ2008	37BI1617	
	37CA0225	37CB4066	37CH0534	37CI0841	37CI1922	37CJ2206	37DA9505	37DC0121	37DJ2213	
	37ED9119	37EE5692	37EH2210	37EI2183	37FA3634	37FD0997	37FH1803	37FH2923	37GC4812	
	37GH0203	37HA6440	37HB8974	37HI1652	37IA7274	37ID4857	37IH1355	37IH2710	37IJ0315	
	37JH1037	37JI2175	37KD7609	37KH1325	37KI1906	37KJ0705	37LB4006	37LC1504	37LH0624	37LH2736
Material #	Product Description						GTIN			
NN-2006	13.5x95 BP, 13x13 MHP/PB - MIPS Pack						00816744022068			
	Lot numbers									
	37AA1903	37AC5628	37BI1618	37CA0229	37CB4072	37CH0535	37CI0842	37CI1921	37CJ2207	
	37DA9506	37DC0122	37EI2185	37EJ2021	37FD0998	37FH1804	37FH2926	37GC4825	37GH0208	
	37GJ1618	37HA6443	37HB8985	37HI1646	37IH1356	37IJ0318	37JH1116	37JI2176	37KD7610	
	37KI1907	37LB4007	37LC1508	37LH0626						

Material #	Product Description						GTIN			
NN-2024	13.5x50 BP, 11x13 MHP/PB - MIPS Pack						00816744022242			
	Lot numbers									
	37AA1904	37DA2762	37ED9122	37EI2186	37GA4303	37GJ0602	37HD2739			
Material #	Product Description						GTIN			
NN-2025	13.5x60 BP, 11x13 MHP/PB - MIPS Pack						00816744022259			
	Lot numbers									
	37AA1905	37AC5350	37AE1265	37DA2764	37EI2191	37GA4304	37GJ0603	37HB8986	37ID4859	
Material #	Product Description						GTIN			
NN-2026	13.5x75 BP, 11x13 MHP/PB - MIPS Pack						00816744022266			
	Lot numbers									
	37AA1906	37CD5330	37DA2768	37ED9124	37EI2197	37GA4305	37GJ0604			
Material #	Product Description						GTIN			
NN-2027	13.5x95 BP, 11x13 MHP/PB - MIPS Pack						00816744022273			
	Lot numbers									
	37AC5351	37EI2198	37ID4860							
Material #	Product Description						GTIN			
NN-2500	11x50 BP, 13x13 MHP/PB - MIPS Pack						00816744022501			
	Lot numbers									
	37HB8995	37HB8995	37KA9437	37KA9437						
Material #	Product Description						GTIN			
NN-2501	11x60 BP, 13x13 MHP/PB - MIPS Pack						00816744022518			
	Lot numbers									
	37HB8996	37JC8651	37KA9382	37KA9438	37LC1512					
Material #	Product Description						GTIN			
NN-2502	11x75 BP, 13x13 MHP/PB - MIPS Pack						00816744022525			
	Lot numbers									
	37HB8997	37II1122	37JC8652	37KA9383	37KA9439	37LC1623				
Material #	Product Description						GTIN			
NN-2503	13.5x50 BP, 13x13 MHP/PB - MIPS Pack						00816744022532			
	Lot numbers									
	37CD5599	37CJ2208	37DI2635	37EC1887	37FD1003	37IA7276	37II1124	37JC8277	37LA0160	37LC1624
Material #	Product Description						GTIN			
NN-2504	13.5x60 BP, 13x13 MHP/PB - MIPS Pack						00816744022549			
	Lot numbers									
	37CD5600	37CJ2209	37DI2636	37ED9126	37FD1005	37GA4555	37HA5487	37HC5380	37HC5573	
	37IA7291	37II1125	37JB1444	37KB3710	37KC0369	37LC1625				
Material #	Product Description						GTIN			
NN-2505	13.5x75 BP, 13x13 MHP/PB - MIPS Pack						00816744022556			
	Lot numbers									
	37CD5604	37CJ2210	37DI2637	37EC1889	37FD1008	37GA4556	37HA5488	37IA7292	37II1126	37JB1449
Material #	Product Description						GTIN			
NN-2506	13.5x95 BP, 13x13 MHP/PB - MIPS Pack						00816744022563			
	Lot numbers									
	37DI2638	37HB8998	37HC5383	37IA7293	37II1127	37KA9952	37KC0372			
Material #	Product Description						GTIN			
NN-3000	11x50 BP, 13x13 MHP/PB - MIPS Illum Pack						00816744023003			
	Lot numbers									
	37EC1890	37AA0405	37BA1865	37DA2218	37HA6021	37BB2278	37AJ2002	37CJ2202	37IJ0209	
Material #	Product Description						GTIN			
NN-3001	11x60 BP, 13x13 MHP/PB - MIPS Illum Pack						00816744023010			
	Lot numbers									
	37LC1628	37KD8021	37AA0407	37BA1866	37DA2219	37HA6023	37BB2279	37AJ2003	37IJ0209	
	37AJ1439	37KI1902	37EJ2017	37II0123	37GI1607					

Material #	Product Description					GTIN			
NN-3002	11x75 BP, 13x13 MHP/PB - MIPS Illum Pack					00816744023027			
	Lot numbers								
	37LD8610	37LD8778	37AA0408	37DA2220	37HA6024	37BB2280	37CJ2203	37IJ0209	37BJ0516
	37KI1903	37AJ2004	37GJ1607	37AJ1439	37JJ0123				
Material #	Product Description					GTIN			
NN-3003	13.5x50 BP, 13x13 MHP/PB-MIPS Illum Pack					00816744023034			
	Lot numbers								
	37FD1010	37ID4861	37CE2573	37AA0409	37DA2221	37HA6025	37BB2283	37EJ2020	37GJ1607
	37IJ0209	37DJ1036	37FJ2203	37CJ2204	37AJ2006	37BJ0516	37CJ0204		
Material #	Product Description					GTIN			
NN-3004	13.5x60 BP, 13x13 MHP/PB-MIPS Illum Pack					00816744023041			
	Lot numbers								
	37AA0410	37AJ2007	37BB2305	37BJ0516	37CE2576	37DA2222	37DJ2212	37FD1014	37GJ1607
	37HA6027	37ID4862	37IJ0209	37KC0073	37LC1629				
Material #	Product Description					GTIN			
NN-3005	13.5x75 BP, 13x13 MHP/PB-MIPS Illum Pack					00816744023065			
	Lot numbers								
	37AA0411	37AC4930	37AJ2008	37BB2307	37BJ0516	37CE2579	37CJ2206	37DA2224	37DJ2213
	37EC1891	37GJ1607	37HA6028	37IJ0209					
Material #	Product Description					GTIN			
NN-3006	13.5x95 BP, 13x13 MHP/PB-MIPS Illum Pack					00816744023065			
	Lot numbers								
	37AA0412	37AC4931	37BB2308	37CJ2207	37DA2226	37EJ2021	37FD1016	37GJ1607	37HA6029
	37IJ0209	37JJ0123	37KD8023						
Material #	Product Description					GTIN			
NN-3007	13.5x50 BP, 11x13 MHP/PB-MIPS Illum Pack					00816744023072			
	Lot numbers								
	37AA0413	37AB1800	37BB2309	37EB5748	37EB6231	37EC1892	37FD1019		
Material #	Product Description					GTIN			
NN-3008	13.5x60 BP, 11x13 MHP/PB-MIPS Illum Pack					00816744023089			
	Lot numbers								
	37AA0414	37AC4928	37BB2310	37EB5757	37EB6232	37EC1895	37GA5424	37LA0663	
Material #	Product Description					GTIN			
NN-3009	13.5x75 BP, 11x13 MHP/PB-MIPS Illum Pack					00816744023096			
	Lot numbers								
	37AA0415	37EB5759	37EB6235	37GA5433	37LA0670				
Material #	Product Description					GTIN			
NN-3010	13.5x95 BP, 11x13 MHP/PB-MIPS Illum Pack					00816744023102			
	Lot numbers								
	37AA0416	37CB3332	37DB4837	37EB5761	37EB6236	37EC1899			
Material #	Product Description					GTIN			
NN-8000	15x10 MHP					00816744028008			
	Lot numbers								
	37BH1415	37CB3768	37DC0517	37EB6074	37GA4892	37GB7928	37IJ2310	37JJ0604	37KB3653
	37KI1915	37LA0162							
Material #	Product Description					GTIN			
NN-8001	15x 13 MHP					00816744028015			
	Lot numbers								
	37BC7385	37DJ1035	37EB5920	37EB6075	37ED9464	37FC3028	37GD1188	37HB0020	37HD3586
	37ID3702	37IA8380	37ID5540	37IJ2312	37KH3023	37KI1916	37LC1191	37LI1727	

Material #	Product Description					GTIN				
NN-8002	15x 25 MHP - MINOP					00816744028022				
	Lot numbers									
	37AC5372	37AC5631	37AE1273	37AH1830	37BA1102	37BB2319	37BC7386	37BD3731	37BJ1508	
	37CB4081	37CD5349	37CD5637	37DA1503	37DA9507	37DC0519	37DH3011	37DI2303	37DJ1036	
	37EC1900	37ED8566	37ED9571	37EE5791	37FC3030	37FI2611	37FJ2203	37GA4306	37GA5370	
	37GC4834	37GI1763	37GJ1619	37HB9007	37HD2741	37HH0317	37HI1617	37HJ1814	37IA7295	
	37ID4143	37ID5031	37II1827	37IJ0319	37JB1412	37JD5541	37JI2101	37KD7643	37KG1725	
	37KH0104	37KH3004	37KJ0706	37LA0163	37LA0680	37LC1192	37LD8928	37LG1840	37LI1731	
	37LJ0817									

Material #	Product Description					GTIN				
NN-8003	13x10 MHP					00816744028039				
	Lot numbers									
	37AJ2009	37BC7387	37BH1419	37BJ1509	37CD5353	37DC0522	37EB5596	37EB5984	37ED8570	
	37EI2141	37FA3637	37FB6706	37FC3037	37FD0078	37GJ1620	37HB0021	37HH0314	37IB0874	
	37ID5058	37IG2927	37II1824	37JA8394	37JG1103	37JH0115	37JI2102	37KB3657	37KD7651	
	37KI1917	37LC1674	37LH1805	37LI1733						

Material #	Product Description					GTIN				
NN-8004	13x13 MHP					00816744028046				
	Lot numbers									
	37AA0522	37AC5387	37AC5632	37AE9800	37AH1052	37AH1741	37AI2352	37AJ2011	37BA0230	
	37BB2320	37BC7388	37BD3754	37BH1422	37BI1926	37CA0227	37CB4082	37CD5553	37CD5638	
	37CH0723	37CJ2211	37DA2200	37DA9508	37DB4672	37DC0524	37DC9362	37DH0251	37DH3014	
	37DI2304	37DJ1037	37EA2247	37EB5985	37EC1901	37ED8572	37EI2145	37EJ2022	37FA3639	
	37FD0084	37FH0606	37FI2612	37FJ2205	37GA4308	37GA5371	37GC4336	37GC4835	37GD1193	
	37GI1764	37GJ1622	37HA6682	37HB0022	37HB9008	37HD2742	37HH0313	37HH2305	37HI1618	
	37HJ1815	37IA7297	37IB0875	37IB0876	37IC6867	37ID4150	37ID5059	37II1833	37JA8395	
	37JC8278	37JD5597	37JG0201	37JG1107	37JH0116	37JI2103	37KB3658	37KB3731	37KC9435	
	37KD7675	37KH0161	37KH3005	37KJ0506	37LA0164	37LA0703	37LC1194	37LC1675	37LD8936	
	37LH1809									

Material #	Product Description						GTIN		
NN-8005	15x25 MHP GAAB						00816744028053		
	Lot numbers								
	37BI1922	37FI2613	37HA6683	37HI1100	37IB0877	37II1835	37LI0818		

Material #	Product Description					GTIN				
NN-8006	11X10 MHP					00816744028060				
	Lot numbers									
	37AC5388	37DB4673	37DC9363	37EA2818	37EC1903	37ED8458	37EI2146	37FD9828	37FI2614	
	37GB7926	37HI1621	37IC6868	37ID4154	37IJ0320	37JB1404	37JG0203	37KI1918	37LA0166	
	37LD8946									

Material #	Product Description					GTIN				
NN-8007	11x13 MHP					00816744028077				
	Lot numbers									
	37AB1587	37AC5634	37AH1825	37AI2353	37AJ2012	37BA0231	37BB2321	37BC7389	37BD3756	
	37BI1924	37BJ1828	37CA1600	37CD5554	37CD5641	37CI1910	37CJ2212	37DA9509	37DB4676	
	37DC9509	37DH0253	37DH3009	37DI2305	37DJ1038	37EA0600	37EA2723	37EA2820	37EB5994	
	37EC1904	37EI2149	37EJ2024	37FA3651	37FH0608	37FI2615	37FJ2206	37GA4309	37GA5372	
	37GC4839	37GD1198	37GI1765	37HA6691	37HD2750	37HH0318	37HH2308	37IA7299	37IB0878	
	37IC6871	37ID5090	37ID5092	37JA8457	37JC8279	37JD5598	37JG0205	37JH0117	37JI2104	
	37JJ0513	37KB3664	37KC9467	37KD7677	37KH0165	37KH3009	37KI1919	37LD8949	37LG1837	
	37LH1810	37LI1800	37LJ0819							

Material #	Product Description						GTIN			
NN-8008	19x28 MHP OI						00816744028084			
	Lot numbers									
	37AB1597	37AJ3013	37BJ1829	37CA0230	37CJ2213	37DJ1039	37FI2616	37GA4310	37GC4841	37HB9009
	37IA7302	37IC6875	37JB1405	37JG0207	37JJ0605	37KG1418	37KG1510	37LI1804	75HG2942	75KI2015

Material #	Product Description					GTIN			
NN-8010	17x31.5 MHP DECQ					00816744028107			
	Lot numbers								
	37BJ1834	37DH0303	37DJ1040	37FA3669	37IJ1015	37LA0704	37LC1195		

Material #	Product Description						GTIN			
NN-8015	NICO Fluid System (NFS)						00816744028152			
	Lot numbers									
	37AE9793	37BD3764	37DA1557	37DA1558	37DC9510	37EA2474	37EB5685	37FC3045	37GA4316	
	37GD1216	37ID5098	37JA8398	37JB1410	37JB1411	37KB3726				
Material #	Product Description						GTIN			
NN-8017	11x13 MHP/PB						00816744028176			
	Lot numbers									
	37AC5635	37AE0103	37AI2354	37BD3759	37BH1426	37BI1930	37CB4085	37CD5642	37CH0720	
	37CI1911	37DA2001	37DC0534	37DH0306	37DH3004	37DI2306	37DJ1041	37EB5686	37EB5997	
	37ED8576	37EE5825	37EI2156	37EJ2026	37FA2828	37FD0144	37FH0611	37FI2617	37FJ2207	
	37GC4847	37GI1767	37HA6692	37HB0024	37HH0321	37HH2317	37HJ1816	37ID4155	37II1841	
	37JG0208	37JG1105	37JH0122	37JJ0608	37KB3666	37KG2107	37KH3020	37LA0705	37LC1200	
	37LD8950	37LG1845	37LH1811	37LI1805						
Material #	Product Description						GTIN			
NN-8018	13X13 MHP/PB						00816744028183			
	Lot numbers									
	37AA0524	37BC7397	37BE1491	37EB6000	37HA6697	37HD2753	37HH0320	37IB0880	37LA0170	37LD8952
Material #	Product Description						GTIN			
NN-8019	19x28 MHP Little LOTTA						00816744028190			
	Lot numbers									
	37AB1601	37AC5636	37AE1285	37AE9804	37AH1036	37BA0232	37BD3760	37BJ1830	37CD5557	
	37DA2202	37DB4677	37DC0535	37DC9514	37DI2633	37DJ1042	37ED8577	37EI2140	37FC3055	
	37FH0602	37FI2609	37GC4851	37GD1223	37GI1768	37HB9026	37HD2754	37HJ1817	37ID5100	
	37JB1422	37JI0713	37KB3674	37KC9470	37KD7723	37KG1416	37KH0163	37LC1210	37LC1681	
	37LJ1527	75HG2936	75HG2937	75HG2940	75HG2941					
Material #	Product Description						GTIN			
NN-8020	19x21.5 MHP PaediScope						00816744028206			
	Lot numbers									
	37BC7399	37CD5558	37EC1905	37EI2139	37EJ2027	37FC3092	37GI1769	37HB0025	37HH0322	
	37HI0130	37IA7303	37LI1807							
Material #	Product Description						GTIN			
NN-8021	15x26.5 MHP LOTTA						00816744028213			
	Lot numbers									
	37AA0526	37AC5637	37AD2692	37BA0400	37BB2322	37BH1429	37BJ1831	37CA0232	37CB4087	
	37CD5559	37CD5644	37CE2596	37CH0715	37CI1909	37CI1912	37DA2239	37DA9512	37DC0536	
	37DD8119	37DH3012	37DI2623	37DJ1043	37EB6002	37EC0830	37EC1908	37ED8582	37EE6617	
	37EI2160	37EJ2028	37FB6637	37FC3093	37FD0150	37FH0613	37FI2618	37GA4345	37GC4852	
	37GD1447	37GI1772	37HB0028	37HB9027	37HD2765	37HI1623	37HJ1101	37HJ1818	37ID4158	
	37ID5104	37II1843	37IJ1005	37JD5600	37JD5782	37JH0125	37JI2165	37JJ0609	37KB3675	
	37KC9471	37KG1412	37KH0162	37KH3022	37LA0707	37LC1683	37LD8956	37LI1809	37LJ1528	
Material #	Product Description						GTIN			
NN-8024	BrainPath Shephard's Hook - Greenberg						00816744028244			
	Lot numbers									
	37AI2944	37AJ2013	37BD3753	37BD3778	37CA1605	37CB3738	37CH0636	37CI1919	37EH0405	
	37FC3059	37FD0154	37GA4512	37GH1715	37GI0110	37GJ1623	37HB0029	37HH2313	37HI1536	
	37HI1631	37IA7694	37ID4162	37IG2900	37JB1413	37JH0131	37JH2903	37KG1406	37KH2953	
	37LA0708	37LG1834								
Material #	Product Description						GTIN			
NN-8025	BrainPath Shephard's Hook - Sugita						00816744028251			
	Lot numbers									
	37BB2326	37BJ1836	37CH0640	37DC0545	37HD2766	37KB3677	37KJ1139			
Material #	Product Description						GTIN			
NN-8026	BrainPath Shephard's Hook - Budde						00816744028268			
	Lot numbers									
	37BJ1837	37CC7852	37CH0643	37DJ1056						
Material #	Product Description						GTIN			
NN-8042	BrainPath Shephard's Hook - Greenberg LK						07613327663280			
	Lot numbers									
	37CC9220	37ED9408	37FJ2202	37II2308	37LD8666					

Material #	Product Description						GTIN			
NN-8500	15X10 MHP						00816744025007			
	Lot numbers									
	37AC5392	37DC9515	37DJ1044	37HI1634	37IJ1007	37KI1920				
Material #	Product Description						GTIN			
NN-8501	15x13 MHP						00816744025014			
	Lot numbers									
	37CC7849	37DC9516	37DI2627	37DJ1045	37FD0155	37HI1638	37JA8443	37KI1921		
Material #	Product Description						GTIN			
NN-8502	15x25 MINOP MHP						00816744025021			
	Lot numbers									
	37BA1106	37CD5646	37DI2628	37DJ1046	37GD1474	37HC4974	37HI1639	37HJ1102	37IA7304	
	37JB1455	37KB3678	37KB3706	37KD6865	37LA0710	37LI1810				

Material #	Product Description					GTIN				
NN-8503	13x10 MHP					00816744025038				
	Lot numbers									
	37AC5393	37BE1488	37BE1638	37CJ2214	37ED9004	37FC3095	37GD2039	37GI1773	37IJ1008	
	37KD6867	37LI1811								
Material #	Product Description					GTIN				
NN-8504	13X13 MHP					00816744025045				
	Lot numbers									
	37AC5646	37AH0600	37AH1829	37AJ2015	37BE1489	37CA0233	37CB4111	37CH0724	37DA9514	
	37EB6004	37EE6465	37EJ2031	37GA5374	37GE8331	37HH0315	37HI1641	37IA7306	37IH2628	
	37KD6878	37LA0729	37LH1812	75AH0216						
Material #	Product Description					GTIN				
NN-8505	15x25 GAAB MHP					00816744025052				
	Lot numbers									
	37AC4878	37DI2629	37DJ1048	37HJ1103	37IA7312					
Material #	Product Description					GTIN				
NN-8506	11X10 MHP					00816744025069				
	Lot numbers									
	37CI1836	37DJ1049	37GA4513	37IJ1009	37JH2433	37KB3681	37LH1813	37LI1814		
Material #	Product Description					GTIN				
NN-8507	11x13 MHP					00816744025076				
	Lot numbers									
	37BA0401	37BB2328	37BJ1838	37CI1916	37DJ1051	37GI1778	37IJ1010	37JB1456	37JG0214	37LA0174
Material #	Product Description					GTIN				
NN-8508	19x28 OI MHP					00816744025083				
	Lot numbers									
	37CI1917	37CJ2215	37HB9031	37IC6878	37JJ0610					
Material #	Product Description					GTIN				
NN-8510	17x31.5 DECQ MHP					00816744025106				
	Lot numbers									
	37BJ1835	37CB4112	37FC3097	75KI2225						
Material #	Product Description					GTIN				
NN-8517	13x13 MHP/PB					00816744025175				
	Lot numbers									
	37AI2946	37CB4118	37DA2701	37EH2906	37GA5375	37GC4337	37HH2903	37JB1457	37KH2955	
Material #	Product Description					GTIN				
NN-8518	11X13 MHP/PB					00816744025182				
	Lot numbers									
	37GA4637	37HB9032	37KJ1717	37LJ1532						
Material #	Product Description					GTIN				
NN-8519	19x28 Little LOTTA MHP					00816744025199				
	Lot numbers									
	37BC7418	37FD0159	37HC4981	37JI1556	37LC1691					

Material #	Product Description					GTIN				
NN-8520	19x21.5 PaediScope MHP					00816744025205				
	Lot numbers									
	37CC7851	37HD2770	37IC6879							
Material #	Product Description					GTIN				
NN-8521	15X26.5 LOTTA MHP					00816744025212				
	Lot numbers									
	37Bj1260	37Bj1832	37CB4144	37CD5653	37DA9528	37HC4984	37HD2772	37Hj1104	37Jl1820	
	37Jj0611	37LC1231								

Business Reply Form

NICO Myriad™ Handpiece | NICO Myriad™ Illumination Pack | NICO BrainPath

October 22, 2025

Please check one of the options below and complete this form. Email the completed form to Instruments.Recalls@stryker.com. **RESPONSE IS REQUIRED.** If an option below is not selected, a follow up from a Stryker representative will occur.

- ☐ I, the customer, confirm my facility has no affected products onsite.
- ☐ I, the customer, choose to keep the impacted product(s). I confirm this notification is placed where the product is stored at my facility and product is marked with "CONTAINS LATEX."
- ☐ I, the customer, destroyed my affected products locally with my Stryker representative present. Refer to attached Certificate of Destruction.
- ☐ I, the customer, choose to return the following product(s) to Stryker for replacement:

Catalog Number	Product Description	Affected Lots	Qty of product (EA)

Form completed by:

Facility Name			
Facility Address			
Printed Name		Title	
Email		Phone	
Signature		Date	

Kindly confirm that you have selected one of the options above and completed the table, if applicable.

Note: Your signature indicates that you have received and understand the enclosed notification and that you have performed all actions requested.

If you have further distributed any affected product, please indicate to whom:

Facility Name		Contact Person	
Full Address			

For reporting any serious adverse events or product quality problems contact Stryker at 1-800-253-3210, ext. 5555 8:00 A.M. to 7:00 P.M. (Eastern Time), Monday – Friday or by email at inst.stryker.cs@stryker.com. In addition to contacting Stryker, healthcare professionals and consumers may report any potential quality problems or adverse reactions, or events associated with the use of this product may be reported to the FDA's MedWatch Adverse Event Reporting program either online, or by fax or phone. Online: www.fda.gov/Safety/MedWatch/HowToReport/default.htm Fax: (800) FDA-0178 Phone: (800) FDA-1088.

Stryker Corporation or its affiliates own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

Certificate of destruction

If your facility chooses to dispose affected products on site, please complete the form in its entirety prior to submitting it via email to Instruments.Recalls@stryker.com. It is required that a Stryker representative is present when destroying affected product. **To receive replacements for products destroyed, this form must be fully completed.**

Method of Destruction:

To destroy impacted products, please do the following for every impacted product

1. Remove the label and carefully open the product packaging.
2. Deactivate the device to prevent further use.
3. Dispose of both the product and its packaging in accordance with local regulations.

Products Destroyed:

Catalog Number Provide part number	Lot Number Provide lot number, one lot number per line	Total EA destroyed

Location of destruction:

Facility Name	
Facility Address	

Destruction completed by:

Printed Name		Title	
Email		Phone	
Signature		Date	

Destruction witnessed by (Stryker representative):

Printed Name		Title	
Email		Phone	
Signature		Date	